136045-12

- Español
- Português
- Nederlands
- English

For QuikRead® 101 Instrument

QuikRead® CRP

with prefilled cuvettes

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- Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. Clin Chem Lab Med 1999, 37, 109–113.
- van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. Baillière's Clinical Rheumatology 1994, 531–552.
- Olaison L, Hogevik H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. Arch Intern Med 1997, 157, 885–892.
- Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. Pediatrics, 1997, 99, 846–850.
- Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. Pediatrics, 2000, 106, 1–5.
- Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. Oxford Textbook of Medicine, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
- Seamark DA, Backhouse SN, Powell R. Field-testing and validation in a primary care setting of a point-of-care test for C-reactive protein. Ann Clin Biochem 2003;40:178–180.
- Thomas L. C-reaktives Protein (CRP) [C-reactive protein (CRP). In German.]. In: Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft; 1992. p. 780–790.

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





Referenties • References

- Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med* 1999, 37, 109–113.
- van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillière's Clinical Rheumatology* 1994, 531–552.
- Olaison L, Hogevik H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med* 1997, 157, 885–892.
- Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics*, 1997, 99, 846–850.
- Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics*, 2000, 106, 1–5.
- Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
- Seamark DA, Backhouse SN, Powell R. Field-testing and validation in a primary care setting of a point-of-care test for C-reactive protein. *Ann Clin Biochem* 2003;40:178–180.
- Thomas L. C-reaktives Protein (CRP) [C-reactive protein (CRP)]. In German. In: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft; 1992. p. 780–790.

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1 Intended use

For quantitative determination of CRP (C-reactive protein) in whole blood, serum or plasma, using the QuikRead® 101 instrument. For *in vitro* diagnostic use.

2 Summary and explanation of the test

CRP is an acute phase protein present in low concentrations in healthy individuals¹. Any pathological condition associated with invasive bacterial infection, inflammation or tissue destruction is accompanied by elevation of the CRP level in the patient's serum. The rise in CRP levels is rapid, and increased levels can be detected within 6 to 12 hours of the onset of the inflammatory process². Quantitative measurement of the CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections, as well as an effective tool in controlling and monitoring postoperative infections²⁻⁶.

3 Principles of the procedure

QuikRead CRP is an immunoturbidimetric test based on microparticles coated with anti-human CRP F(ab)₂ fragments. The CRP present in the sample reacts with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead 101 instrument. The sample is added to the assay buffer. The blood cells in whole blood sample are haemolysed. The assay is performed in the same cuvette. The reagents are precalibrated, and the calibration curve, which is lot-specific, is coded on a magnetic card provided with each kit. QuikRead CRP correlates well with results obtained using immunoturbidimetric and other standard laboratory methods⁷.

4 Reagents

Content of the kit

Component name and origin	Symbol	QuikRead® CRP, Cat. no. 134191 50 assays	QuikRead® CRP, Cat. no. 134192 50 assays
QuikRead® CRP Reagent Caps ORIG SHP	REAG CPS	2 x 25	2 x 25
Buffer in prefilled cuvettes	BUF	2 x 25 x 1 ml	2 x 25 x 1 ml
Capillaries (20 µl)	CAPIL HEP	50	–
Plungers	PLUN	50	–
Magnetic card	MGN CRD	1	1
Instructions for use			

Reagents contain preservatives, please see Section 5 "Warnings and precautions".

Storage

Kit component	Storage at 2...8°C	Storage at 18...25°C
Reagent caps (in opened and unopened aluminium tube)	Until the expiry date of the kit	24 h per day – 1 month 7.5 h per day – 3 months
Prefilled cuvettes in unopened foil pouches	Until the expiry date of the kit	Until the expiry date of the kit
Prefilled cuvettes without the foil pouch	6 months	3 months
Opened prefilled cuvette	2 h	2 h

Mark the date of opening the foil pouch on the cuvette rack.

Reagent preparation and storage conditions

All reagents are ready to use. Keep the CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

5 Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water.
- All patient samples and controls should be handled as potentially infectious material.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagents contain < 1% sodium azide, which may be harmful (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412), avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Chapter 14

Warning

Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches (2 pcs) protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using and individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.

- QuikRead CRP Reagent Caps are colour-coded blue to distinguish them from other QuikRead analytes.
- Keep the QuikRead CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Make sure that the correct magnetic card for the QuikRead CRP test is used. The magnetic card must not be bent or exposed to liquids or strong magnetic fields.
- Do not splash liquid into the measuring well of the instrument.

6 Sample collection and preparation

6.1 Fingertip whole blood

Prick the finger with a lancet and discard the first drop. Wipe the finger and collect 20 µl of blood from the second drop in a glass capillary. Perform the test as outlined in Section 7 "Procedure".

Use of capillaries and plungers

The QuikRead CRP kit Cat. No. 134191 contains 20 µl glass capillaries and plungers, which are recommended for adding the sample. Insert the plunger into the capillary at the end with the blue stripe. Fill the capillary with sample up to the white stopper. Ensure that there are no air bubbles in the capillary. Use a soft tissue to wipe excess sample from the outer surface of the capillary. Make sure the soft tissue does not draw any sample from inside the capillary. Put the end of the capillary in the buffer and empty the capillary into the buffer by pressing down the plunger. Make sure the capillary is emptied completely.

6.2 Alternative sample materials

Anticoagulated whole blood: Whole blood collected in a tube containing heparin or EDTA can be used. Mix the whole blood by inverting the tube. Perform the test as outlined in Section 7 "Procedure". The sample volume is 20 µl. The 20 µl capillaries are recommended for adding the sample.

Serum: Take a normal venous blood sample and separate the serum. Perform the test as outlined in Section "Procedure". Sample volumes of 12 or 20 µl can be used, noting the effect of sample volume as stated below.

Plasma: Collect whole blood in a tube containing an anticoagulant such as heparin or EDTA. Separate the plasma from the blood cells as soon as possible to avoid haemolysis. Mild haemolysis in samples will not affect the test results. Perform the test as outlined in Section 7 "Procedure". Sample volumes of 12 or 20 µl can be used, noting the effect of sample volume as stated below.

Sample volume for serum or plasma samples

- 1) Pipette 12 µl of the plasma or serum sample. The final result appears in the display.
OR
- 2) Pipette 20 µl of the plasma or serum sample and multiply the result shown in the display by 0.6.

Sample storage

Sample material	Short-term storage	Long-term storage
Anticoagulated whole blood	2...8°C for 3 days	Separate plasma and store below -20°C
Serum	2...8°C for 7 days	Below -20°C
Plasma	2...8°C for 7 days	Below -20°C

Samples should not be repeatedly frozen and thawed. Allow the samples to reach room temperature (18...25°C) before testing. Frozen samples must be thawed completely, mixed thoroughly, brought to room temperature and centrifuged prior to testing.

6.3 Control materials

QuikRead CRP Control: QuikRead CRP Control is ready-to-use. Perform the test as outlined in Section 7 "Procedure". The sample volume is 20 µl. The 20 µl capillaries are recommended for adding the sample.

Other commercially available controls: Handle and use according to the instructions for each control substance. Serum- and plasma-based controls should be used in the same way as plasma and serum samples, see Section 6.2.

7 Procedure

Materials required but not supplied

	Cat. no.	
QuikRead 101 Instrument	06078	
Capillaries (20 µl) 50 pcs	67692	Supplied with 134191
Plungers 50 pcs	67966	Supplied with 134191
QuikRead CRP Control	154812	Recommended for quality control
Fingertip lancets for the whole blood procedure		
A pipette for dispensing 12 µl volumes		For plasma and serum samples (see Section 6.2)

Assay procedure

All reagents should reach room temperature (18...25°C) before use. Upon removal from refrigerated storage, discard the foil pouch of cuvettes and allow to stay at room temperature before use. It will take 15 minutes for an individual refrigerated prefilled cuvette to reach room temperature.

The measurement time of the assay depends on the instrument software version. With the software version 7.0 or newer, the assay reaction time is max. 1 minute. With older software versions the assay reaction time is max. 2 minutes.

The subheadings correspond to messages displayed by the instrument.

Read the card

Read the magnetic card by pulling it through the reader slot with the magnetic band facing you.

Ready for use CRP

- 1) Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on results. The test must be performed within two hours from opening the cuvette.
- 2) Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
Add 20 µl of undiluted whole blood sample to the cuvette, using a capillary. For plasma and serum sample volumes and the use of capillaries, see Section 6 "Sample collection and preparation".
- 3) Close the cuvette tightly with a CRP Reagent Cap. Put the cap properly on the cuvette to prevent leakage from the cuvette. Do not press down the blue inner part of the cap. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Mix by shaking gently (do not turn the cuvette upside down). Wait until the whole blood sample solution becomes bright red upon haemolysis. Insufficient haemolysis may cause the error message "Unstable sample. Please run again" to be displayed. If not measured immediately, the sample may be stored in the buffer for two (2) hours. Keep the cuvette in an upright position, do not shake. Before measuring the blank, mix again.
- 4) Put the cuvette in the measurement well of the instrument.


Measuring blank

The instrument measures the sample blank. This will take a maximum of 40 seconds. Do not touch the cuvette in the measurement well during the measurement. In case the error message "Unstable sample. Please run again" is displayed, repeat the blank measurement with the same cuvette.

**Blank measured
Add reagent**

OR

**Add reagent
Lift the cuvette**

- 1) Press down the blue inner part of the cap (using, for example, your finger or a pen). This will release the dried CRP reagent into the cuvette.
- 2) Take the cuvette from the measurement well and mix the contents very rapidly by shaking the cuvette vigorously. The whole blood sample solution in the cuvette becomes opaque red. Foaming of the mixture does not disturb the measurement.

Shake the cuvette

The display shows you how long you should shake the cuvette (approx. 6 seconds). The moving dots indicate the pace of vigorous shaking. Too slow shaking may cause the error message "Faulty reagent addition. Please run again" to be displayed. In that case a new test needs to be performed.

Insert the cuvette for measurement

Place the cuvette in the measurement well of the instrument. If the cuvette is inserted back too early, the error message "Test cancelled" or "Cuvette inserted too early" will be displayed. If the cuvette is reinserted too late, the error message "Cuvette inserted too late" will be displayed. In both cases, a new test needs to be performed.

Measuring 55 s

The instrument measures the CRP concentration within the time shown on the display. The display shows the time left before completion. Do not remove the cuvette from the measurement well before the result is displayed and recorded. Cuvette removal will terminate the measurement immediately.

Sample: <XXX> / date
Result: <XXX> CRP

The result is shown in the display as mg/l.

Performing self test Wait....

When the cuvette is removed, the QuikRead 101 instrument automatically performs a self-check. If the self-check is not successful, discard the previous test result. After a successful self-check the LCD will display:

Ready for use CRP

To begin a new test, place a cuvette in the measurement well.

Brief instructions for use

1	Remove the foil cover from the cuvette. Beware of splashing the liquid.
2	Add 20 µl of whole blood. Put the cap on and mix gently (not upside down).
3	Blank measurement.
4	Add reagent by pressing down the inner part of the cap.
5	Take out the cuvette and mix thoroughly. Put the cuvette back in the measurement well.
6	The instrument measures CRP concentration.
7	Read the CRP result.

Sample dilution

QuikRead CRP test gives numerical value up to 160 mg/l. To obtain numerical value for higher concentrations, dilute the sample with 0.9 % NaCl. Dilution ratio is 1+1 (1 part of volume sample and 1 part of volume 0.9 % NaCl). Take 20 µl of the diluted sample into a new cuvette, repeat the assay and multiply the result by 2.

8 Quality control

Regular use of QuikRead CRP Control (Cat. No. 154812) is recommended. The control value is specified for a volume of 20 µl with the same instructions for use as for whole blood samples.

9 Results' interpretation

Increases in CRP are nonspecific and should not be interpreted without a complete clinical history.

Elevated CRP range

Concentrations of < 10 mg/l exclude many acute inflammatory diseases but do not specifically exclude inflammatory processes. Elevated concentrations of < 50 mg/l in acute disease occur in the presence of slight to moderate inflammatory processes. Values > 50 mg/l indicate high and extensive inflammatory activity⁸.

Effect of haematocrit for whole blood samples

The calibration of the reagents for whole blood has been made assuming that the value for haematocrit is approximately 40 % and a plasma proportion is 60 %. If the actual haematocrit value deviates from 40 %, the corrected CRP value can be calculated by multiplying the whole blood result by the factor given below:

Haematocrit (%)	Factor	Haematocrit (%)	Factor	Haematocrit (%)	Factor
20-29	0.8	59-61	1.5	73	2.2
30-36	0.9	62-63	1.6	74	2.3
37-42	1.0	64-65	1.7	75	2.4
43-47	1.1	66-67	1.8	76-77	2.6
48-51	1.2	68-69	1.9	78	2.7
52-55	1.3	70	2.0	79	2.9
56-58	1.4	71-72	2.1	80	3.0

10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

11 Expected values

Measurement range: 8–160 mg/l

Reference values: < 10 mg/l

It is recommended that each laboratory determine their own reference values according to the health and age characteristics of their target populations.

12 Performance characteristics

Precision

Precision study was performed by the manufacturer according to Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2.

Within-run, between-day and total precision

Sample	Number of days	Mean CRP (mg/l)	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Sample 1	20	14	5.1	2.4	6.3
Sample 2	20	77	4.9	0.8	5.6
Sample 3	20	142	5.0	1.3	5.2
Sample 4 (whole blood)	20	51	1.7	2.5	3.5

Interference

Interfering substance	Concentration	Interference
Rheumatoid factors (RF)	≤ 525 IU/ml	none
Bilirubin	≤ 400 μmol/l	none
Triglycerides**	≤ 10 mmol/l (Intralipid)	none

** The instrument displays the error message "Sample blank too high. Please run again" if the sample is too lipaemic. It is recommended that highly lipaemic or precipitated samples be clarified by centrifuging (e.g. 10 min at 15000 x g) before assaying.

EDTA and heparin as anticoagulants in the samples do not interfere with the test.

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test, as the assay antibodies lack the Fc-part. In rare cases, interference by IgM myeloma protein has been observed.

Antigen excess

CRP concentrations less than 600 mg/l do not give falsely low results.

Whole blood versus plasma comparability

In a comparison involving 180 patient samples, the results for whole blood and plasma were found comparable.

$$y \text{ (whole blood)} = 1.04 \times \text{(plasma)} - 0.32$$
$$r = 0.99$$

13 Traceability

The calibrators used to calibrate the CRP assay of the QuikRead CRP test are traceable to the ERM®-DA474 reference material.

14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, used caps, cuvettes, capillaries and plungers should be handled and disposed of as potentially infectious material.
- Materials of the components:
 - Paper: Instructions for use
 - Cardboard: Kit box
 - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, plunger and capillary tubes
 - Glass: Capillaries
 - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
 - Several (not to be recycled): Lids of reagent cap tubes
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

15 Troubleshooting

Problem	Possible cause	Corrective action
Unexpected low results	Too small sample volume.	Retest the sample. Make sure that the capillary is completely filled. Avoid air bubbles.
	Incorrect reagent storage.	Retest the sample. Make sure that all the reagents are stored according to the instructions for use.
	Components from different kit lots or tests are used.	Retest the sample. Make sure that all the reagents are from the same kit lot.
Unexpected high results	Too large sample volume.	Retest the sample. Wipe excess sample from the outer surface of the capillary. Make sure that the sample is taken at the end with the white stopper and the plunger is inserted at the end with the blue stripe.
	A sample volume of 20 μl was used for plasma or serum sample but the result was not multiplied by the correction factor.	Multiply the result on the display by 0.6.
	The cuvette is dirty.	Retest the sample. Do not touch the clear flat surfaces in the lower part of the cuvette.
	Incorrect reagent storage.	Retest the sample. Make sure that the reagents are stored according to the instructions for use.
	Components from different kit lots or tests are used.	Retest the sample. Make sure that all the reagents are from the same kit lot.

For troubleshooting, see also the QuikRead 101 instrument instructions for use.

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QuikRead® is a registered trademark of Aidian Oy.