132515-21

- English
- Deutsch
- Français
- Suomi
- 中國

For QuikRead® 101 Instrument

# QuikRead® U-ALB

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- 1) Viberti GC, Mogensen CE, Passa P et al. St Vincent declaration, 1994: Guideline for the prevention of diabetic renal failure.
- 2) Bigazzi R, Bianchi S, Baldari D and Campese V. Microalbuminuria predicts cardiovascular events and renal insufficiency in patients with essential hypertension. J Hypertension, 1998; 16, 1325–1333.
- 3) Mogensen CE, Keane WF, Bennett PH et al. Prevention of diabetic renal disease with special reference to microalbuminuria. Lancet, 1995; 346, 1080–1084.

QuikRead® is a registered trademark of Aidian Oy.  
 QuikRead® ist eingetragenes Warenzeichen von Aidian Oy.  
 QuikRead® est une marque déposée de Aidian Oy.  
 QuikRead® on Aidian Oy:n rekisteröimä tavaramerkki.

# AIDIAN

Aidian Oy  
 Koivu-Mankkaan tie 6 B  
 P.O. Box 83, FI-02101 Espoo, Finland  
[www.aidian.eu](http://www.aidian.eu)  
[www.quikread.com](http://www.quikread.com)

## 1 Intended use

For the quantitative determination of albumin in urine (U-ALB), using the QuikRead® 101 instrument. For *in vitro* diagnostic use.

## 2 Summary and explanation of the test

Nephropathy is one of the most common complications of diabetes, and it is internationally recommended that diabetics should be monitored for incipient nephropathy<sup>1</sup>. Since microalbuminuria is suggestive of incipient nephropathy, diabetics should have urinary albumin (U-ALB) measurements a few times a year. The prompt detection of microalbuminuria is important because early intervention can halt and even reverse the progression of nephropathy. As hypertensive individuals may also develop microalbuminuria, they too benefit from regular monitoring of urinary albumin concentration<sup>2</sup>.

## 3 Principles of the procedure

QuikRead U-ALB is an immunoturbidimetric test based on microparticles coated with anti-human albumin. The albumin present in the sample reacts with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead 101 instrument. Urine sample is added to the buffer and the assay is performed in the same cuvette. The reagents are precalibrated, and the calibration curve, which is batch specific, is coded onto a magnetic card provided with each kit.

QuikRead U-ALB correlates well with results obtained using immunoturbidimetric methods.

## 4 Reagents

### Contents

Component name and origin	Symbol	QuikRead® U-ALB Cat. No. 06014 50 assays
QuikRead® U-ALB Albumin-Reagent Caps ORIG   PIG	REAG   CPS	50
Buffer	BUF	1 x 120 ml
Cuvettes	CUVET	50
Magnetic card	MGN   CRD	1
Instructions for use	-	1

Reagents contain sodium azide, please see Section 5 "Warnings and precautions".

### Reagent storage

Kit component	Storage at 2...8°C	Storage at 18...25°C
Reagent caps	Until the expiry date of the kit	24 h per day – 1 month 7.5 h per day – 3 months
Buffer	Until the expiry date of the kit	Until the expiry date of the kit

### Reagent preparation and storage conditions

All reagents are ready to use. Keep the Albumin-Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

### Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

## 5 Warnings and precautions

### Health and safety information

- For *in vitro* diagnostic use.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water.
- All patient samples and controls should be handled as potentially infectious material.
- Lyophilised reagents contain <1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501).
- Reconstituted and liquid reagents contain <0.1% sodium azide, which is not considered a harmful concentration.
- Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Chapter 14.

### Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- QuikRead U-ALB Reagent Caps are colour-coded green to distinguish them from other QuikRead analytes.
- Keep the QuikRead U-ALB Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Make sure that the correct magnetic card for the QuikRead U-ALB test is used. The magnetic card must not be bent or exposed to liquids or strong magnetic fields.
- Make sure that the dispenser is properly cleaned when taking a new reagent kit into use.
- Always check before use that there are no air bubbles in the dispenser tubes. Remove any air bubbles by dispensing the buffer once before actual use.
- Do not splash liquid into the measuring well of the instrument.

## 6 Sample collection and preparation


The QuikRead U-ALB test is performed with urine. Adding preservatives is not recommended. If the test cannot be performed immediately, the urine may be stored at 2...8°C for 14 days. Urine samples should not be frozen. Turbid samples can be centrifuged before the assay is performed.

<b>Urine collection overnight</b>	Urine is collected overnight, preferably for at least six hours. The duration of collection and the volume of urine collected are recorded. The albumin excretion rate is calculated according to the formula in Section 12.
<b>Spot urine sample</b>	For screening of microalbuminuria a spot sample of urine can be used. With spot samples of urine, the most reliable results are obtained using the first morning urine.

Allow the samples to reach room temperature (18...25°C) before testing.

### Use of capillaries and plungers

Insert the plunger into the capillary from the end with the blue stripe. Fill the capillary with sample up to the white stopper. Ensure that there are no air bubbles in the capillary. The volume of the capillary is 20 µl. Use a soft tissue to wipe excess sample from the outer surface of the capillary. Make sure the soft tissue does not suck any sample from the inside of the capillary. Put the end of the capillary in the buffer and empty the capillary into the buffer by pressing down the plunger. Make sure the capillary is emptied completely.

### Use of QuikRead 1 ml Dispenser

The dispenser is a 1 ml accuracy dispenser. Put the dispenser on the buffer bottle and close it tightly. Ensure the functioning of the dispenser by dispensing a few times and discarding the dispensed buffer. The buffer bottle contains an excess of buffer in order that the functioning of the dispenser may be ensured. Press the plunger button to the bottom in one smooth stroke. Avoid letting the button to "snap" back. One stroke corresponds to 1 ml of buffer. To avoid contamination, use the dispenser only with QuikRead U-ALB. The dispenser is delivered with labels indicating the tests for which the dispenser can be used. After unpacking the dispenser, stick the label indicating the test you are going to use the dispenser for, onto the dispenser.

See the instructions manual of the dispenser for more detailed information.

## 7 Procedure

### Materials required but not supplied

	Cat. No.	
QuikRead 101 Instrument incl. 1 ml Dispenser or	06040	Recommended for quality control
QuikRead 101 Instrument without Dispenser +	06078	
Buffer dispenser (1 ml)	05560	
Capillaries (20 µl), 50 pcs	67962	
Plungers, 50 pcs	67966	
Urine Albumin Control	67881	

### Assay procedure

All reagents should be allowed to reach room temperature (18...25°C) before use. Use the samples undiluted.


**The subheadings correspond to messages displayed by the instrument.**

#### Read the card

Read the magnetic card by pulling it through the reader slot with the magnetic band facing you.

#### Ready for use U-ALB

- 1) Dispense 1 ml of buffer into the cuvette, using the QuikRead Dispenser. You can make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
- 2) Add 20 µl of undiluted urine sample into the cuvette, using a capillary.
- 3) Close the cuvette tightly with the Albumin-Reagent cap and mix by shaking gently (not upside-down). Do not press the inner part of the cap. Do not touch the clear flat surfaces on the lower part of the cuvette (the optical part).
- 4) Put the cuvette into the measurement well of the instrument.

#### Measuring blank


The instrument measures the sample blank. This will take a maximum of 30 seconds. Do not touch the cuvette in the measurement well during the measurement. In case the error message "Unstable sample. Please run again" is displayed, repeat the blank measurement with the same cuvette.

#### Blank measured Add reagent

OR

#### Add reagent Lift the cuvette

- 1) Press down the inner part of the cap (using, for example your finger or a pen). This will release the dried Albumin-Reagent into the cuvette.
- 2) Take the cuvette from the reading well and mix very rapidly by shaking the cuvette vigorously. The solution in the cuvette becomes milky white.

#### Shake the cuvette .....

The display shows you how long you should mix the cuvette (approx. 5 seconds). The moving dots indicate the pace of vigorous shaking. Too slow shaking may cause the error message "Faulty reagent addition. Please run again" to be displayed. In that case a new test needs to be performed.

#### Insert the cuvette for measurement

Place the cuvette in the measurement well of the instrument. If the cuvette is inserted back too early, the error message "Test cancelled" or "Cuvette inserted too early" will be displayed. If the cuvette is reinserted too late, the error message "Cuvette inserted too late" will be displayed. In both cases, a new test needs to be performed.

#### Measuring 120 s

The instrument measures the albumin concentration in 2 minutes. The display shows the time left before completion. Do not remove the cuvette from the measurement well before the result is given and recorded. Cuvette removal will terminate the measurement immediately.

Sample: <XXX> / date  
Result: <XXX> U-ALB

The result is shown in the display as mg/l.

#### Performing self test Wait...

When the cuvette is removed, the QuikRead 101 instrument automatically performs a self-check. If the self-check is not successful, discard the previous test result. After a successful self-check the LCD will display:

#### Ready for use U-ALB

To begin a new test, place a cuvette in the measurement well.

### Brief instructions for use

1	Dispense 1 ml of buffer into the cuvette.
2	Add 20 µl of urine. Put the cap on and mix gently (not upside-down).
3	Blank measurement max. 30 seconds.
4	Add reagent by pressing down the inner part of the cap.
5	Take out the cuvette and mix thoroughly. Put the cuvette back in the measurement well.
6	Instrument measures albumin concentration during 120 seconds.
7	Read the U-ALB result.

### Sample dilution

QuikRead yields a numerical reading up to 150 mg/l. To obtain numerical exact readings for higher concentrations, the sample must be diluted and reanalysed.

For instance, predilution 1:5 (1+4) or 1:25 (1+24) with 0.9 % NaCl or QuikRead U-ALB buffer can be used. Analyse the diluted sample and multiply the result with the dilution factor.

Examples:

Dilution 1:5	Dilution 1:25
1) Take 250 µl of undiluted urine.	1) Take 100 µl of undiluted urine.
2) Add 1 ml of 0.9 % NaCl or buffer.	2) Add 2400 µl of 0.9 % NaCl or buffer.
3) Mix the diluted sample.	3) Mix the diluted sample.
4) Perform the assay in the same way as with undiluted samples.	4) Perform the assay in the same way as with undiluted samples.
5) Multiply the albumin result by 5.	5) Multiply the albumin result by 25.

## 8 Quality control

We recommend regular use of Urine Albumin Control (Cat. No. 67881). The target value for QuikRead U-ALB is indicated on the label of the control vial.

## 9 Results' interpretation

Microalbuminuria is defined as an overnight urinary excretion of 20–200 µg albumin/min or an albumin concentration of 20–200 mg/l in a spot urine sample. The diagnosis of microalbuminuria is confirmed if at least two out of three samples are positive over a period of three to six months<sup>3</sup>.

### Calculation of albumin excretion rate with timed urine collection

The QuikRead 101 instrument displays albumin concentration as mg/l (= µg/ml). The albumin excretion rate for a timed urine sample (µg/min) can be calculated using the following formula:

$$\frac{\text{total urine volume (ml)}}{\text{collection period (min)}} \times \text{QuikRead value of U-ALB (mg/l)} = \text{albumin excretion rate (µg/min)}$$

## 10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics". U-ALB results should never be used alone, without a complete clinical evaluation. Intra- individual variations in Urine Albumin are significant and should be taken into account – e.g., by means of serial measurements – in interpretation of the values.

## 11 Expected values

- Measurement range: 5–150 mg/l
- Reference values:
  - Overnight collection: < 20 µg/min
  - Spot urine: < 20 mg/l

## 12 Performance Characteristics

The following performance characteristics have been determined for the QuikRead 101 instrument.

### Precision and Reproducibility

An example of typical results:

Intra-assay precision (n=10)			Inter-assay precision (n=10)		
Urine samples	Mean (mg/l)	CV (%)	Urine samples	Mean (mg/l)	CV (%)
Sample 1	10	4.7	Sample 1	11	6.6
Sample 2	56	1.2	Sample 2	57	1.5
Sample 3	136	0.5	Sample 3	135	2.3

Detection limit: < 5 mg/l (2.1 mg/l)

Linearity: 89–100 % (at 132 mg/l)

### Interferences and limitations

No interference with the test has been found with the following substances:

Acetone	< 2.5 g/l	Glucose	< 50 mmol/l
Ascorbic acid	< 10 g/l	Hemoglobin	< 500 mg/l
β-2-microglobulin	< 250 mg/l	IgG	< 500 mg/l
Bilirubin	< 100 µmol/l	IgA	< 180 mg/l
Citric acid	< 6 mmol/l	Myoglobin	< 100 mg/l
Creatinine	< 60 mmol/l	Sodium chloride	< 250 mmol/l
Fibronectin	< 4 mg/l	Urea	< 500 mmol/l

As human serum contains large amounts of albumin, the presence of even small quantities of whole blood in urine samples give elevated values in U-albumin tests.

### Antigen excess

The instrument has a built-in antigen-excess detection system. As a result, albumin concentrations exceeding 150 mg/l and up to 40 000 mg/l are displayed as ">150 mg/l".

## 13 Traceability

The calibrators used to calibrate the QuikRead U-ALB test are traceable to the ERM®-DA470 reference material.

## 14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, used cuvettes, capillaries and plungers should be handled and disposed of as potentially infectious material.
- Materials of the components
  - Paper: Instructions for use
  - Cardboard: Kit box including its inner parts
  - Plastic: Cuvettes, reagent caps, plastic covering cuvette rack, cuvette rack, buffer bottle, plungers, plunger and capillary tubes
  - PVC: Magnetic card
  - Glass: Capillaries
  - Metal: Reagent cap tubes, plunger and capillary tube caps
  - Several (not to be recycled): Lids of reagent cap tubes
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

## 15 Troubleshooting

For troubleshooting, see the QuikRead® 101 instrument instructions for use.

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