

# QuikRead go Plus Instrument

**155448-1**

- English
- Deutsch
- Français
- Nederlands

	Explanation of symbols	Zeichenerklärung	Explication des symboles	Verklaring van symbolen
	<i>In vitro</i> diagnostic medical device	<i>In-vitro</i> -Diagnostikum	Dispositif médical de diagnostic <i>in vitro</i>	Medisch hulpmiddel voor <i>in-vitro</i> -diagnostiek
	Device for near-patient testing	Produkt für patientennahe Tests	Dispositif de diagnostic près du patient	Hulpmiddel voor patiëntnabije tests (hulpmiddel voor near-patient testing)
	Serial number	Seriennummer	Numéro de série	Serienummer
	Catalogue number	Bestellnummer	Référence du catalogue	Catalogus nummer
	Manufacturer	Hersteller	Fabricant	Fabrikant
	Date of manufacture	Herstellungsdatum	Date de fabrication	Productie datum
	Contents	Inhalt	Contenu	Inhoud
	Instrument	Gerät	Appareil	Instrument
	Power supply	Netzteil	Alimentation	Stroomadapter
	Mains cable	Netzkabel	Câble d'alimentation	Netkabel
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter la notice d'utilisation	Raadpleeg de gebruiksaanwijzing
	Caution	Achtung	Attention	Voorzichtig
	Temperature limitation	Temperaturbegrenzung	Limites de température	Temperatuurlimiet
	Fragile, handle with care	Zerbrechlich, mit Sorgfalt behandeln	Fragile, manipuler avec précautions	Breekbaar, voorzichtig behandelen
	Keep dry	Vor Nässe schützen	Conserver au sec	Droog houden
	China RoHS GB/T 26572	China RoHS GB/T 26572	China RoHS GB/T 26572	China RoHS GB/T 26572
	Waste electrical and electronic equipment	Elektro- und Elektronik-Altgeräte	Déchets d'équipements électriques et électroniques	Afgedankte elektrische en elektronische apparatuur

	Explanation of symbols	Zeichenerklärung	Explication des symboles	Verklaring van symbolen
	This product fulfills the requirements of the European Parliament and of the Council on <i>in vitro</i> diagnostic medical devices	Dieses Produkt erfüllt die Anforderungen des Europäischen Parlaments und des Rates für <i>In-vitro</i> -Diagnostika	Ce produit répond aux exigences du Parlement européen et du Conseil sur les dispositifs médicaux de diagnostic <i>in vitro</i>	Dit product voldoet aan de eisen van het Europees Parlement en de Raad voor medische hulpmiddelen voor <i>in-vitro</i> diagnostiek
	Authorised representative in Switzerland	Autorisierter Vertreter in der Schweiz	Représentant autorisé en Suisse	Gemachtigde vertegenwoordiger in Zwitserland
	USB	USB	USB	USB
	RJ-45	RJ-45	RJ-45	RJ-45
	Power	Strom	Alimentation	Spanning
	On-Off switch	Ein-/Ausschalter	Interrupteur marche-arrêt	Aan/uit-schakelaar
	Restriction	Einschränkung	Restriction	Beperking

	Abbreviation	Abkürzung	Abréviation	Afkorting
HIS	Hospital information system	Krankenhausinformationssystem (Hospital information system)	Système d'information hospitalier	Informatiesysteem van het ziekenhuis
ID	Identifier	Kennung (Identifier)	Identifiant	Identificatienummer
IP	Internet protocol	Internet-Protokoll	Protocole Internet	Internetprotocol
LIS	Laboratory information system	Laborinformationssystem	Système d'information du laboratoire	Informatiesysteem van het laboratorium
POCT01	Point-of-Care Connectivity	Point-of-Care-Konnektivität	Connectivité point of care (POC)	Point-of-care-verbinding
QC	Quality Control	Qualitätskontrolle	Contrôle de qualité	Kwaliteitscontrole
SVHC	Substances of Very High Concern	Besonders besorgniserregende Stoffe	Substances extrêmement préoccupantes	Zeer zorgwekkende stoffen
USB	Universal Serial Bus	Universeller serieller Bus	Bus universel en série	Universele Seriele Bus

## About this Instructions of use

QuikRead go Plus Instrument instructions for use will guide you in installation, use, and maintenance of your instrument. The instructions for use also explains the testing principle, the quality control functions, and aids in troubleshooting.

Before analysing patient samples or quality control samples, please read the instructions for use carefully and all QuikRead go reagent instructions for use documents supplied with each reagent kit.

## Symbols

QuikRead go Plus Instrument packaging and device shield symbols, please see page 2.

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
# CONTENTS

<b>1 INTRODUCTION .....</b>	<b>7</b>	Quality control sample measurement .....	16	Self test .....	23
Intended purpose .....	7	QC result view .....	16	Patient result management .....	23
QuikRead go Plus Instrument .....	7	Detailed QC result view .....	16	Delete patient result list .....	23
Safety information .....	7	Results .....	17	Factory reset .....	23
Precautions and limitations .....	7	Deleting results .....	17	<b>7 ABOUT VIEW .....</b>	<b>24</b>
<b>2 GETTING STARTED .....</b>	<b>8</b>	LIS offline list .....	17	<b>8 INSTRUMENT MAINTENANCE ...</b>	<b>24</b>
Unpacking .....	8	<b>4 SETTINGS .....</b>	<b>18</b>	Instrument calibration .....	24
Accessories .....	8	Language, display, sound .....	19	Cleaning the instrument .....	24
QuikRead go Plus Instrument parts .....	8	Date and time .....	19	Software update .....	24
Lifting or transporting the QuikRead go Plus Instrument .....	10	Power saving .....	19	Changing the clock battery .....	24
Location and environment .....	10	Measurement .....	19	<b>9 INSTRUMENT SPECIFICATION ..</b>	<b>24</b>
During use .....	10	Test parameters .....	19	Declaration of conformity .....	24
During transport and storage .....	10	Accessory .....	20	Technical specification .....	24
Power supply cable and battery .....	10	User management .....	20	Photometer .....	24
Connectors and cables .....	10	Operator settings .....	20	Touchscreen display .....	24
Plugging the power supply cable .....	10	Patient settings .....	20	Dimensions & power requirements .....	24
Inserting the rechargeable battery .....	10	Additional ID settings .....	20	Instrument identifier .....	24
Power ON, OFF or sleep mode .....	10	Connectivity .....	21	Memory .....	24
Turning power on .....	10	LIS .....	21	Power supply .....	25
Turning power off .....	10	TCP/IP .....	21	Radio module .....	25
Use of the touch screen .....	10	WLAN .....	21	LIS connection .....	25
First-time setup .....	11	<b>5 QC MANAGEMENT .....</b>	<b>22</b>	USB connection .....	25
User interface .....	12	QC settings .....	22	Service .....	25
Status area icons .....	12	Add new QC .....	22	Warranty .....	25
Menu .....	13	QC list .....	22	Disposal .....	25
Notifications .....	13	Add new QC vial in use .....	23	Revision history .....	25
Language selection .....	13	<b>6 MAINTENANCE MENU .....</b>	<b>23</b>	<b>10 TROUBLESHOOTING .....</b>	<b>26</b>
General settings .....	14	Software update .....	23	User role rights .....	29
<b>3 INSTRUMENT OPERATION .....</b>	<b>14</b>	Logs .....	23		
Performing a measurement .....	14	Error log .....	23		
Patient sample measurement .....	14	Audit logs to USB .....	23		
Patient sample result view .....	15	Measurement logs to USB .....	23		
Detailed sample result view .....	15	Connectivity logs .....	23		

# 1. INTRODUCTION

## Intended purpose

QuikRead go Plus Instrument is an easy-to-use *in vitro* diagnostic test system. It has been designed to measure various analytes from patient samples needed as an aid in diagnosis and treatment monitoring. The system consists of QuikRead go Plus instrument and QuikRead go reagent kits.

QuikRead go Plus Instrument is an automated instrument designed and calibrated for both photometric and turbidimetric measurements. The instrument is intended for quantitative and qualitative determination of various QuikRead go reagent kit analytes from human samples such as whole blood, serum, plasma, throat swabs and faecal samples to be used as an aid in diagnosis and treatment monitoring. QuikRead go Plus Instrument is intended to be used by healthcare professionals in the clinical laboratory and near-patient testing setting.

## QuikRead go Plus Instrument

The instrument will guide you through the assay procedure by using a series of messages and animations shown on the screen. During start up, the instrument performs a self test to ensure proper functionality of the instrument's internal components such as the optical system and mechanical parts.

QuikRead go Plus Instrument measures the absorbance or turbidity of the contents of a cuvette and converts the value into an analyte concentration value or into a qualitative result based on preset test calibration data. The calibration data for each test is encoded on the label of the cuvette. This information is transferred to QuikRead go Plus Instrument automatically during a measurement.

The assays are performed according to the instructions for use accompanying each QuikRead go reagent kit. The results are available in minutes. The instrument can be operated on mains power or with a battery, it has USB connections for an external printer, a keyboard, and a barcode reader.

QuikRead go Plus Instrument can be connected to a remote Laboratory and Hospital Information System (LIS/HIS). The instrument uses a standardized data transfer protocol. Contact your local supplier for more details.

## Safety information

For your safety, comply with all warning and caution statements. To alert you to potential electrical or operational hazards, warning and caution statements are provided where applicable. QuikRead go Plus Instrument contains Substances of Very High Concern (SVHC) as specified in REACH regulation (EU 1907/2006), for more information [aidian.eu](http://aidian.eu). The SVHC substances are bound to the inside components of QuikRead go Instrument. Specific precautions in handling are not required. Before taking QuikRead go Plus Instrument into use, please read the following precautions and limitations carefully. In case of a serious incident, please report it to the manufacturer or its representative and/or national authority.

## Precautions and limitations

- For *in vitro* diagnostic use only.
- Do not spill any liquids or drop any objects onto or into the instrument.
- Spillage of potentially infectious material should be wiped off immediately with an absorbent paper tissue and the contaminated areas swabbed with a standard disinfectant or 70% ethyl alcohol (see the Section "Cleaning the instrument").
- Material used to clean spills, including gloves, should be disposed of as biohazardous waste.
- Follow the QuikRead go reagent instructions for use supplied with each reagent kit.
- Only QuikRead go reagents can be used.
- Materials required but not provided are listed in QuikRead go reagent kit instructions for use.
- Do not mix components with different lot numbers or different assays.
- Never place a cuvette without a tightly placed reagent cap into the instrument.
- Do not insert fingers or any external devices into the instrument when the lid is open.

- Ensure that the cuvette sealing foil is completely removed.
- Pressing hard or using sharp items may damage the screen.
- Only use the power supply supplied with the instrument and make sure that the plug is placed so that it is removable.
- Only use the official QuikRead go Plus Instrument battery supplied by Aidian.
- Do not remove or shut down a USB storage during data transfer.
- Use USB cables shorter than 3 m to connect accessories to the instrument.
- Do not open the instrument's covers by opening any screws. If the warranty seal is broken, the instrument warranty is not valid (see **Picture 3**).
- Use a secure internal network or Virtual Private Network (VPN) when connecting the QuikRead go to a LIS/HIS using LAN.
- Do not use or connect the instrument to LAN if the warranty seal is broken.
- Do not leave an instrument unlocked in publicly accessible areas.
- Delete user accounts no longer in use from the instrument (including their related local data).
- The electromagnetic environment should be evaluated before the operation of the instrument.
- This equipment has been designed and tested to CISPR 11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

## 2. GETTING STARTED

### Unpacking

Open the package and check that it contains all of the following items:

- Instrument
- Instructions for use
- Power supply
- Power cable
- Certificate of analysis

Carefully examine the instrument to ascertain that it has not been damaged during shipment. If damage has occurred or any parts are missing, immediately notify your local supplier.

The instrument package can be stored to be used for transporting i.e., if sending to service.

### Accessories

QuikRead go system accessories make the system more effective, more secure, and easier to use. More information can be found at [aidian.eu](http://aidian.eu).

Available accessories:

#### Printer

- An external printer can be connected to the instrument. A list of compatible printers and configuration parameters can be found at [aidian.eu](http://aidian.eu).
- Connect the printer to a USB port. When a printer is connected the printer icon is shown in the status bar.

#### Barcode reader

- An external barcode reader can be connected to the instrument. A list of compatible barcode readers can be found at [aidian.eu](http://aidian.eu).
- Connect the barcode reader to a USB port. When a barcode reader is connected the barcode reader icon is shown in the status bar.


#### Keyboard

- An external keyboard can be connected to the instrument.
- Connect the keyboard to a USB port. When a keyboard is connected the keyboard icon is shown in the status bar.


**Note!** Do not use USB cords 3 m in length or longer to connect accessories to the instrument.

### QuikRead go Plus Instrument parts


The parts of the instrument are shown in **Picture 1**, **Picture 2**, and **Picture 3**.

**Picture 1** Front view  
1. Measurement well  
2. Touch screen  
3. Power button  
4. USB port (type A)

**Picture 2** Back view

1. RJ-45 port
2. Connector for power supply
3. USB port (type C)
4. USB port (type C)
5. USB port (type A)
6. USB port (type A)
7. Handle

RJ-45



Power

USB

Consult instructions for use

**Picture 3** Underside

1. Battery unit cover
2. Battery connector
3. Battery unit
4. Clock battery

## Lifting or transporting QuikRead go Plus Instrument

When lifting or transporting the instrument, always handle it carefully. To lift the instrument use the handle at the back of the instrument (**Picture 2**).

## Location and environment

### During use

QuikRead go Plus Instrument should be placed on a flat clean, horizontal surface, and the following should be observed:

- Indoor use
- Altitude up to 2000 m
- The ambient temperature must be between 15°C and 35°C.
- Maximum relative humidity 80 % for temperatures up to 31°C decreasing linearly to 67 % relative humidity at 35°C (non-condensing).
- Do not place the instrument in direct sunlight.
- Do not move the instrument or expose it to shaking during measurements.
- Do not measure in a moving vehicle.
- Mains supply voltage fluctuations up to  $\pm 10$  % of the nominal voltage.
- Place the instrument in such a way that it is easy to turn the power off and unplug the power supply cable.
- Do not charge the instrument with a power bank.
- Do not connect external devices not authorized by Aidian. EMC may be exceeding
- Do not insert LTE mobile data USB to the instrument
- All external circuits of devices connected to the equipment shall be provided with at least double insulation to mains.
- Do not place the instrument into a strong magnetic or electric field.
- Do not use this instrument in proximity to sources of strong electromagnetic radiation (e.g., unshielded intentional RF sources), as these can interfere with proper operation.
- Pollution degree 2.
- Installation category II (2500 V transient).

- Restrictions: This device WiFi or Bluetooth is not allowed to be used within a 20 km radius of the centre of Ny-Ålesund at Svalbard, Norway.

### During transport and storage

- The ambient temperature must be between 2°C and 35°C.
- Protect from rain and humidity.
- Handle the instrument with care.

## Power supply cable and battery

QuikRead go Plus Instrument can be used either with the power supply cable or the battery. The battery is charged automatically when the power cable is plugged in. List of compatible batteries can be found at [aidian.eu](http://aidian.eu).

### Connectors and cables

The instrument has a USB-A port on its side (**Picture 1**) and five connectors: two USB-A ports, two USB-C ports, and one LAN connector at the back of the instrument (**Picture 2**), and one USB-A port is placed on the right side of the instrument. The wiring diagram is described at [aidian.eu](http://aidian.eu).

### Plugging the power supply cable

Plug the power supply cable into the back of the instrument (see **Picture 2**). Plug the power supply into a power outlet.

### Inserting the rechargeable battery

Please, follow the steps below carefully when inserting the rechargeable battery into the instrument (see **Picture 3**).

1. Turn off the instrument (if on) and unplug the power supply cable.
2. Turn the instrument onto its side on an even surface and open the battery unit cover.
3. Attach the battery connector to the battery unit.

4. Press the battery into its place and ensure that it is placed properly.
5. Close the battery unit cover and turn the instrument back to an upright position.

## Power ON and OFF

### Turning power on

To turn the instrument on, press the power button on the side panel (see **Picture 1**) until the backlight of the screen is lit up. Light on the power button indicates that the instrument is powered on. If nothing happens, ensure that the power cable is connected. If the instrument is on battery usage and doesn't start by pressing the power button the battery might have run out. To start the instrument with an empty battery attach the power cable and press the power button.

### Turning power off

To turn the instrument off, press the power button. A pop-up to confirm shutdown will appear. If confirmed the instrument will turn off. Pressing the power button for multiple seconds will confirm the shutdown even without touching the pop-up on the screen.

If a cuvette is inside the instrument when shutting down, the cuvette is raised and the instrument will ask you to remove it.

## Use of the touchscreen

The touchscreen uses capacitive touch technology and supports gestures such as sliding, where applicable in the user interface. The screen can be used with bare fingers, with disposable gloves, or a touch pen. The touchscreen does not require much pressure when interacting with it. There is multisensorial feedback to touching a button: The button will indicate a touch both visually, by changing appearance, and with a sound. A command is registered when the finger is released from the virtual button touched. If the pressure is released outside the initial button area, no command is registered.


## First-time setup

When starting the QuikRead go Plus Instrument for the first time, the instrument will require the user to set the Admin password. Next, the instrument will prompt guided setup that the user can skip to use factory default settings (**Picture 4**). After completing the first-time setup or skipping it the screen will enter the home screen (**Picture 5**).

To change the settings later, go to the *Menu* ☰ → *Settings*.

**Note!** To find more information about a setting tap the *Info* ⓘ icon next to the setting on QuikRead go Plus Instrument.

**Note!** Settings sent from POCT01-A2 middleware will override settings made on the instrument unless otherwise noted.


**Picture 4**  
Set Admin password

**Note!** For connections use only the strongest available security policy

### First-time setup:

1. Language
  - Select the system language.
2. Display, sound
  - Adjust the screen brightness, alert volume, and keypad volume using the sliders.
3. Date and time
  - Select the date format and set the current date.
  - Select the clock format and set the current time.
4. Power saving
  - Using the slider select the period of inactivity after which the instrument enters sleep mode.
  - Select sleep mode functionality.

**Picture 5**  
First-time setup

5. Measurement
  - You can choose to enable serial measurement.
    - » When enabled the instrument will automatically start a new patient sample measurement after accepting or rejecting a previous result.
  - You can choose to enable emergency measurement.
    - » When enabled a user can perform a test using the emergency measurement without the need to log in. If enabled, set the number of allowed emergency measurements.
6. Test parameters
  - The view shows the currently available tests and currently selected unit selections. Select the test you wish to adjust the test parameters for.
    - » For detailed description of test parameter settings see Section "Test parameters".
7. Printing
  - You can choose to enable printing. If printing is enabled you can choose to enable automatic printing.
8. Operator settings
  - You can choose to enable operator ID. When enabled, the instrument will prompt the user to enter the operator ID when performing measurements.
    - » If enabled, you can choose to enable that the instrument proposes the previously entered operator ID each time.
  - You can choose to set user log in mandatory. If enabled the user must log in using their username and password.
    - » If enabled, using the slider you can select the period of inactivity after which the user will be logged out automatically. You can also select to automatically log the user out after each measurement.
    - » If enabled, you can choose to allow remote operator list.
      - If enabled locally saved log in accounts will be disabled and LIS middleware accounts will be active. Requires an active POCT01-A2 connection to function properly. See Section "Connectivity to setup a LIS connection".

9. General user settings
  - Select whether general users can view result list.
  - Select whether general users can view QC list.
    - » **Note!** When logged in as an admin the result list and QC list are always available for viewing.
10. Patient ID settings and additional ID settings
  - You can choose to enable patient ID. When enabled, the instrument will prompt the user to enter the patient ID when performing measurements.
  - You can choose to enable additional ID. When the additional ID is enabled, the instrument will prompt the user to enter the additional ID when performing measurements.
    - » Read more about Manage additional IDs from Section "Manage additional IDs".
11. Result storage
  - Select for how long the instrument keeps patient results in storage. After the selected time the instrument will automatically delete older results from storage. If left empty instrument will not delete results based on time.
    - » Select years, days, or hours.
  - Select the maximum number of patient results to be stored. When the limit is reached the instrument will automatically delete the oldest result from storage.


## User interface

QuikRead go Plus Instrument is used by a graphical user interface designed to guide its use. In this section the main principles of the user interface are explained. The screen displays virtual buttons and information about the operation and status of the instrument. The screen allows for adjusting sliders by pressing on one and moving it by sliding on the screen. Similarly for certain views with more content than can fit the screen at a time you can slide the view up and down by pressing on the screen outside an active selection region and moving it.

Each view shows the status area and its symbols, the menu bar, the main view, and at the bottom shortcuts for language selection and general settings for display brightness and volume setting.

All users can change the language, display, and sound settings. These changes when made by a general user are reverted when the instrument is turned off or the user logs out.

Touch the icons on the screen to interact with the instrument.

**Picture 6** Home screen

1. Status area
2. Menu bar
3. View
4. Language selection
5. General settings

## QuikRead go Plus Instrument user interface symbols


	Measure sample		List actions
	Quality Control		Home screen
	Menu		Notification
	Print		Language selection
	Comment		General settings
	Reject /Delete		More information

## Status area icons

	Battery level indicator		Power cord connected (battery recharging if connected)
	Battery level low warning (symbol red)		WiFi status
	USB storage connected		Connected to LIS
	Keyboard connected		LIS connection interrupted
	Barcode reader connected		LIS offline results waiting to be sent
	Printer connected		

## Menu

All the features of the instrument can be accessed from the main menu . The menu has the *Measure sample* button  in the top bar. Press this icon for a shortcut to sample measurement. The top bar also has a button to close the menu and return to the previous view and the home button to return to the home screen.


Picture 7 Menu

## Notifications

The instrument notifications can be viewed in the *Notification list*  in the menu bar. The icon change when notifications are available. Selecting a notification from the list will give more information and guide troubleshooting. Notifications are categorised into three classes

- Non-time critical notifications are blue i.e., measurement ready.
- Important errors are yellow i.e., LIS connection error.
- System errors are red i.e., self test failed.


For more information on errors and troubleshooting see Section 10 “Troubleshooting”.

## Language selection

Interface language can be temporarily changed by all users from the *Language selection*  visible in the lower right corner of the screen. The language selection affects the menu and instruction language as well as the keypad layout. This selection is reverted to system language when the users logs out or the instrument is turned off.

Admin user can change the system language from *Menu*  → *Settings* → *Language, display, sound*.


## General settings

General settings can be temporarily changed by all users from the *General settings*  shortcut visible in the lower right corner of the screen. The user can change display brightness, alert volume, and keypad volume. These selections are reverted to system default settings when the user logs out or the instrument is turned off.

- Display brightness can be adjusted by moving the slider.
- Alert volume can be adjusted by moving the slider. Alerts can be muted by setting the slider to the far left.
- Keypad volume can be adjusted by moving the slider. Keypad can be muted by setting the slider to the far left.

Admin user can change these general settings from *Menu*  → *Settings* → *Language, display, sound*.

## 3. INSTRUMENT OPERATION

The use of QuikRead go Plus Instrument can be divided into three main operations:

- Performing a measurement
- Viewing results
- Changing instrument settings

### Performing a measurement

The instrument can be used in two measurement modes: patient sample measurement and quality control sample measurement. Both modes can be modified through the settings.

**Note!** Only QuikRead go reagent kits can be used to perform an assay. Read the instructions for use of the corresponding QuikRead go reagent kit before use. The instructions provide more detailed information on carrying out tests and handling samples.


### Patient sample measurement

When patient sample is to be analysed select *Measure sample*  on the screen. The instrument will guide the user to carry out a QuikRead go test on the patient sample that the instrument will analyse and display the result on the screen once completed.

To perform a measurement:

1. Choose *Measure sample*  on the screen and follow the instructions on the screen.
2. Insert a cuvette into the measurement well with the barcode of the cuvette facing the front (see **Picture 8**)
3. The lid closes and the instrument starts the measurement. The progress of the analysis is shown on the screen.
4. The measurement can be aborted by selecting cancel. This will raise the cuvette and return the view to the home screen.
5. During a measurement the user can view previous patient results and QC results and access the about view through the menu. A progress icon will be shown on the screen when in other views and the instrument will notify the user when the measurement is ready.

6. After the analysis is complete, the result is displayed, and the cuvette is raised.
7. Remove the cuvette, *Accept*  or *Reject*  the result and the screen returns to the home screen. You can view previous results in the patient result list from the *Menu* .

**Picture 8**




Insert a cuvette into the measurement well with the barcode facing the front.


- The admin can enable multiple functionalities to the patient sample measurement (see the Section "Measurement").
- Operator ID, Patient ID, and Additional ID can be enabled and the instrument will require the user to enter these selected IDs to complete the measurement. Read more about Additional ID in Section "Additional ID settings".
- The IDs can be entered with the touchscreen or with the help of a barcode reader or an external keyboard. It is possible to enable the instrument to propose the last used Operator ID in the prompt when log in is disabled.
- Log in can also be enabled to require the user to log onto the instrument using their username and pass-

word. Log in automatically enters the Operator ID for each measurement.



- Emergency measurement can be enabled. When enabled a user can perform a test using the emergency measurement without the need to log in. If enabled, set the number of allowed emergency measurements. The number of acceptable emergency measurements can be set.
- Serial measurement mode can be enabled. When enabled the instrument will automatically start a new patient sample measurement after the user accepts or rejects the previous result.
- Results printing and sending results to a LIS system can be enabled through the settings.

## Patient sample result view

The patient sample result view is available after a completed analysis. The analysis result and the test information are displayed. The user can accept the result and it will be stored in patient result list using the  icon. The result can be printed by selecting  if printing is enabled and a printer is connected. A comment can be attached to the result. To write a comment select . The user can decide to reject the result by selecting . When rejecting the user must leave a comment. The result will be stored in the patient result list as rejected and the detailed information can still be viewed.


## Detailed sample result view

Detailed sample result view is available through *Menu*  → *Patient results* by selecting a result line. The detailed sample result view displays additionally the instrument serial ID, comments, and additional patient information. For CRP test results sample type will be displayed. It is possible to add comments, by selecting , and enter the additional patient information.

## Quality control sample measurement

The instrument has a separate measurement mode for quality control samples that are measured like patient samples through the *Quality Control*  on the home screen. The results are stored separately from patient results and can be viewed from QC results in the *Menu* .

To perform quality control choose *Quality Control*  on the home screen and follow the instructions on the screen.

The admin can enable multiple functionalities to the quality control sample measurement (see Section "QC settings").

Control ID can be enabled and the instrument will require the user to enter it to perform the quality control. QC checking can be enabled to indicate whether a QC result is within set range for the QC sample.

When enabled QC lock can be enabled to restrict use of the instrument in case the last QC measurement failed or the QC schedule has expired. Once expired the user can't measure patient samples before performing a quality control that passes.


## QC result view

The QC result view is available after a completed QC measurement. The result and the test information is displayed. The user can accept the QC result and it will be stored in QC result list using the  icon. The result can be printed by selecting  if printing is enabled and a printer is connected. A comment can be attached to the result. To write a comment select .

## Detailed QC result view

Detailed QC result view is available through *Menu*  → *QC results* by selecting a result line. The detailed QC result view displays additionally the QC concentration, additional QC information, and comments. It is possible to add comments by selecting .


## Results

Results can be viewed through the *Menu* ☰. Results are stored in separate lists for *Patient results*, *QC results*, and *LIS offline results*. Results view can be scrolled by swiping. Selecting a result displays its detailed information. The result lists can be sorted by touching a column header.

The patient list view can be searched by *Patient ID* and sorted by *Date*, *Test*, and *Operator ID*. The QC result list view can be filtered by *Test*, *Control ID*, *QC date*, and *Operator ID*. To filter a list access the *List actions* ⓘ from the menu bar.

Individual results and result lists can be printed. The lists can be printed with or without filters and transferred to USB storage. Result list needs to be filtered before USB transfer is enabled. To print an individual result select the result line and select print on the results view. To print a list access the *List actions* ⓘ from the menu bar.

The screenshot shows a mobile application interface for 'Patient results'. At the top, there is a status bar with the date '2023-07-06 08:05' and a home indicator. Below the status bar is a navigation bar with a menu icon (☰), a home icon, a notification icon, and a profile icon. The main title is 'Patient results' with a filter icon (⋮) on the right. Below the title is a table with columns: Patient ID, Time, Test, and Result. The table contains several rows of data, including patient IDs like C10190-KKXX, 2F0212-YYYY, C10656-AAAA, and test names like CRP, Strept A, FOB, and HbA1c. The results are displayed in a clean, modern font with alternating row colors.

Patient ID	Time	Test	Result
C10190-KKXX	17:00 05-19-2024	CRP	< 1.0 mg/l
2F0212-YYYY	09:05 04-19-2024	Strept A	1 Pos +
C10656-AAAA	19:12 04-18-2024	FOB	Neg -
Jimmy Johnny Doe	18:32 04-18-2024	C RP Hb	12 mg/ 122 g/l
--	18:05 04-18-2024	HbA1c	Rejected
--	17:59 04-18-2024	C RP	2.4 mg/l

## Deleting results

An admin can delete individual results from the instrument memory. Select the result from the list and select delete result from the bottom of the view. The instrument will require an additional confirmation.

## LIS offline list

*Menu* ☰ → *LIS offline list*

Results that have been generated by an instrument with LIS connection enabled while offline are saved in *LIS offline results*. Results are saved in case the connection is interrupted e.g., due to a momentary network problem. The results can be sent to LIS after the connection is restored. Once the results are transferred to LIS they are removed from the LIS offline list.

To send offline results to LIS, select the results by tapping the check box on the result line and selecting *Send*. To delete results, select the results to be deleted and select delete result. All results can be selected at once by tapping the check box on the title row of the results table.


When automatic offline result transfer is enabled in general LIS settings and the LIS connection is restored a pop-up is displayed for sending the offline results to LIS. The instrument checks for the LIS connection during start-up and after every measurement.

## 4. SETTINGS

Settings can be accessed from the *Menu* ☰ by selecting *Settings*. In the settings, the admin user can set permanent settings for the following. General users can change some settings and view those that require admin authorisation.

- Language, display, sound
- Date and time
- Power saving
- Measurement
- Accessory
- User management
- Connectivity

**Note!** The system uses an *Auto-save feature* **Saved** ✓. Any changes in the changes will be saved automatically. When auto-save is active **Saved** ✓ will be shown in the menu bar.

## Language, display, sound

Menu  → Settings → Language, display, sound

- System language can be changed by selecting a language from the selection. The language selection affects the menu and instruction language as well as the keypad layout.
- Display brightness can be adjusted by moving the slider.
- Alert volume can be adjusted by moving the slider. You can mute alerts by setting the slider to the far left.
- Keypad volume can be adjusted by moving the slider. You can mute the keypad by setting the slider to the far left.

## Date and time

Menu  → Settings → Date and time

- The date format can be changed from the selection. The date can be changed by setting a new date.
- The clock format can be changed between 24 h and 12 h format. The current time can be changed by setting a new time.

**Note!** The instrument checks the expiry date of QuikRead go reagents before a measurement. Setting the wrong date might restrict the use of reagents due to failed expiry date check.

## Power saving

Menu  → Settings → Power saving

- The sleep mode delay can be changed by moving the slider. After this time the instrument enters sleep mode.
- The sleep mode functionality can be changed from the selection.
  - » Full stand-by: In sleep mode the instrument closes the lid and dims the screen before turning it off. Screen can be awakened by touching the screen or pressing the power button quickly.
  - » Close lid only: In sleep mode the instrument closes the lid and the screen is dimmed.

## Measurement

Menu  → Settings → Measurement


- Change test parameters by selecting a test. Each analyte has specific settings for that test type.
- Serial measurement mode can be enabled or disabled. When enabled the instrument will automatically start a new patient sample measurement after accepting or rejecting a previous result.
- Emergency measurement mode can be enabled or disabled. When enabled a user can perform a test using the emergency measurement without the need to log in. If enabled, set the number of allowed emergency measurements.

## Test parameters

Menu  → Settings → Measurement → Test parameters

- Change the unit selection of in which unit(s) the results are shown for the test. For some tests it is possible to switch between quantitative and qualitative results.
- If multiple units are selected for displaying results and the unit in which the result is sent to LIS can be selected.


- Indicator limits can be set for each test. Limits can be set for low, decreased, increased, and/or high results. The test result view will display the indicator limits set and the result. If the result is below the limit for decreased a single down arrow will be displayed next to the result. Similarly, if the result is below the limit for low two down arrows will be displayed and if the result is above the limit for increased a single up arrow will be displayed and if the result is above the limit for high two up arrows will be displayed.

- For qualitative results the indicator can be set for "Pos +" or "Neg -". When the sample result is the according to the selection an exclamation mark will be displayed next to the result. When enabled the results view will display the text "Result indicator on" for all results.

Sample type flagging can be enabled for some tests. When enabled the instrument will give an error when it detects the wrong type of sample during measuring.

## Accessory

Menu  → Settings → Accessory

Printing can be enabled or disabled. When printing is enabled and a compatible printer is connected automatic printing can be enabled. With automatic printing the results of each sample will be automatically sent to the printer.

<b>CRP</b> <b>Hb</b>	<b>21 mg/l</b> <b>100 g/l</b>	<b>QuikRead go Plus</b> 38/A22064512288
Patient ID Operator ID Time Sample type	Joe Smith Nurse 1 2021-01-29 10:18 Whole blood	 JE45  JF08  2025-01-29

Picture 9 QuikRead go Plus result when printed

## User management

Menu  → Settings → User management

### Operator settings

- Operator ID can be enabled or disabled. When enabled, the instrument will prompt the user to enter an operator ID when performing measurements.
  - When enabled, you can choose to enable that the previously entered operator ID will be proposed each time.
  - Log in can be enabled or disabled. If enabled the user must log in using their password.
  - When enabled, you can choose after what period of inactivity the user will be logged out automatically using the slider. You can also select to automatically log the user out after each measurement.
- Note!** It is strongly recommended to enable password login to ensure information security, particularly in environments with low physical security measures, such as shared or unlocked working spaces or when visiting patients.
- When log in is enabled, you can choose to enable remote operator list. This option can be enabled when the instrument has POCT01-A2 connection enabled. When enabled the operator list is managed by the

LIS middleware and cannot be edited through the instrument. Enabling log in will disable the propose previous selection.

- Patient result list viewing can be enabled or disabled for general users.
- QC result list viewing can be enabled or disabled for general users.

**Note!** When logged in as an admin the patient result list and QC list are always available for viewing.

### Manage operators

Menu  → Settings → User management → Manage operators

- An admin can add new operators, view complete operator list and roles, and edit accounts.
- When adding a new operator ID:
  - Select *Operator role* (General/Admin)
    - » Enter Operator ID
    - » Enter Password
    - » Enter operator name
  - When editing accounts, the user password can be changed.

See user role rights table on page 29.


### Patient settings

Patient ID can be enabled or disabled. When enabled, the instrument will prompt the user to enter the patient ID when performing measurements.


### Additional ID settings

Additional ID can be used as an additional text field or to select from a predefined list of additional ID items for each patient sample measurement. When the additional ID is enabled, the instrument will require the user to enter an additional ID when performing measurements. To enter the additional ID list items or customise the list name select *Manage additional IDs*.

### Manage additional IDs

Menu  → Settings → User management → Manage additional IDs

Enter additional ID list items by selecting *Add list item*. Adding a list item will activate a list button next to the text field when performing a measurement (**Picture 10**). Customise the list name by selecting *Edit*  next to the list name. This will replace the *Additional ID* title when performing a measurement (**Picture 10**).

Picture 10 Additional ID

When additional ID list items are predefined a list button is activated next to the text field (right). When the additional ID list name has been edited the name will be shown above the text field (right).

## Connectivity

Menu  → Settings → Connectivity

The instrument can be connected to a Laboratory or Hospital Information System (LIS/HIS) with a LAN or WLAN connection. These connections can be used to automatically transfer data between the instrument and the system for remote management and result viewing. Contact your supplier for more information about connectivity.

**Note!** Contact your system administration for your connectivity settings.

### LIS

Menu  → Settings → Connectivity → LIS

- LIS connection can be enabled or disabled. When enabled automatic result transfer can be selected. If connection is interrupted results that have been generated but not sent to LIS will be visible in the LIS offline results view (see Section "LIS offline list")
- Connection type can be selected. Available connection types are:
  - » POCT1-A2
  - » LIS01-A2 (TCP/IP)
  - » LIS01-A2 (RS232)
- LIS server address and TCP port can be changed.
- Conversational interval can be changed using the slider.
- TLS encryption can be enabled. The TLS certificate can be uploaded from a USB storage.
- Patient identifier validation can be enabled. When enabled the instrument will check the patient ID from the POCT01-A2 middleware. If the patient ID entered is invalid the instrument will give an error message.
- Patient data shown can be enabled. When enabled the instrument will fetch patient data from POCT01-A2 middleware based on the entered patient ID.

- Offline result transfer can be toggled between automatic and manual mode. When automatic is selected LIS offline results will be automatically transferred when connection is re-established. When manual is selected the user will need to choose the results to be transferred to LIS.

### TCP/IP

Menu  → Settings → Connectivity → TCP/IP

- IP address type can be toggled between Static IP and Dynamic IP.
- IP address can be changed.
- Subnet mask can be changed.
- Default gateway can be changed.
- DNS server can be changed.

### WLAN

Menu  → Settings → Connectivity → WLAN

WLAN connection can be enabled or disabled. When enabled connection can be set by scanning the available WLAN networks by choosing Scan or through manual setup. Scan will show available WLAN networks list. Select the desired network from the list. If the network is password protected a password pop-up will be shown.

## 5. QC MANAGEMENT

QC management can be accessed from the *Menu*  by selecting *QC management*. QC management settings control the Quality Control measurement functionality intended for monitoring the performance of the instrument using QuikRead go reagent kits and controls.

### QC settings

*Menu*  → *QC management* → *QC settings*

Control ID can be enabled and the instrument will require the user to enter it to perform the quality control. If enabled, QC checking can be enabled to indicate whether a QC result is within set range for the QC sample. When QC checking is enabled and QC result is within the accepted range *QC passed* will be displayed with the result. If the QC result is out of range *QC failed*. *Result out of range* will be displayed with the result.

When enabled QC lock can be enabled to restrict use of the instrument in case the last QC measurement failed or the QC schedule has expired. Once expired the user can't measure patient samples before performing a quality control that passes. Set the schedule and select daily, after selected number of tests, or selected days. The instrument will require an input for the latter two. If schedule is set for daily the QC lock will become active on the following day if QC measurement wasn't done. If schedule is set for selected days the QC lock will become active on the following day if QC measurement wasn't done on time.

### Add new QC

*Menu*  → *QC management* → *Add new QC*


Add new QC is used to input information of a new quality control (QC) when Control ID setting is enabled in QC settings. The added QC will be available for selection when performing QC measurements. The instrument will use the control information when determining whether a QC is passed or failed.

QC name can be entered to set a name visible in QC list selection. QC unit displayed with the results can be changed. QC range is mandatory and is used to determine whether a QC measurement using the control ID is passed or failed. QC concentration is mandatory and is the target value of the control measurement. QC lot expiry is mandatory and is used to check if expired reagents are used for QC.

QC vial in use expiry date can be set. QuikRead go reagent kit control vials have an expiry date and a separate in-use expiry period once opened. To assist in managing the QC vial in use expiry periods the instrument offers the option to enter the last use date for an opened QC vial. When performing QC and selecting a control ID if the open vial expiry date has passed the instrument displays a pop-up. The pop-up suggests opening a new vial and entering a new open vial in use expiry date.

Additional QC information can be entered for each Control ID.

**Note!** Refer to the label and the instructions for use of the control when entering information for a new control.

### QC list

*Menu*  → *QC management* → *QC list*

Control IDs' information can be viewed and changed from QC list. The list view shows name, expiry date, test, and range for each control.

The list can be sorted using the column headers.

## Add new QC vial in use

Menu  → QC management → Add new QC vial in use


A new QC vial in use expiry date can be entered for an existing control ID. Select the control ID you want to add a new QC vial in use expiry date from the list. A pop-up will appear guiding with the setup.

The date can be added as a date or as number of days from the date of input. A previously entered open vial expiry date can be deleted.

**Note!** Check the in-use expiry period from the instructions for use of the control.

Exp.ry date		
2025	June	30
2024	July	31
2025	August	1
2026	September	2
2027	October	3
2028	November	4

The QC has an open vial with expiry date  
2024-05-15

Do you want to replace this vial with a new one?

No Yes

## 6. MAINTENANCE MENU

Maintenance menu can be accessed from the Menu  by selecting *Maintenance*.

### Software update

Menu  → Maintenance → Software update

The instrument software version can be updated with the help of a USB storage. Contact your local supplier for more details on software updates.

### Logs

Menu  → Maintenance → Logs


The instrument records measurement, error, and connectivity and audit logs in its memory. Read the following sections for further information.

### Error log

Menu  → Maintenance → Logs → Error log

The instrument errors are stored in its memory. Previous errors can be viewed in the error log. The view can be sorted by touching a column header. Error codes can be transferred to a USB storage by selecting the transfer to USB in the menu bar.

### Audit logs to USB

Menu  → Maintenance → Logs → Audit logs to USB

Audit logs contain login information, records of sample measurement actions, rejections of results, result deletions, and log deletions. The logs can be used to monitor the use of the instrument.


The instrument audit logs can be transferred to USB storage.

### Measurement logs to USB

Menu  → Maintenance → Logs → Measurement logs to USB

Measurement logs are used for troubleshooting by Aidian for example in case of a customer complaint case. Measurement logs are encrypted and only available for Aidian personnel. They contain instrument performance logs without patient or user information.

### Connectivity logs

Menu  → Maintenance → Logs → Connectivity logs

Connectivity logs can be stored in the memory. Change LIS log collection.


Connectivity logs can be transferred to a USB storage.

### Self test

Menu  → Maintenance → Self test

The instrument performs operations check to ensure proper functioning. To manually perform a self test select Yes. The system will return to the home screen after complete.


### Patient result management

Menu  → Maintenance → Patient result management

The duration of patient result storage can be set or changed. Set the storage time in years, days, or hours. Older patient results will be automatically deleted after the selected time. If left empty instrument will not delete results based on time.


The maximum number of patient results in storage can be changed. When the limit is reached the instrument will automatically delete the oldest result from storage.

### Delete patient result list

Menu  → Maintenance → Delete patient result list

Patient result list can be deleted entirely. This option can be used for data protection purposes.

### Factory reset

Menu  → Maintenance → Factory reset

The instrument can be reset to factory settings.

## 7. ABOUT VIEW

The *About* view shows the serial number of the instrument, current software version, and connection information. In addition, the instrument can be named. Once named, the instrument name is displayed in the lower left corner of the screen in all views.

About can be accessed from the *Menu*  by selecting *About*.

## 8. INSTRUMENT MAINTENANCE

QuikRead go Plus Instrument has been designed to be user-friendly with no need for regular maintenance. For repair and service requirements, contact your local supplier.

### Instrument calibration

The instrument is factory calibrated. The proper functioning of the instrument is checked by the self-check procedure during start-up and with every measurement. In case of malfunction, an error message is displayed. The calibration data defining each test is encoded on the cuvette labels. This information is read automatically by the instrument with every measurement.


### Cleaning the instrument

Periodically clean the exterior of the instrument using a lint-free cloth dampened with water. Pay special attention to cleaning the display. Take care that no liquid enters the instrument through the edges of the display, the measurement well opening, or the connectors. If necessary, a mild detergent may be used. Do not use organic solvents or corrosive substances. Spillage of potentially infectious material should be wiped off immediately with an absorbent paper tissue and the contaminated areas swabbed with a standard disinfectant or 70% ethyl alcohol. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

Acceptable disinfectants:

- 70% ethyl alcohol
- 70% isopropanol
- 0.5% sodium hypochlorite
- 2% glutaraldehyde

### Software update

The instrument software can be updated. Access the software update through *Menu*  → *Maintenance* → *Software update* and follow the instructions on screen. Ask your local supplier for more information.

### Changing the clock battery

The instrument has a battery for the internal clock to keep time when powered off. If the clock battery charge is running low a warning is displayed. The clock battery can be replaced with a type CR 2032 3V battery (**Picture 3**).

## 9. INSTRUMENT SPECIFICATION

### Declaration of conformity

QuikRead go Plus Instrument complies with the regulation (EU) 2017/746 on *in vitro* diagnostic medical devices, the directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment together with the delegated directive (EU) 2015/863 amending Annex II to Directive 2011/65/EU and the directive 2012/19/EU on waste electrical and electronic equipment (WEEE). QuikRead go Plus Instrument complies with the electromagnetic emission and immunity requirements described in the standard IEC 61326-2-6:2012. The instrument fulfills the FCC class A requirements. QuikRead go Plus Instrument complies with the regulation (EU) 1907/2006 on the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH). The instrument also complies with the radio equipment directive 2014/53/EU.

For further information, please contact your local supplier.

## Technical specification

The instrument has a pre-programmed microprocessor that controls the assay steps and data processing. Test identification, timing, and calibration curve or cut-off value data are contained on a barcode on each cuvette. Once activated by the cuvette label, the microprocessor controls and guides all assay steps and converts the absorbance values of the samples into concentration units or cut-off values.

### Photometer

QuikRead go Plus Instrument photometer consists of a measurement well, three LEDs, and light detectors. The photometer has been designed and calibrated both for photometric and turbidimetric measurements.

### Touchscreen display

The user interface is based on an easy-to-use capacitive touchscreen. It provides the user with messages, prompts, and animations for performing each assay step, and displays test results and error messages.

### Dimensions & power requirements

- Weight: 2.0 kg without power supply
- Size: 26.9 x 16.6 x 13.9 cm
- Power requirements
  - Voltage: 100-240 V AC
  - Frequency: 50-60 Hz
  - Power consumption:
    - Idle 4 W
    - Stand-by 1 W
    - Max 35 W (without external accessories)
  - Input 18V --- 3.3A

### Instrument identifier

Each QuikRead go Plus Instrument has a unique serial number which can be found on the label at the bottom of the instrument and in the *About* view accessible from the *Menu* .

### Memory

QuikRead go Plus Instrument has an internal memory for results history and logs. The instrument can store 6,000 patient results and 6,000 QC results.

## Power supply

The instrument is powered by a power supply supplied with the instrument. In addition to the power supply, the instrument can use a battery as a power source. An internal switch inside the cable connector will switch automatically from battery-powered use to mains-powered use. For instructions on installing a battery, see Section “Inserting the rechargeable battery”.

## Radio module

Wi-Fi Dual Band 2.4 GHz/5 GHz and Bluetooth Module  
Panasonic model: PAN9028.

## LIS connection

The connection can be made using:

- An RJ-45 connector and a supported 10BASE-T/100BASE-TX/1000BASE-TX Ethernet connection. Unshielded twisted pair cable should be used.
- A WLAN connection.
- Power over Ethernet (PoE) is not supported.

Contact your supplier for more details.

## USB connection

The instrument has three USB type A connectors and two USB type C connectors. These connectors can be used for printers, barcode readers, keyboards, and USB storages.

## Service

QuikRead go Plus Instrument is designed to be free of regular maintenance with built-in self-check operations. In case of instrument malfunction, please contact your local supplier. Before shipping the instrument to service, delete all patient results from the maintenance view by selecting the delete patient result list and clean the exterior of the instrument. See the Section “Cleaning the instrument”.

## Warranty

The manufacturer’s warranty for QuikRead go Plus Instrument covers defects in materials or manufacturing for a period of two years from the date of purchase. For the warranty to be valid, the warranty seal (**Picture 3**) needs to be intact. The manufacturer agrees to repair

or replace the instrument if it becomes inoperative due to the failure of any internal part of the instrument. The warranty does not cover damage caused by inappropriate use. The manufacturer is under no obligation to modify or update the instrument once it has been manufactured unless a manufacturing defect is identified. In the event of an instrument malfunction, please contact your local supplier.

## Disposal

QuikRead go Plus Instrument is a low-voltage electronic device. A used QuikRead go Plus Instrument must be treated as potentially biohazardous waste. The instrument should be disposed of as electrical and electronic equipment waste (WEEE 2012/19/EU) if local and national law do not require the instrument to be collected and disposed as potentially infective clinical waste. You can find more information about instrument disposing on our website [aidian.eu](http://aidian.eu).

The packaging materials are recyclable materials. The battery unit and the internal clock battery should be disposed of according to national and local battery collection rules based on Directive 2006/66/EC or Batteries Regulation 2023/1542/EU.

## Revision history

Revision history can be found from [aidian.eu](http://aidian.eu).

Frequency bands and their maximum radio-frequency powers:

Supported technologies	Operation mode	Band	Frequency range	Maximum conducted average output power	Antenna gain (integrated chip antenna)	Maximum radiated average output power
			[MHz]	[dBm]	[dBi]	[dBm]
Bluetooth	Basic Rate (BR)	2.4 GHz	2 400... 2 483.5	4	2.1	<10
Bluetooth	Enhanced Data Rate (EDR)	2.4 GHz	2 400... 2 483.5	4	2.1	<10
Bluetooth	Low energy (LE)	2.4 GHz	2 400... 2 483.5	4	2.1	<10
WLAN	IEEE 802.11b/g/n	2.4 GHz	2 400... 2 483.5	16	2.1	<19
WLAN	IEEE 802.11a/n/ac	5 GHz	5 150... 5 250	16	1.5	<18
WLAN	IEEE 802.11a/n/ac	5 GHz	5 250... 5 350	16	1.5	<18
WLAN	IEEE 802.11a/n/ac	5 GHz	5 470... 5 725	16	1.5	<18
WLAN	IEEE 802.11a/n/ac	5 GHz	5 725... 5 825	16	1.5	<18

## 10. TROUBLESHOOTING

QuikRead go Plus Instrument displays error messages and guides the user in case it detects errors. Follow the instruction displayed and see the troubleshooting table in this instruction for use and QuikRead go kit instructions for use.

**Note!** For more troubleshooting support visit [aidian.eu](http://aidian.eu) or contact your local supplier.

Error message / Troubleshooting	Possible cause	Corrective action
Error code with the message "Please restart QuikRead go Plus" is displayed.	Temporary malfunction of the instrument.	Restart the instrument. If this error message is displayed frequently, contact customer service.
Error code with the message "Please contact customer service" is displayed.	Permanent malfunction of the instrument.	Contact customer service.
"Battery level is low. Please connect to the mains cable in order to continue operation" an error message is displayed.	The charge in the battery is low.	Connect the power unit into QuikRead go Plus Instrument power connector.
"Cuvette's position is not correct. Remove cuvette." error message is displayed.	Remnants of cuvette sealing foil is left on the cuvette collar.	Remove cuvette when the instrument has lifted it. Make sure that all remnants are removed when doing the next measurement.
	The instrument has a mechanical malfunction.	Check the item above. If it does not apply, restart the instrument. If the problem persists, contact customer service.

Error message / Troubleshooting	Possible cause	Corrective action
"Measurement prohibited."	A reagent cap is missing or a cuvette is used.	Check that the cuvette has a reagent cap and the inner coloured part of the cap is not pressed down.
	Reading lot data from the barcode failed.	Try again. If the problem persists, cancel the test.
	Kit lot expired.	Discard the expired kit lot. Use a new one.
	The cuvette temperature is too low.	Let the cuvette warm up to room temperature. Test the same cuvette again.
	The cuvette temperature too high.	Let the cuvette cool down to room temperature. Test the same cuvette again.
"Test cancelled."	The blank too high.	Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed.
	Unstable blank.	
	Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. Make sure that the cap is properly closed.
	Instrument failure.	Perform a new test. If this message appears often, contact customer service.
QuikRead go Plus Instrument does not start.	The power unit is not connected.	Connect the power unit and retry.
	Battery has emptied.	Connect the power cable and retry.
	The instrument has an electronic malfunction.	Contact customer service.

Error message / Troubleshooting	Possible cause	Corrective action
The touch screen is not functioning properly.	Touch panel is not responding at all.	Contact customer service.
Instrument alarm sounds cannot be heard.	Volume is set at a low level.	Set volume according to the procedure described in the Section "Language, display, sound".
	Instrument sound system has a malfunction.	Restart QuikRead go Plus Instrument. If the problem persists, please contact customer service.
The printer is not printing.	The printer is off or the printer cable is not connected or the printer has malfunctioned or the settings are not correct.	Make sure that the printer is connected and its power is on. Check the settings. If the problem continues start-up the instrument and printer and retries printing from the Results menu. If the problem continues, contact customer service.
Barcode reader is not functioning.	The barcode reader is not connected or the barcode reader has a malfunction or the settings are not correct.	Make sure that the barcode reader is connected. Check the settings. If the problem continues start-up the instrument and retries the barcode reading. If the problem continues, contact customer service.
The battery must be recharged frequently.	The battery storage capacity decreases during its lifetime.	Replace the old battery with a new one according to the procedure described in Section "Inserting the rechargeable battery".
A clock battery warning is displayed.	The internal clock battery is empty.	Replace the clock battery according to the procedure described in Section "Changing the clock battery".

## User role rights

Action	General user	Admin
Patient measurement	X	X
View patient result list	- / X*	X
View LIS offline result list	- / X*	X
QC measurement	X	X
View QC results	X	X
Add new Control ID	X	X
Clear all patient information	-	X
Transfer results to USB	- / X*	X
Change password	X**	X
Set Admin password	-	X
Settings	X***	X
Power saving	-	X
Temporary change: Power saving	X	-
Language, screen brightness, volume	-	X
Temporary change: language, screen brightness, volume	X	X
Maintenance	-	X
Maintenance: error log	X***	X
About	X***	X

\*If enabled by Admin    \*\* Only local user. Restricted when remote operator list is active.    \*\*\* Viewing only

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QuikRead and QuikRead go are trademarks or registered trademarks of Aidian Oy.
Aidian Oy