








154470-4

- English
- Deutsch
- Français
- Nederlands

# QuikRead go<sup>®</sup> HbA1c

## Assay Procedure • Testverfahren

Measure		RESULT
HbA1c:		42 mmol/mol
Patient ID:	XXXXXXXXX	
Measurement time:		2019-02-21 12:19
Test:	HbA1c	Result info
Choose <b>Result info</b> to view result information. Remove cuvette to perform a new measurement.		
Exit	Print	New measurement

## 1 Intended use

QuikRead go<sup>®</sup> HbA1c is an *in vitro* diagnostic test for quantitative measurement of glycosylated hemoglobin (HbA1c) from a finger prick capillary blood or anticoagulated (EDTA or heparin) venous whole blood sample. The test is carried out using the automated QuikRead go<sup>®</sup> instrument. The assay is intended to be used by trained healthcare professionals in clinical laboratories and near-patient testing (NPT) settings. The test can be used in monitoring the long-term blood glucose control in individuals with diabetes mellitus, as an aid in diagnosis of diabetes and in identifying patients at risk of developing diabetes mellitus.

**For *in vitro* diagnostic use only.**

## 2 Summary and explanation of the test

Hemoglobin (Hb) is an iron-containing protein located in red blood cells. The most prevalent variant of Hb in adults is HbA covering over 95% of total Hb. Glycation of HbA is a two-step process where glucose binds chemically to the N-terminal end of hemoglobin beta chain<sup>1</sup>. This process happens in every red blood cell. The relative amount of HbA1c increases with glucose levels in blood during the life span of red blood cells (2–3 months)<sup>2,3</sup>. As a result, HbA1c concentration rises in individuals who have chronic elevated blood sugar levels, meaning diabetic patients. HbA1c is a measure of long-term glycemic status for diabetic patients whereas blood glucose levels represent the daily variations. Quantitative measurement of HbA1c concentration is an established method used for estimating the status of diabetes treatment<sup>4–6</sup>. HbA1c can also be used as an aid in diagnosis of diabetes mellitus<sup>7</sup>.

## 3 Principles of the procedure

QuikRead go HbA1c is an immunoturbidimetric test based on agglutination reaction. The blood sample is added into the cuvette containing a hemolysing solution with latex particles. The released hemoglobin is bound to the latex particles. The addition of a second reagent containing anti-HbA1c antibodies induces an agglutination reaction together with the latex-bound HbA1c. Hemoglobin and HbA1c bind to the particles in the same ratio that is present in the blood.

The change in turbidity is measured photometrically. The turbidity is thus proportional to the amount of HbA1c molecules bound to the latex and, moreover, to the proportion of HbA1c in total hemoglobin. The calibration data of the test is in the barcode label of the cuvette, which QuikRead go instrument reads automatically before the start of the measurement. The HbA1c measurement range of the assay is 20–140 mmol/mol.

## 4 Reagents

### Content of the kit

Component name and origin	Symbol	QuikRead go <sup>®</sup> HbA1c Cat. no. 151058 25 tests
QuikRead go HbA1c Reagent Caps ORIG MOU	REAG CPS	25
Prefilled cuvettes	CUVET	25 x 0.9 ml
QuikRead go Sample Collector 1µl (sample collector)	SAMPL COL EDTA	25
Instructions for use		

Reagent caps contain methylchloroisothiazolinone and methylisothiazolinone, please see Section 5. "Warnings and precautions".

### Storage

Kit component	Opened at 2...8°C	Opened at 18...25°C	Unopened at 2...8°C	Unopened at 18...25°C
Reagent Caps	6 months	2 months	Until the expiry date of the kit	2 months
Prefilled cuvettes in foil pouch	6 months Single prefilled cuvette: 1 hour	2 months		
Sample collector	Storage at 2...25°C until the expiry date of the QuikRead go Sample Collector 1µl package			

Mark the date of opening the foil pouch on to the cuvette rack.

### Reagent preparation and storage conditions

All reagents are ready to use. Let the reagents reach room temperature. Keep the HbA1c Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

### Reagent deterioration

Do not use cuvettes with visible dirt in the solution. The solution shall be milky and homogenous.

## 5 Warnings and precautions

### Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled.
- Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- Lyophilised reagent inside the caps contains 0.0064% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3). May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing dust (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Reconstituted and liquid reagents contain < 0.0015% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1), which is not considered as a harmful concentration.
- Disposal: see Section 15.

**Warning**

### Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Do not freeze the QuikRead go HbA1c kit.
- The QuikRead go HbA1c kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components for performing a test.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. Always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- QuikRead go HbA1c Reagent Caps are colour-coded light violet to distinguish them from other QuikRead analytes. QuikRead go Sample Collector 1µl has a colorless plastic frame. Do not use broken sample collectors.
- Keep the QuikRead go HbA1c Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps. Do not use the reagent caps if they have not been stored according to the instructions.
- Do not splash liquid into the measuring well of the instrument.
- Do not perform the HbA1c test on an inclined surface.

## 6 Sample collection and preparation

### Sample material, collection and volume

QuikRead go Sample Collector 1µl (Cat. no. 154457) included in the kit (Cat. no. 151058) is required for transferring the sample into a cuvette. The capillary part of the sample collector is covered with EDTA and the appearance might vary depending on the crystallization of the EDTA.

Sample material	Sample volume	Sample collection
Finger prick capillary blood	1 µl	Prick the clean and dry finger with a lancet. Wipe off the first drop and collect 1 µl of blood from the second drop into EDTA covered capillary part of the sample collector.
Anticoagulated whole blood	1 µl	Use a venous blood sample collected in a tube containing Li-heparin, K2-EDTA or K3-EDTA. Mix the whole blood by inverting the tube several times and collect 1 µl using the sample collector.

### Sample dilution

Do not dilute samples.

### Sample storage

Sample material	Short-term storage	Long-term storage
Finger prick capillary blood	Max. 1 minute in the sample collector.	Do not store
Anticoagulated whole blood	2...8°C for 3 days. Do not use hemolyzed samples.	Do not store
Sample in the cuvette	18...25°C for max. 5 minutes in the sample collector inside the cuvette with a closed cap.	Do not store

Allow the samples to reach room temperature (18...25°C) before testing.

## 7 Procedure

### Materials required but not provided

The following materials are required to perform the test but are not provided or they are available also separately. The materials provided are listed in Section 4. "Reagents".

Material	Cat. no.	Additional info
QuikRead go® Instrument with software version 8.1.1 or newer	133893, 149915	–
QuikRead go® Sample Collector 1µl	154457	–
QuikRead go® HbA1c Control Set	154520	For quality control
Fingertip lancets	–	–

### Assay procedure

Before starting an assay, read and follow the QuikRead go instrument and QuikRead go HbA1c test instructions for use. Open the foil pouch of the cuvette rack, and mark the opening date on the cuvette rack.

### Sample collection (see figures 1–6)

QuikRead go Sample Collector 1µl is required for transferring the sample into a prefilled cuvette. The QuikRead go HbA1c kit contains 25 pieces of the sample collectors (QuikRead go Sample Collector 1µl).

1. A prefilled cuvette must reach room temperature (18...25°C) before use. This will take 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed solution on the foil cover has no effect on results. The test must be performed within 1 hour of opening the cuvette.
2. Open the sample collector box slightly to get one single sample collector out of the box. The package filling can be discarded.
3. Hold the sample collector almost horizontally and take the sample by placing the tip just beneath the blood drop. Fill the capillary part of the sample collector completely.
4. Check that the capillary part is completely filled. Ensure that there are no air bubbles in the capillary part. Do not wipe the sample collector. If there is excess of sample on the sample collector, discard the sample collector and take a new sample with a new sample collector.
5. Insert the sample collector into the cuvette within 1 minute. The sample collector does not touch the solution in the cuvette. The sample collector remains in the cuvette during the whole test measurement procedure.
6. Close the cuvette tightly with the HbA1c Reagent Cap. Do not press down the violet inner part of the reagent cap. Perform the test within 5 minutes. Keep the cuvette upright and do not shake.

**Note!** If a venous sample is used, it is recommended to take the sample directly from the tube. Tilt the tube slightly and collect the sample just beneath the surface of the sample with the sample collector. Do not wipe the sample collector.

### Analysing the sample (see figures 7–9)

Begin the test by turning on the instrument. For full instructions, see the QuikRead go instrument instructions for use. The display will guide you through the test.

7. Choose **Measure** on the display of the QuikRead go instrument. Insert user ID and patient ID if needed.
8. Place the cuvette in the measurement well of the instrument. The barcode should be facing you as shown in the Figure 8. The display shows how the measurement is proceeding.
9. The result appears on the display and the cuvette rises up automatically when the measurement is completed.

### Analysing the control

Please see the instructions for use provided with each control set. Analyse control samples as patient samples but choose **Quality Control** on the display of the QuikRead go instrument. The result will be stored as a control measurement.

## 8 Quality control

For quality control use QuikRead go HbA1c Control Set (Cat. no. 154520). Regular use of controls is recommended. Please consult Section 16. "Troubleshooting" and contact Aidian or your local distributor if you observe any changes in test functionality.

QuikRead go HbA1c Control Set features two ready-to-use liquid controls: QuikRead go HbA1c Control and QuikRead go HbA1c Control High. The control values have been determined for the QuikRead go instrument. The control is applied using the same procedure as for venous blood sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (1 µl). Analyse control samples by choosing Quality Control on the display of the QuikRead go instrument.

Quality control testing must be performed in accordance with applicable national or local regulations or accreditation requirements and laboratory standard Quality Control Procedures.

It is recommended to perform the control at least with each new shipment and new reagent lot and every time an unexpected result is obtained.

The control should be used each time new QuikRead go system users are trained.

## 9 Results' interpretation

The HbA1c result should not be interpreted without patient's medical history, clinical examination and other laboratory findings.

### HbA1c measuring units

Three different measuring units are used for reporting the results of the QuikRead go HbA1c test. The result will be displayed in mmol/mol HbA1c (aligned to the IFCC Reference Procedure<sup>8,9</sup>) or in % HbA1c (aligned to the DCCT study assay, also known as NGSP-HbA1c<sup>9</sup>) or in eAG (estimated average blood glucose), where a HbA1c value is converted to an equivalent glucose value. Linear relationship between these HbA1c results and estimated average glucose (eAG) have been established<sup>10</sup>. Patient's HbA1c values should be reported according to the national recommendations. Choose the right units from the instrument display: **Settings** -> **Measurement Flow** -> **Test Parameters**. Please see the instrument manual.

NGSP	=	0.0915 x IFCC + 2.15%
IFCC	=	10.93 x NGSP – 23.5 mmol/mol
eAG (mmol/l)	=	1.59 x DCCT-HbA1c (%) – 2.59
eAG (mg/dl)	=	28.7 x HbA1c (%) – 46.7 <sup>9</sup>

If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

## 10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics". Test results should never be used alone, without a complete clinical evaluation.

In conditions where the life span of red blood cells is shortened, the HbA1c results are decreased regardless of the method used. The life span is reduced for example in conditions such as hemolytic anemia or other hemolytic diseases, pregnancy, homozygous sickle cell trait, iron deficiency, blood loss, polycythemia, etc. Caution with interpreting HbA1c results in patients with these conditions is needed.

The HbA1c test should not be used to diagnose patients during pregnancy, patients that have received cancer chemotherapy within the past 3 weeks, patients that have received a blood transfusion within the past 3 weeks, patients with an elevated fetal hemoglobin such as hereditary persistence of fetal hemoglobin (HPFH) or patients with a hemoglobinopathy but normal red cell turnover (e.g. sickle cell trait). In case of rapidly evolving type 1 diabetes, the increase in HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions HbA1c value cannot be used for diagnosis of diabetes mellitus.

The test cannot be completed if hemoglobin in the sample is out of the measuring range (< 60 g/l).

## 11 Expected values

HbA1c test result	Interpretation of result <sup>7</sup>
≥ 48 mmol/mol (or 6.5%)	Indicates increased value of glycosylated hemoglobin (increased blood sugar) during the previous 3 months. Cut-off used for diagnosis of diabetes mellitus. Country-specific recommendations may vary.
39–47 mmol/mol (or 5.7–6.4%)	Values for intermediate hyperglycaemia.
< 39 mmol/mol (or 5.7%)	Indicates normal value of glycosylated hemoglobin. Does not exclude the possibility of patient having diabetes.

## 12 Performance characteristics

### Precision and reproducibility

A precision study was performed in line with the Clinical and Laboratory Standards Institute (CLSI) guideline EP05-A3:2014. QuikRead go HbA1c Control and Control High, and four venous EDTA whole blood samples were analyzed for 5 days, 2 runs per day and 2 replicates per run with 4 instruments. Precision data is summarized in Table 1. Multisite precision was performed at 3 independent NPT clinical sites using quality control material. The two controls were run in 5 replicates on 5 days. Multisite precision data is summarized in Table 2.

**Table 1.** Precision data. CV = Coefficient of Variation. N (Number of replicates) = 80

Sample	Mean Value (mmol/mol)	N	Repeatability		Within-Laboratory Precision	
			SD	%CV	SD	%CV
Low	35	80	0.7	2.1	0.9	2.6
Threshold	49	80	1.0	2.1	1.3	2.6
Medium	63	80	1.1	1.7	1.5	2.4
High	84	80	2.1	2.5	2.3	2.7
Control	37	80	0.6	1.7	0.7	1.9
Control High	112	80	1.8	1.6	2.2	1.9

**Table 2.** Multisite precision data from 3 independent NPT clinical sites. CV = Coefficient of Variation. N (number of replicates) = 75

Sample	Mean Value (mmol/mol)	N	Repeatability		Between-Site		Reproducibility	
			SD	%CV	SD	%CV	SD	%CV
Control	46	75	1.4	2.9	0.5	1.1	1.6	3.4
Control High	119	75	2.1	1.7	1.8	1.5	2.7	2.3

A precision study with finger prick samples was performed at 3 independent NPT clinical sites. Between-site precision was analyzed from 36 patients. Three (3) operators collected a finger prick sample from the same patient. The data is summarized in Table 3.

**Table 3.** Precision data of finger prick samples. CV = Coefficient of Variation. N (Number of replicates) = 9.

Sample	Mean Value (mmol/mol)	CV (%)
Low	35	2.3
Threshold	46	2.7
Medium	59	2.7
High	80	3.7

### Analytical specificity

Analytical specificity	Significant interference (≥ 10%)
HbA2	No
HbAC	Yes
HbAD	No
HbAE	No
HbAS	No
HbF	Yes (> 7%)
Acetylated Hb	No
Labile (pre-glycosylated) Hb	No

## Interference

Interfering substance	No significant interference ( $\leq 10\%$ ) found up to concentration
Biotin	3.51 mg/l
Bilirubin	200 mg/l
Glucose	10 g/l
Intralipids	20 g/l
Rheumatoid factor	780 000 IU/l
Albumin	150 g/l
Glycated albumin	7.7 g/l
Vitamin C (Ascorbic acid)	52.5 mg/l
Acetaminophen	156 mg/l
Acetylcysteine	150 mg/l
Acetylsalicylic acid	30 mg/l
Ampicillin	75 mg/l
Cefoxitin	6600 mg/l
Cyclosporine A	1.8 mg/l
Cyclosporine C	5.0 mg/l
Doxycycline	18 mg/l
Glyburide	0.72 mg/l
Heparin	3300 U/l
Ibuprofen	219 mg/l
Levodopa	7.5 mg/l
Metformin	12 mg/l
Methyldopa	22.5 mg/l
Metronidazole	123 mg/l
Phenylbutazone	321 mg/l
Rifampicin	48 mg/l
Salicylic acid	286 mg/l
Theophylline	60 mg/l

## Measuring range

The QuikRead go instrument displays the HbA1c results in mmol/mol and/or percentages (%). The calculated estimated average glucose (eAG) can be displayed with mmol/mol or % -results.

	HbA1c	HbA1c	eAG
Unit	mmol/mol	%	mmol/l
Measuring range	20–140	4.0–15.0	3.8–21.3

## 13 Traceability

The calibrators used to calibrate the HbA1c assay of the QuikRead go HbA1c test are traceable to the IFCC Reference Measurement Procedure<sup>8</sup>.

## 14 Method comparison

Finger prick capillary blood and anticoagulated whole blood samples from 170 patients were measured with the QuikRead go HbA1c method at a NPT clinical site. The reference values were determined with an HPLC method at clinical laboratory. Summary of correlation study is presented in Tables 4 and 5.

Table 4. Method comparison. QuikRead go HbA1c compared to a laboratory HPLC method.

Passing-Bablok analysis					
HPLC method					
Capillary sample			Venous blood sample		
$y = 1.02x - 1.8$			$y = 1.02x - 2.1$		
N = 170	Range = 29–130 mmol/mol	Correlation coefficient $r = 0.990$	N = 170	Range = 29–130 mmol/mol	Correlation coefficient $r = 0.990$

Table 5. Estimated difference at clinically significant HbA1c concentrations.

Expected value (mmol/mol)	Observed value (mmol/mol)	Difference (%)
42	41	-2.5
48	47	-1.9
53	52	-1.6

## 15 Disposal

- Dispose of contents according to national and local law.
- All patient samples and sample collectors, controls, used caps, cuvettes and tubes should be handled and disposed of as potentially infectious material.
- Materials used in various QuikRead go kit components:
  - Paper: Instructions for use
  - Cardboard: Kit box including its inner parts and sample collector box
  - Plastic: Cuvettes, reagent caps, cuvette rack, plungers, plunger and capillary tubes, filling inside capillary tube and the filling of sample collector box, swabs and extraction vials and tubes
  - Glass: Capillaries
  - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
  - Several (not to be recycled): Lids of reagent cap tubes, foil covering cuvette rack, unopened cuvettes and unused sample collectors
- When used in accordance with Good Laboratory Practice, good occupational hygiene, and the instructions for use, the reagents supplied should not present a hazard to health.

## 16 Troubleshooting

### Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the QuikRead go instrument instructions for use.

Error messages	Corrective actions
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner light violet part of the cap is not pressed down.
Reagent and cap lot do not match.	Discard the test. Perform a new test with new components. Make sure that cuvettes and caps are from the same kit lot.
Unable to read cap label.	Check from the cap tube that the cap is from the same kit than the cuvette. If the cap and cuvette are from the same lot, and the sample collector has not been over 5 minutes in the cuvette, you can continue the measurement. If lots are not the same, cancel the measurement, discard the test and perform a new test. Make sure that all the reagents are from the same kit lot.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25°C). Perform a new test with new components.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18...25°C). Perform a new test with new components.
Test cancelled. Temperature error.	Perform a new test. There has been some problem with the reaction temperature.
Test cancelled. Reagent error.	Perform a new test. There has been some problem with the reagents. If this error message appears again, please make sure that the HbA1c reagent caps have not been exposed to humidity. Also ensure that the reagent cap tube is properly closed.
No result. Sample error.	Perform a new test. There has been a problem with the sample. Ensure that the test has been done according to the instructions. The test cannot be completed if hemoglobin in the sample is out of the range (< 60 g/l).

### Unexpected low and high results







Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected high or low result.	Sample too long in the sample collector.	Perform a new test. Sample can be stored max. 1 minute in the sample collector.
	Too large or small sample volume.	Perform a new test. Take the sample by placing the tip of the sample collector just beneath the surface of the blood and fill the capillary part completely. Check that the sample collector is completely filled and there are no air bubbles in the capillary part. Ensure that there is no blood on the plastic frame of the sample collector and the sample has been taken just beneath the surface of the blood. Do not wipe the sample collector.
	Use of incorrect sample collector.	Perform a new test. Use QuikRead go Sample Collector 1µl for transferring the sample into the cuvette. Do not use pipette or other sampling apparatus.
	Use of components from different kit lots or tests.	Perform a new test. Make sure that all the reagents are from the same kit lot.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored in the manner described in the instructions for use.
Unexpected high result.	The cuvette is dirty.	Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette.
	Sample released in the solution before the analysis.	Perform a new test. Make sure that the sample is not released from the sample collector to the solution before inserting the cuvette in the instrument.

## References • Referenzen

1. Bunn F et al. The Biosynthesis of Human Hemoglobin A1c. Slow glycosylation of hemoglobin in vivo. J Clin Invest 1976; 57:1652-1659.
2. Beltran Del Rio, M., Tiwari, M., Amodu, L.I., Cagliani, J. & Rodriguez Rilo, H.L. 2016, "Glycated Hemoglobin, Plasma Glucose, and Erythrocyte Aging", Journal of diabetes science and technology, vol. 10, no. 6, pp. 1303-1307.
3. Sherwani SI, Khan HA, Ekhzaimy A, Masood A, Sakharkar MK. Significance of HbA1c Test in Diagnosis and Prognosis of Diabetic Patients. Biomark Insights 2016; 11: 95-104
4. Weykamp, C. 2013, "HbA1c: a review of analytical and clinical aspects", Annals of laboratory medicine, vol. 33, no. 6, pp. 393-400.
5. Sacks, D.B., Arnold, M., Bakris, G.L., Bruns, D.E., Horvath, A.R., Kirkman, M.S., Lernmark, A., Metzger, B.E. & Nathan, D.M. 2011, "Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus", Clinical chemistry, vol. 57, no. 6, pp. e1-e47.
6. European Diabetes Policy Group (1999) A desk-top guide to type 2 diabetes mellitus. Diabet Med 16:716–730.
7. International Expert Committee, International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. Diabetes Care 2009; 32:1327-1334.
8. Jeppsson, JO et al., Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, Clin Chem Lab Med 2002; 40(1):78-89.
9. Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine and International Diabetes Federation Consensus Committee. Diabetes Care 2007;30:2399-2400.
10. Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1C assay into estimated average glucose values. Diabetes Care 2008; 31:1473-1478.

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





- Bunn F et al. The Biosynthesis of Human Hemoglobin A1c. Slow glycosylation of hemoglobin in vivo. J Clin Invest 1976; 57:1652-1659.
- Beltran Del Rio, M., Tiwari, M., Amodu, L.I., Cagliani, J. & Rodriguez Rilo, H.L. 2016, "Glycated Hemoglobin, Plasma Glucose, and Erythrocyte Aging", Journal of diabetes science and technology, vol. 10, no. 6, pp. 1303-1307.
- Sherwani SI, Khan HA, Ekhzaimy A, Masood A, Sakharkar MK. Significance of HbA1c Test in Diagnosis and Prognosis of Diabetic Patients. Biomark Insights 2016; 11: 95-104
- Weykamp, C. 2013, "HbA1c: a review of analytical and clinical aspects", Annals of laboratory medicine, vol. 33, no. 6, pp. 393-400.
- Sacks, D.B., Arnold, M., Bakris, G.L., Bruns, D.E., Horvath, A.R., Kirkman, M.S., Lernmark, A., Metzger, B.E. & Nathan, D.M. 2011, "Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus", Clinical chemistry, vol. 57, no. 6, pp. e1-e47.
- European Diabetes Policy Group (1999) A desk-top guide to type 2 diabetes mellitus. Diabet Med 16:716-730.
- International Expert Committee, International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. Diabetes Care 2009; 32:1327-1334.
- Jeppsson, JO et al., Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, Clin Chem Lab Med 2002; 40(1):78-89.
- Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine and International Diabetes Federation Consensus Committee. Diabetes Care 2007;30:2399-2400.
- Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1C assay into estimated average glucose values. Diabetes Care 2008; 31:1473-1478.

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QuikRead go® is a registered trademark of Aidian Oy.