

## Disposal instructions of QuikRead go Instrument

**Product: 133893 QR Go Instrument**

**Date: 7.3.2025**

Aidian.eu website instructions under instrument disposal information and country specific additional info.

### Generic instruction

The following is generic instruction to be added to the Aidian.eu website regarding disposal and disinfection.

#### **How should a QuikRead go Instrument be disposed of?**

QuikRead go Instrument is a low-voltage electronic device. A used QuikRead go Instrument must be treated as potentially infectious waste. The instrument should be disposed of as electrical and electronic equipment (WEEE 2012/19/EU) if local and national laws do not require the instrument to be collected and disposed of as potentially infective clinical waste.

The packaging materials are recyclable.

Be sure to disinfect the device and to delete all patient data before disposing of the instrument.

#### **Disinfection**

Before sending QuikRead go Instrument to recycling clean the external surfaces and decontaminate the interior of the instrument.

Clean the exterior of the instrument using a lint-free cloth dampened with water. If necessary, a mild detergent may be used. Do not use organic solvents or corrosive substances. Spillage of potentially infectious material should be wiped off immediately with absorbent paper tissue and the contaminated areas swabbed with 70% ethyl alcohol or standard disinfectant. For safety, use chemical resistant gloves and follow the instructions in the safety data sheet. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

To decontaminate the inside of the device, open the lid and spray 5% sodium hypochlorite into the instrument through the cuvette hole or if not able to open lid, access interior of the instrument by removing the accumulator unit cover at the bottom of the instrument and opening the bottom hatch. Do not spray disinfectant on the accumulator or battery.

Let the disinfectant dry for 24 h and in addition, the device must be aged for two weeks after the last use the device before sending QuikRead go Instrument for recycling. For your safety, use chemical-resistant gloves, and follow the instructions in the safety data sheet of the disinfectant.



**Data protection**

Used QuikRead go Instrument may contain sensitive patient data. Please note, that the end user of the device is responsible for deleting all the data on the old devices to be disposed of.

To delete the data, go to Results > Delete result history.

**Please see country specific disposal information to correctly dispose your QuikRead go Instrument**