
153741-6

- Hrvatski
- Srpski
- Română
- English

QuikRead go[®] easy CRP

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Reference • Reference

1. Price CP et al. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med.* 1999; 37(2):109–13.
2. Pepys MB, Hirschfield GM. C-reactive protein: a critical update. *J Clin Invest* 2003; 111:1805-1812.
3. Bruns AH et al. Usefulness of consecutive C-reactive protein measurements in follow-up of severe community-acquired pneumonia. *Eur Respir J* 2008; 32:726-732.
4. Coelho L et al. Usefulness of C-reactive protein in monitoring the severe community-acquired pneumonia clinical course. *Crit Care* 2007; 11:R92.
5. Coelho LM et al. Patterns of C-reactive protein RATIO response in severe community-acquired pneumonia: a cohort study. *Crit Care* 2012; 16:R53.
6. Povoja P et al. C-reactive protein, an early marker of community-acquired sepsis resolution: a multi-center prospective observational study. *Crit Care* 2011; 15:R169.
7. Gans SL et al. Diagnostic value of C-reactive protein to rule out infectious complications after major abdominal surgery: a systematic review and meta-analysis. *Int J Colorectal Dis* 2015; 30:861-873.
8. Straatman J et al. Predictive Value of C-Reactive Protein for Major Complications after Major Abdominal Surgery: A Systematic Review and Pooled-Analysis. *Plos One* 2015; 10:e0132995.
9. Benitz WE et al. Serial serum C-reactive protein levels in the diagnosis of neonatal infection. *Pediatrics* 1998; 102:E41.
10. Ehl S et al. C-reactive protein is a useful marker for guiding duration of antibiotic therapy in suspected neonatal bacterial infection. *Pediatrics* 1997; 99:216-221.
11. Meem M et al. Biomarkers for diagnosis of neonatal infections: A systematic analysis of their potential as a point-of-care diagnostics. *J Glob Health* 2011; 1:201-209.
12. Aksu F et al. C-reactive protein levels are raised in stable Chronic obstructive pulmonary disease patients independent of smoking behavior and biomass exposure. *J Thorac Dis* 2013; 5:414-421.
13. Zhang Y et al. Association between C-reactive protein concentration and chronic obstructive pulmonary disease: a systematic review and meta-analysis. *J Int Med Res* 2012; 40:1629-1635.
14. van Durme YM et al. C-reactive protein levels, haplotypes, and the risk of incident chronic obstructive pulmonary disease. *Am J Respir Crit Care Med* 2009; 179:375-382.
15. Dessein PH et al. High sensitivity C-reactive protein as a disease activity marker in rheumatoid arthritis. *J Rheumatol.* 2004; 31 (6): 1095–7.
16. Poole CD et al. The association between C-reactive protein and the likelihood of progression to joint replacement in people with rheumatoid arthritis: a retrospective observational study. *BMC Musculoskelet Disord.* 2008; 9: 146.
17. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples (WHO/DIL/LAB/99.1 Rev.2), 2002. Available: http://whqlibdoc.who.int/hq/2002/WHO_DIL_LAB_99.1_Rev.2.pdf. Accessed 21.8.2014.
18. Thomas L. C-reaktives Protein (CRP) in: Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780–90.
19. Chiesa C et al. C reactive protein and procalcitonin: Reference intervals for preterm and term newborns during the early neonatal period. *Clin Chim Acta.* 2011; 412 (11-12): 1053–9.
20. Shine et al. Solid phase immunoassays for human C-reactive protein. *Clin Chim Acta.* 1981; 117: 13-23.
21. Schlebusch H et al. High Sensitive CRP and Creatinine: Reference intervals from Infancy to Childhood. *J Lab Med* 2002; 26 (5/6): 341-6.

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
1. Price CP et al. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med*. 1999; 37(2):109–13.
2. Pepys MB, Hirschfield GM. C-reactive protein: a critical update. *J Clin Invest* 2003; 111:1805-1812.
3. Bruns AH et al. Usefulness of consecutive C-reactive protein measurements in follow-up of severe community-acquired pneumonia. *Eur Respir J* 2008; 32:726-732.
4. Coelho L et al. Usefulness of C-reactive protein in monitoring the severe community-acquired pneumonia clinical course. *Crit Care* 2007; 11:R92.
5. Coelho LM et al. Patterns of C-reactive protein RATIO response in severe community-acquired pneumonia: a cohort study. *Crit Care* 2012; 16:R53.
6. Povoia P et al. C-reactive protein, an early marker of community-acquired sepsis resolution: a multi-center prospective observational study. *Crit Care* 2011; 15:R169.
7. Gans SL et al. Diagnostic value of C-reactive protein to rule out infectious complications after major abdominal surgery: a systematic review and meta-analysis. *Int J Colorectal Dis* 2015; 30:861-873.
8. Straatman J et al. Predictive Value of C-Reactive Protein for Major Complications after Major Abdominal Surgery: A Systematic Review and Pooled-Analysis. *Plos One* 2015; 10:e0132995.
9. Benitz WE et al. Serial serum C-reactive protein levels in the diagnosis of neonatal infection. *Pediatrics* 1998; 102:E41.
10. Ehl S et al. C-reactive protein is a useful marker for guiding duration of antibiotic therapy in suspected neonatal bacterial infection. *Pediatrics* 1997; 99:216-221.
11. Meem M et al. Biomarkers for diagnosis of neonatal infections: A systematic analysis of their potential as a point-of-care diagnostics. *J Glob Health* 2011; 1:201-209.
12. Aksu F et al. C-reactive protein levels are raised in stable Chronic obstructive pulmonary disease patients independent of smoking behavior and biomass exposure. *J Thorac Dis* 2013; 5:414-421.
13. Zhang Y et al. Association between C-reactive protein concentration and chronic obstructive pulmonary disease: a systematic review and meta-analysis. *J Int Med Res* 2012; 40:1629-1635.
14. van Durme YM et al. C-reactive protein levels, haplotypes, and the risk of incident chronic obstructive pulmonary disease. *Am J Respir Crit Care Med* 2009; 179:375-382.
15. Dessein PH et al. High sensitivity C-reactive protein as a disease activity marker in rheumatoid arthritis. *J Rheumatol*. 2004; 31 (6): 1095–7.
16. Poole CD et al. The association between C-reactive protein and the likelihood of progression to joint replacement in people with rheumatoid arthritis: a retrospective observational study. *BMC Musculoskelet Disord*. 2008; 9: 146.
17. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples (WHO/DIL/LAB/99.1 Rev.2), 2002. Available: http://whqlibdoc.who.int/hq/2002/WHO_DIL_LAB_99.1_Rev.2.pdf. Accessed 21.8.2014.
18. Thomas L. C-reaktives Protein (CRP) in: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780–90.
19. Chiesa C et al. C reactive protein and procalcitonin: Reference intervals for preterm and term newborns during the early neonatal period. *Clin Chim Acta*. 2011; 412 (11-12): 1053–9.
20. Shine et al. Solid phase immunoassays for human C-reactive protein. *Clin Chim Acta*. 1981; 117: 13-23.
21. Schlebusch H et al. High Sensitive CRP and Creatinine: Reference intervals from Infancy to Childhood. *J Lab Med* 2002; 26 (5/6): 341-6.

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1 Intended use

The QuikRead go easy CRP test is an immunoturbidimetric assay for quantitative determination of C-reactive protein (CRP) values in whole blood, serum and plasma. The test is carried out using the QuikRead go instrument. The assay is intended to be used by trained health care professionals in the clinical laboratory and Point-of-Care (POC) settings.

CRP is an acute-phase protein present in low concentrations in the blood of healthy individuals.¹ Measurement of CRP aids in the detection and evaluation of infection, tissue injury and inflammatory disorders. It provides information for the diagnosis, therapy and monitoring of inflammatory disorders.

For *in vitro* diagnostic use only.

2 Summary and explanation of the test

Infections, tissue injuries and inflammatory disorders may be accompanied by a nonspecific elevation of the CRP level in the patient's blood. The rise in CRP levels is rapid, and increased levels can be detected within 6 hours of the onset of the inflammatory process.²

Quantitative measurement of CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections³⁻⁶, as well as an effective tool in controlling and monitoring postoperative infections^{7,8}. Modestly elevated CRP values have been shown to be clinically relevant e.g. in neonatal sepsis⁹⁻¹¹, and conditions that include systemic inflammation (such as COPD¹²⁻¹⁴ and Rheumatoid arthritis¹⁵⁻¹⁶).

The advantages of measuring CRP lie in it being a sensitive, early and reliable reflector of the acute phase response. It can also be used for monitoring of the patient's condition and response to therapy. However, the rise in CRP level is non-specific.

3 Principles of the procedure

The CRP measurement of the QuikRead go easy CRP test is an immunoturbidimetric test based on the agglutination reaction. Microparticles are coated with anti-human CRP F(ab)₂ fragments, and the CRP present in the sample reacts with microparticles. The resultant change in the turbidity of the solution is measured by the QuikRead go instrument. The CRP assay correlates well with results obtained from several other immunoturbidimetric methods.

The sample is added to the buffer in prefilled cuvettes and the blood cells in a whole blood sample are haemolysed. The calibration data of the test are found in the barcode label of the cuvette, which the QuikRead go unit reads automatically before the start of the test. The CRP value is automatically corrected on the basis of the sample's hematocrit level. Hematocrit results are not displayed but are used in the calculations. The CRP measurement range of the assay is 1–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 1–120 mg/l with serum or plasma samples.

4 Reagents

Content of the kit

| Component name and origin | Symbol | QuikRead go® easy CRP Cat. no. 153287 50 tests |
|---|-------------------|--|
| QuikRead go easy CRP Reagent Caps ORIG SHP | REAG CPS | 2 x 25 |
| Buffer | BUF | 2 x 25 x 1 ml |
| QuikRead go Sample Collector 10µl (sample collector) | SAMPL COL EDTA | 50 |
| Instructions for use | | |

The reagents contain preservatives, please see Section 5. "Warnings and precautions".

Storage

| Kit component | Opened at 2...8°C | Opened at 18...25°C | Unopened at 2...25°C |
|------------------|--|-------------------------|--------------------------------------|
| Reagent Caps | 6 months | 6 months | Until the expiration date of the kit |
| Buffer | Foil pouch: 6 months Single prefilled cuvette: 2 hours | Foil pouch: 3 months | |
| Sample collector | Storage at 2...25°C until the expiry date of QuikRead go Sample Collector 10µl box | | |

Reagent preparation and storage conditions

All reagents are ready to use. Keep the CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

Reagent deterioration

Do not use a cuvette with visible dirt in the buffer.

5 Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink on premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagents inside a reagent cap contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Section 15.

Warning

Analytical precautions

- Do not use the product after the expiration date marked on the outer package.
- Do not exceed the stated stability periods for opened reagents.
- The QuikRead go easy CRP kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (the optical part). Discard any cuvettes with fingerprints.
- The QuikRead go easy CRP Reagent Caps have turquoise color.
- Keep the QuikRead go easy CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Do not splash liquid into the measuring well of the instrument.

6 Sample collection and preparation

Sample material, collection and volume

QuikRead go Sample Collector 10µl (Cat. no. 153211) included in the kit should be used for adding the sample into the cuvette. The capillary part of the sample collector is covered with EDTA and the appearance might vary depending on crystallization of the EDTA.

| Sample material | Sample volume | Sample collection |
|--------------------------------|---------------|---|
| Fingertip sample (whole blood) | 10 µl | Prick the clean and dry finger with a lancet. Wipe off the first drop and collect 10 µl of blood from the second drop into EDTA covered capillary part of the sample collector. |
| Anticoagulated whole blood | 10 µl | Use a venous blood sample collected in a tube containing heparin or EDTA. Mix whole blood by inverting the tube several times and take 10 µl with the sample collector. |
| Plasma | 10 µl | Use EDTA/heparin plasma. Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. |
| Serum | 10 µl | Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. |

Sample dilution

| Sample material | Instructions |
|-----------------|--|
| Whole blood | Do not dilute whole blood samples. |
| Plasma/serum | Plasma or serum samples may be diluted with a 0.9% NaCl solution before the sample is added to the cuvette. The recommended dilution ratio is 1+3 (by volume, 1 part sample + 3 parts 0.9% NaCl). Take 10 µl of the diluted sample into a cuvette with the sample collector, analyse the sample and multiply the obtained result by 4. |

Sample storage

| Sample material | Short-term storage | Long-term storage |
|---|--|--|
| Fingertip blood (whole blood) | Max. 2 minutes in the sample collector | Do not store |
| Anticoagulated whole blood | 2...8°C for 3 days | Separate plasma and store as indicated below |
| Plasma | 2...8°C for 7 days | Below -20°C for 3 years ¹⁷ |
| Serum | 2...8°C for 7 days | Below -20°C for 3 years ¹⁷ |
| Sample (whole blood, plasma, serum) in buffer | 18...25°C for 2 hours | Do not store |

Allow the samples to reach room temperature (18...25°C) before testing. Frozen samples must be thawed completely, mixed thoroughly, and brought to room temperature prior to testing. If the thawed plasma samples contain clots, the samples should be centrifuged. Samples should not be repeatedly frozen and thawed.

7 Procedure

Materials required but not provided

The following materials are required to perform the test but are not provided or they are available also separately. The materials provided are listed in Section 4. "Reagents".

| Material | Cat. no | Additional info |
|--|----------------|---------------------------------|
| QuikRead go® Instrument with software version 7.6.2 or newer | 133893, 149915 | – |
| QuikRead go® CRP Control Low | 153765 | Recommended for quality control |
| QuikRead go® CRP Control | 153764 | Recommended for quality control |
| QuikRead go® CRP Control High | 153763 | Recommended for quality control |
| QuikRead go® Sample Collector 10µl | 153211 | Supplied with 153287 |
| Fingertip lancets | – | – |

Assay procedure

Before starting an assay, read and follow the QuikRead go instrument and QuikRead go easy CRP test instructions for use.

Open the foil pouch of the cuvette rack, and mark the opening date on the cuvette rack.

Sample collection (see figures 1–6)

Sample collectors (10 µl) for QuikRead go easy CRP should be used for adding the sample to a prefilled cuvette. The QuikRead go easy CRP kit (Cat. no. 153287) contains the sample collectors.

1. A prefilled cuvette must reach room temperature (18...25°C) before use. This takes 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces of the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on the results. The test must be performed within two (2) hours of opening the cuvette. The liquid surface should lie between the two lines marked on the cuvette.
2. Open the sample collector box slightly to get one single the sample collector out of the box. The package filling can be discarded.
3. Hold the sample collector almost horizontally and take the sample by placing the tip of the sample collector just beneath the surface of the blood drop. Fill the capillary part of the sample collector completely.
4. Check that the sample collector is completely filled and there are no air bubbles in the capillary. Do not wipe the sample collector. If there is excess of sample on the sample collector, discard the sample collector and take a new sample with a new sample collector.
5. Within 2 minutes insert the sample collector into the cuvette. The sample collector remains in the cuvette during the whole test measurement procedure.
6. Close the cuvette with the reagent cap. Do not press down the turquoise inner part of the reagent cap. Once the sample has been added to the buffer the assay should be run within two hours. Keep the cuvette upright and do not shake.

Note! If a sample tube is used, it is recommended to take the sample directly from the tube. Tilt the tube slightly and collect the sample from just beneath the surface of the sample with the sample collector. Do not wipe the sample collector.

Analysing the sample (see figures 7–9)

Begin the test by turning on the instrument. For full instructions, see the QuikRead go instrument instructions for use. The display will guide you through the test.

7. Choose **Measure** on the display of the QuikRead go instrument.
8. Place the cuvette in the instrument with the barcode facing you. The instrument detects the sample type. The display shows how the measurement is proceeding.
9. The result appears on the display when the measurement is completed and the cuvette automatically rises from the measurement well.

Analysing the control

Please see the instructions for use provided with each control vial. Analyse the control samples as patient samples but select **Quality Control** on the display of the QuikRead go instrument. The result will be stored as a control measurement.

8 Quality control

Regular use of QuikRead go CRP Control Low (Cat. no 153765), QuikRead go CRP Control (Cat. no 153764) and/or QuikRead go CRP Control High (Cat. no 153763) is recommended. Please consult Section 16. "Troubleshooting" and contact Aidian or your local distributor if you observe any changes in test functionality.

• QuikRead CRP controls:

The QuikRead go CRP Control Low, QuikRead go CRP Control and QuikRead go CRP Control High are ready-to-use controls. The control values have been determined for the QuikRead go instrument and the control is applied using the same procedure as for a plasma/serum sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (10 µl). The control is stable for 15 minutes in the cuvette.

• Other commercially available CRP controls:

Handle and use the controls in accordance with the instructions for each control substance. Perform the test as outlined in Section 7. "Procedure". Acceptable limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells, because these might not haemolyse normally.

Quality control testing must be performed in accordance with applicable national or local regulations or accreditation requirements and laboratory standard Quality Control Procedures. It is recommended to perform the control at least with each new shipment and new reagent lot and every time an unexpected result is obtained. The control should be used each time new QuikRead go system users are trained.

9 Results' interpretation

A rise in the CRP level is non-specific and should not be interpreted without a complete clinical history. Serial CRP measurements may be required to monitor a patient's condition or response to therapy and when measuring CRP levels in newborns. Additionally, establishing the patient-specific baseline CRP level is recommended when assessing slightly elevated CRP levels.

| CRP test result | Interpretation of result |
|-----------------|---|
| < 10 mg/l | <ul style="list-style-type: none"> Excludes many acute inflammatory diseases but does not specifically exclude inflammatory processes.¹⁸ Possible inflammations in newborn babies are not necessarily excluded.⁹⁻¹¹ Possible inflammations or exacerbations in conditions with systemic inflammation are not necessarily excluded.¹³⁻¹⁶ |
| 10–50 mg/l | Elevated concentrations are found in acute disease that occur in the presence of a slight to moderate inflammatory process. ¹⁸ |
| > 50 mg/l | Indicates high and extensive inflammatory activity. ¹⁸ |

If diluted plasma or serum samples have been used, the final results must be calculated in accordance with the instructions in Section 6. "Sample Collection and Preparation". If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

11 Expected values

Reference value: <3 mg/l¹⁹⁻²¹, median value in apparently healthy adults: 0.8 mg/l²⁰. CRP levels in newborns can be affected by e.g. prematurity, duration of delivery and complications related to the delivery.

Immediately after birth the CRP levels in term neonates are generally low and will increase to some extent (median below 2 mg/l) over the next few days¹⁹.

12 Performance characteristics

Precision and reproducibility

A precision study was performed in line with the Clinical and Laboratory Standards Institute (CLSI) guidelines, EP5-A3.

| Within-run, between-day and total precision | | | | | | |
|---|---------------|----------------|-----------------|-------------------|--------------------|--------------|
| Sample material | Sample number | Number of days | Mean CRP (mg/l) | Within-run CV (%) | Between-day CV (%) | Total CV (%) |
| Whole blood | Sample 1 | 5 | 14 | 2.9 | 0.9 | 3.2 |
| | Sample 2 | 5 | 60 | 3.8 | 1.5 | 4.7 |
| | Sample 3 | 5 | 175 | 8.0 | 2.2 | 8.3 |
| Plasma | Sample 1 | 20 | 1.1 | 6.8 | 1.7 | 7.8 |
| | Sample 2 | 20 | 44 | 3.4 | 1.4 | 3.7 |
| | Sample 3 | 20 | 96 | 8.8 | 0.5 | 8.9 |
| Control 1 | Sample 1 | 20 | 2.1 | 2.3 | 1.3 | 4.7 |
| Control 2 | Sample 1 | 20 | 24 | 4.8 | 0.9 | 5.3 |
| Control 3 | Sample 1 | 20 | 80 | 5.2 | 4.2 | 8.1 |

Interference

| Interfering substance | No interference found up to concentration |
|-------------------------------------|---|
| Bilirubin | 400 µmol/l |
| Vitamin C | 200 µmol/l |
| Triglycerides | 11.5 mmol/l |
| Cholesterol | 10.0 mmol/l |
| Rheumatoid factors (RF) | 525 IU/ml |
| Leucocytes | 87x10 ⁹ cells/l |
| Anticoagulants (Li-heparin or EDTA) | No interference |

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test because the assay antibodies lack the Fc part. In rare cases, interference by IgM myeloma protein has been observed.

Limits of detection

The detection limit of the assay is < 0.8 mg/l.

Measurement range

For whole blood samples, the measurement range for CRP is 1–200 mg/l at the hematocrit level of 40%. If the hematocrit is lower or higher than 40%, the measurement range changes in the manner shown in the table below. The CRP result is not displayed if the level of hematocrit is outside the range 15–75%.

| Sample type | Hematocrit % | Sample volume | Measurement range (mg/l CRP) | |
|---------------------|--------------|---------------|------------------------------|---|
| Whole blood sample | 15–19 | 10 µl | 1.0–140 | In cases when the result is outside the measurement range, the result will be displayed as "CRP < 1,0 mg/l" or, for instance, as "CRP > 200 mg/l" |
| | 20–24 | | 1.0–150 | |
| | 25–29 | | 1.0–160 | |
| | 30–33 | | 1.0–170 | |
| | 34–36 | | 1.0–180 | |
| | 37–39 | | 1.0–190 | |
| | 40–42 | | 1.0–200 | |
| | 43–45 | | 1.1–210 | |
| | 46–50 | | 1.2–220 | |
| | 51–53 | | 1.3–240 | |
| | 54–57 | | 1.4–260 | |
| | 58–60 | | 1.5–280 | |
| | 61–64 | | 1.7–300 | |
| | 65–68 | | 1.9–340 | |
| 69–71 | 2.1–380 | | | |
| 72–75 | 2.4–420 | | | |
| Plasma/serum sample | – | 10 µl | 1.0–120 | |

It is possible to turn off the sliding result range, in which case the fixed measurement range for whole blood samples is 2.4–140 mg/l. In order to turn off the sliding measurement range, you have to change the settings of the QuikRead go instrument. The change is made to easy CRP test parameters. Please see the instrument manual.

Antigen excess

CRP concentrations of less than 1000 mg/l do not give falsely low results.

Linearity

The recoveries in serial dilutions of CRP plasma samples have ranged from 85% to 104%.

13 Traceability

The calibrators used to calibrate the CRP assay of the QuikRead go easy CRP test are traceable to the ERM[®]-DA474 reference material.

14 Method comparison

Patient plasma samples were analysed via clinical laboratory method, a point-of-care test and the QuikRead go easy CRP method. A summary of correlation studies is presented in the table below.

| Passing-Bablok analysis | |
|--|--|
| Clinical laboratory method | Point-of-care test |
| $y = 1.00x - 0.5$ $r = 0.988$ $n = 92$ | $y = 0.93x + 1.2$ $r = 0.982$ $n = 89$ |

Whole blood versus plasma comparability

In comparison of 95 clinical samples, whole blood (y) and plasma (x) yielded comparable results: $y = 0.99x - 0.4$

15 Disposal

- Dispose of contents according to national and local law.
- All patient samples, sample collectors, controls, used tubes, caps and cuvettes should be handled and disposed of as potentially infectious material.
- Materials of the components of the QuikRead go easy CRP kit:
 - Paper: Instructions for use
 - Cardboard: Kit and sample collectors boxes including inner parts
 - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, sample collector frame, cuvette rack, sample collector package filling
 - Glass: Capillary part of sample collectors
 - Metal: Reagent cap tubes, cuvette lids
 - Not to be recycled: Lids of reagent cap tubes (several)
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

16 Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the QuikRead go instrument instructions for use.

| Error messages | Corrective actions |
|---|---|
| Measurement prohibited. Please check the reagent cap. | Check that the cuvette has a reagent cap and the inner turquoise part of the cap is not pressed down. |
| Reagent and cap lots do not match. | Discard the test. Perform a new test with new components. Make sure that all the reagents are from the same kit lot. |
| Unable to read cap label. | Check from the cap tube that the cap is from the same kit than the cuvette. If the cap and cuvette are from the same lot, you can continue the measurement. If lots are not the same, cancel the measurement, discard the test and perform a new test. Make sure that all the reagents are from the same kit lot. |
| Measurement prohibited. Cuvette temperature too low. | Let the cuvette warm up to room temperature (18...25°C). Test the same cuvette again. |
| Measurement prohibited. Cuvette temperature too high. | Let the cuvette cool down to room temperature (18...25°C). Test the same cuvette again. |
| Test cancelled. Blank too high. | Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed. |
| Test cancelled. Unstable blank. | |
| Test cancelled. Error in reagent addition. | Perform a new test. There has been some problem during reagent addition. If this error message appears again, please make sure that the easy CRP reagent caps have not been exposed to humidity. Also ensure that the reagent cap tube is properly closed. |
| No result. Hematocrit too low. | If the level of hematocrit in a whole blood sample is below the hematocrit measurement range of the instrument (see section 12), no result is shown. |
| No result. Hematocrit too high. | If the level of hematocrit in whole blood samples is above the hematocrit measurement range of the instrument (see section 12), no result is shown. |

Unexpected low and high results

Possible reasons for unexpected low and high results are listed in the table below.

| Problem | Possible cause | Corrective action |
|-----------------------------|---|---|
| Unexpected low result. | Too small sample volume. | Perform a new test. Make sure that the sample collector is completely filled and there are no air bubbles in the capillary. It is not possible to overfill the capillary. Do not wipe the sample collector. |
| | Sample too long in the sample collector. | Perform a new test. Sample can be stored max. 2 minutes in the sample collector. |
| Unexpected high result. | Too large sample volume. | Perform a new test. Take the sample by placing the tip of the sample collector just beneath the surface of the sample and fill the capillary part completely. Do not wipe the sample collector. Make sure that there is no sample on the turquoise part of the sample collector and the sample has been taken just beneath the surface of the sample. |
| | Too low buffer volume. | Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette. |
| | The cuvette is dirty. | Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette. |
| Unexpected low/high result. | Use of components from different kit lots or tests. | Perform a new test. Make sure that all the reagents are from the same kit lot. |
| | Incorrect reagent storage. | Perform a new test. Make sure that the reagents are stored in the manner described in the instructions for use. |

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QuikRead go[®] is a registered trademark of Aidian Oy.