Intended use

QuikRead go iFOBT Control Quantitative is intended for quality control of the QuikRead go iFOBT assay 151051.

For *in vitro* diagnostic use.**Principles of the procedure**

QuikRead go iFOBT Control Quantitative contains human haemoglobin (Hb) in buffer. The solution is suitable for controlling the quantitative Hb determination in concentrations comparable to occult blood in faeces. The Hb present in the solution reacts with the microparticles in FOB reagent, and the resultant change in turbidity is measured by the QuikRead go instrument.

Kit contents

QuikRead go® iFOBT Control Quantitative	Cat. no. 152390
QuikRead go iFOBT Control Quantitative	CONTROL ORIG HUM 2 ml
QuikRead FOB Pastette	18 pcs
Instructions for use	1

Storage

Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 1 month at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

Preparation

QuikRead go iFOBT Control Quantitative is ready-to-use.

Use of the control

- Take one QuikRead FOB pastette.
- Hold the pastette by the upper bulb (A), squeeze the bulb firmly before inserting the tip (C) into the control liquid (picture 1).

- Release the bulb, lift the pastette and make sure that the liquid has filled the entire tip (picture 2).
- Release the liquid from the pastette into the test cuvette by squeezing the upper bulb firmly (picture 3). All of the liquid in the tip (C) should be released, but some liquid may remain in the lower bulb (B). Do not squeeze the lower bulb at any point.
- Discard the used pastette.
- Continue with the assay procedure, as described in the instructions for use of QuikRead go iFOBT.

Expected values

The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.

Warnings and precautions

- Consult the warnings and precautions in the instructions for use of QuikRead go iFOBT.
- QuikRead go iFOBT Control Quantitative contains material of human origin that has been tested and found negative for anti-human immunodeficiency virus (anti-HIV1&2), anti-hepatitis C virus (anti-HCV) antibodies, and hepatitis B surface antigen (HBsAg). Handle the product as if it is capable of transmitting infectious agents.
- The reagent contains <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal of all test materials should be in accordance with national and local law.

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