Revision History/ Revisionshistorie/ Historique des révisions/ Historial de revisión/ Revisiegeschiedenis/ Versiooni ajagulu/ Historie verzí/ História verzí/ Historia zmian/ Revisionshistorik/ Revisionshistorik/ Revisjonshistorikk/ Versiohistoria

<p>Version and publication date/ Version und Veröffentlichungsdatum/ Version et date de publication/ Fecha de versión y publicación/ Versie en releasedatum/ Versioon ja väljalaskekuupäev/ Verze a datum vydání/ Verzia a dátum vydania/ Wersja i data wydania/ Version och utgivningsdatum/ Version og udgivelsesdato/ Versjon og utgivelsesdato/ Versionumero sekä julkaisuajankohta</p>	<p>Change classification/ Klassifizierung ändern/ Classification de la modification/ Clasificación de cambio/ Classificatie wijzigen/ Muuda klassifikatsiooni/ Klasifikace změny/ Klasifikácia zmeny/ Zmiana klasyfikacji/ Ändringsklassificering/ Muutoksen luokitus/ Ændring af klassifikation/ Endringsklassifisering</p>	<p>Section of change/ Abschnitt der Änderung/ Section de la modification/ Sección de cambio/ Sectie van wijziging/ Muutuste jagu/ Změněná část/ Zmenená časť/ Sekcja zmiany/ Endringsavsnitt/ Ändrat avsnitt/ Afsnit af ændring/ Kappale, jossa muutos esiintyy</p>
<p>151054-5 02/2023</p>	<p>Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Cambio normativo/ Regulatieve wijziging/ Regulatiivne muudatus/ Regulační změna/ Regulačná zmena/ Zmiana regulacyjna/ Regulatorisk ändring/ Regulatorisk ændring/ Regulatorisk endring/ Säätelyn muutos</p>	<p>Explanation of symbols/ Zeichenerklärung/ Explication des symboles/ Explicación de los símbolos/ / Verklaring van symbolen/ Sümbolite seletus/ Vysvětlení symbolů/ Vysvetlivky symbolov/ Symbolförklaring/ Symbolforklaring/ Symbolforklaringer/ Symbolien selitykset</p>

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English

Intended purpose

QuikRead iFOB Sampling Set is intended to collect a specified amount of faeces and transfer the diluted and filtered faecal sample to enable QuikRead go iFOBT (a faecal immunochemical test, FIT) to be used as intended. Each sampling tube in QuikRead iFOB Sampling Set consists of a sample collecting stick, dilution buffer, preservative, and a filter. For *in vitro* diagnostic use only.

Sampling instructions

Before starting the sample collection, read and follow these instructions and consult the separate pictorial instructions.

- Obtain a sample of faeces in a clean container or on clean collecting paper (not provided). Avoid contact with the toilet water when you are collecting the sample. If there is visible blood in the faeces sample or in case of diarrhea, do not collect a sample but wait for the next bowel movement and collect a sample from that.
- Unscrew the top of the sampling tube containing the sample collecting stick. Collect the sample by poking the sample collecting stick (attached to the tube cap) in three (3) different places on the faeces sample.
- Insert the sample collecting stick back into the sampling tube and screw the cap on firmly.
- Write patient number and sample collection date on a separate label and attach the label to the sampling tube.
- Return the iFOB sampling tube to the laboratory or doctor as soon as possible after sampling. Protect the sampling tube from breakage or leakage during transportation. The sample in the sampling tube can be stored in a refrigerator (2...8°C) for up to 7 days or at room temperature (max. 25°C) for up to 5 days. Storage in refrigerator is however preferred.

Storage

Store the sampling tubes at 2...25°C.

Warnings and precautions

- The sample collection and transfer to the test are intended to be performed by healthcare professionals in the clinical laboratory and near-patient testing (NPT) setting (also known as point of care setting). The patient may also self-collect the sample, in which case the healthcare professional is expected to provide adequate guidance to the patient.
- For use with QuikRead go iFOBT kit, Ref. 151051.
- Do not drink the liquid in the sampling tube.
- Take care not to bend or break the tip of the tube.
- Do not use the sampling tubes after the expiry time marked on the label as yyyy-mm.
- The reagent contains 0.1% formaldehyde. May cause cancer (Carc. 1B, H350) and may produce an allergic reaction (EUH208). Obtain special instructions before use (P201). Wear protective gloves/protective clothing (P280). IF exposed or concerned: Get medical advice/attention (P308+P313). Dispose of contents according to national and local law (P501). Restricted to professional users.
- Reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates very toxic gas when in contact with acids (EUH032). Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- All patient samples should be handled and disposed of as potentially infectious material.
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, reagents supplied should not present a hazard to health.
- In case of a serious incident, please report it to the manufacturer or its representative and/or national authority.

Danger

Revision history

Revision history can be found from aidian.eu.

Consult instructions for use: QuikRead go iFOBT, Ref. 151051

QuikRead go® is a registered trademark of Aidian Oy.

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