
146364-10







- English
- Deutsch
- Français
- Italiano

# QuikRead go<sup>®</sup>

## wrCRP+Hb

### Assay Procedure • Testverfahren

8	Measure	RESULT
	CRP:	1,5 mg/l
	Hb:	125 g/l
Patient ID:	XXXXXXXX	Measurement time: 2015-03-25 12:19
Test:	wrCRP+Hb	Result info
	Choose <b>Result info</b> to view result information. Remove cuvette to perform a new measurement.	
	Exit	Print
		New measurement

## 1 Intended use

The QuikRead go wrCRP+Hb test is an immunoturbidimetric assay for quantitative determination of C-reactive protein (CRP) values in whole blood, serum and plasma and of haemoglobin (Hb) values in whole blood. The test is carried out using the QuikRead go instrument.

CRP is an acute-phase protein present in low concentrations in the blood of healthy individuals.<sup>1</sup> Measurement of CRP aids in the detection and evaluation of infection, tissue injury and inflammatory disorders. It provides information for the diagnosis, therapy and monitoring of inflammatory disorders.

Hb is the iron-containing protein in red blood cells that carries oxygen. Determining its concentration is a common procedure in both primary health care and acute care.

For *in vitro* diagnostic use only.

## 2 Summary and explanation of the test

### **C-reactive protein (CRP)**

Infections, tissue injuries and inflammatory disorders may be accompanied by a non-specific elevation of the CRP level in the patient's blood. The rise in CRP levels is rapid, and increased levels can be detected within 6–12 hours of the onset of the inflammatory process.<sup>2</sup>

Quantitative measurement of CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections, as well as an effective tool in controlling and monitoring postoperative infections<sup>2-6</sup>. Modestly elevated CRP values as well as patient-specific variations within the range of 0.5–5 mg/l have been shown to be clinically relevant e.g. in neonatal sepsis<sup>7-9</sup>, recurrent coronary events<sup>10-12</sup> and conditions that include systemic inflammation (such as COPD<sup>13-15</sup> and Rheumatoid arthritis<sup>16,17</sup>).

The advantages of measuring CRP lie in it being a sensitive, early and reliable reflector of the acute phase response. It can also be used for accurate monitoring of the patient's condition and response to therapy. However, the rise in CRP level is non-specific.

### **Haemoglobin (Hb)**

A decrease in the Hb concentration may result from iron deficiency or loss of blood. An elevated concentration of Hb in blood may indicate a low oxygen supply, e.g. in those living in mountainous regions or in smokers.

Determining the haemoglobin level rapidly gives a general picture of the individual's health/overall physical condition. However, the reasons for changes in haemoglobin levels alone need to be investigated and measurement of other blood values is recommended. Haemoglobin levels may also be affected by the individual's plasma volume.

## 3 Principles of the procedure

### **C-reactive protein (CRP)**

The CRP measurement of the QuikRead go wrCRP+Hb test is an immunoturbidimetric test based on the agglutination reaction. Microparticles are coated with anti-human CRP F(ab)<sub>2</sub> fragments, and the CRP present in the sample reacts with microparticles. The resultant change in the turbidity of the solution is measured by the QuikRead go instrument. The CRP assay correlates well with results obtained from several other immunoturbidimetric methods.

The sample is added to the buffer in prefilled cuvettes and the blood cells in a whole blood sample are haemolysed. The calibration data of the test are found in the barcode label of the cuvette, which the QuikRead go unit reads automatically before the start of the test. The CRP value is automatically corrected on the basis of the sample's hematocrit level. Hematocrit results are not displayed but are used in the calculations. The CRP measurement range of the assay is 0.5–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 0.5–120 mg/l with serum or plasma samples.

### **Haemoglobin (Hb)**

The Hb measurement using the QuikRead go wrCRP+Hb test is based on measuring oxyhaemoglobin photometrically at two wavelengths. The buffer haemolyses the sample's red blood cells in the cuvette, whereupon the absorption of released Hb can be measured. The instrument corrects for the effect of lipids or other turbidity-inducing factors on the results. The barcode on the cuvette contains the required parameters for calculating the Hb concentration, and the QuikRead go instrument displays the obtained Hb results. The Hb measuring range of the Hb assay is 50–220 g/l (5.0–22.0 g/dl, 3.1–13.7 mmol/l).

## 4 Reagents

### Content of the kit

Component name and origin	Symbol	QuikRead go® wrCRP+Hb	
		Cat. no. 146255 50 tests	Cat. no. 152633 50 tests
QuikRead go wrCRP+Hb Reagent Caps ORIG SHP	REAG CPS	2 x 25	2 x 25
Buffer	BUF	2 x 25 x 1 ml	2 x 25 x 1 ml
Capillaries (10 µl)	CAPIL HEP	50	–
Plungers	PLUN	50	–
Instructions for use			

The reagents contain preservatives, please see Section 5 "Warnings and precautions".

### Stability

Kit component	Opened at 2...8°C	Opened at 18...25°C	Unopened at 2...25°C
Reagent Caps	6 months	6 months	Until the expiration date of the kit
Buffer	Foil pouch: 6 months Single prefilled cuvette: 2 hours	Foil pouch: 3 months	

### Reagent preparation and storage conditions

All reagents are ready to use. Keep the CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

### Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

## 5 Warnings and precautions

### Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink on premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagent inside a reagent cap contains < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Section 15.

### Analytical precautions

- Do not use the product after the expiration date marked on the outer package.
- Do not exceed the stated stability periods for opened reagents.
- The QuikRead go wrCRP+Hb kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (the optical part). Discard any cuvettes with fingerprints.
- The QuikRead go wrCRP+Hb Reagent Caps are colour-coded turquoise to distinguish them from other QuikRead analytes.
- Keep the QuikRead go wrCRP+Hb Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Do not splash liquid into the measuring well of the instrument.

## 6 Sample collection and preparation

### Sample material, collection and volume

CRP results can be obtained from whole blood, plasma and serum samples, whereas Hb results can only be obtained from whole blood samples.

The capillaries included in the kit (Cat. no. 146255) are recommended for adding the sample.

Sample material	Sample volume	Obtainable result	Sample collection
Fingertip sample (whole blood)	10 µl	CRP and Hb	Prick the clean, dry finger with a lancet and discard the first drop. Wipe the finger and collect 10 µl of blood from the second drop in a heparinized capillary.
Anticoagulated whole blood	10 µl	CRP and Hb	Use a venous blood sample collected in a tube containing heparin or EDTA. Mix whole blood by inverting the tube several times and take 10 µl with a capillary.
Plasma	10 µl	CRP	Use EDTA/heparin plasma. Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results.
Serum	10 µl	CRP	Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results.

### Sample dilution

Sample material	Instructions
Whole blood	Do not dilute whole blood samples.
Plasma/serum	Plasma or serum samples may be diluted with a 0.9% NaCl solution before the sample is added to the cuvette. The recommended dilution ratio is 1+3 (by volume, 1 part sample + 3 parts 0.9% NaCl). Take 10 µl of the diluted sample into a cuvette, analyse the sample and multiply the obtained result by 4.

### Sample storage

Sample material	Short-term storage	Long-term storage
Fingertip blood (whole blood)	Max. 5 minutes in the heparinized capillary	Do not store
Anticoagulated whole blood	2...8°C for 3 days	Separate plasma and store as indicated below
Plasma	2...8°C for 7 days	Below -20°C for 3 years <sup>20</sup>
Serum	2...8°C for 7 days	Below -20°C for 3 years <sup>20</sup>
Sample (whole blood, plasma, serum) in buffer	18...25°C for 2 hours	Do not store

Allow the samples to reach room temperature (18...25°C) before testing. Frozen samples must be thawed completely, mixed thoroughly, and brought to room temperature prior to testing. If the thawed plasma samples contain clots, the samples should be centrifuged. Samples should not be repeatedly frozen and thawed.

## 7 Procedure

### Materials required but not provided

The following materials are required to perform the test but are not provided or they are available also separately. The materials provided are listed in Section 4."Reagents".

Material	Cat. no.	Additional info
QuikRead go® Instrument with software version 7.5 or newer	133893, 149915	–
Capillaries (10 µl) 50 pcs	147851	–
Plungers 50 pcs	67966	–
QuikRead go® CRP Control	153764	Recommended for quality control
QuikRead go® CRP Control High	153763	Recommended for quality control
QuikRead go® CRP Control Low	153765	Recommended for quality control
QuikRead go® Hb 10 µl Control	153656	Recommended for quality control
Fingertip lancets	–	–

### Assay procedure

Before starting an assay, read and follow the QuikRead go instrument and QuikRead go wrCRP+Hb test instructions for use.

Open the foil pouch of the cuvette rack, and mark the opening date on the cuvette rack.

### Sample collection (see figures 1–5)

Glass capillaries (10 µl) and plungers are recommended for adding the sample to a prefilled cuvette. The QuikRead go wrCRP+Hb kit (Cat. no.146255) contains capillaries and plungers.

1. A prefilled cuvette must reach room temperature (18...25°C) before use. This takes 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces of the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on the results. The test must be performed within two (2) hours of opening the cuvette. The liquid surface should lie between the two lines marked on the cuvette.
2. Insert the plunger into the capillary at the end with the orange stripe.
3. Fill the capillary with the sample up to the white stopper (10 µl). Ensure that there are no air bubbles in the capillary. Use a soft cloth or tissue to wipe any residual sample from the outside of the capillary.
4. Place the sample capillary into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.
5. Close the cuvette tightly with a wrCRP Reagent Cap. Do not press down the turquoise inner part of the reagent cap. After the sample has been added to the buffer, the solution is stable for two (2) hours. Keep the cuvette upright and do not shake.

### Analysing the sample (see figures 6–8)

Begin the test by turning on the instrument. For full instructions, see the QuikRead go instrument's instructions for use. The display will guide you through the test.

6. Select **Measure** on the display of the QuikRead go instrument.
7. Place the cuvette in the instrument with the barcode facing you, as shown in Figure 7. The instrument detects the sample type. The display shows how the measurement is progressing. The instrument first measures the sample blank and then the CRP and Hb concentrations.

8. The result appears on the display when the measurement is completed and the cuvette automatically rises from the measurement well.

### Analysing the control

Please see the instructions for use provided with each control vial. Analyse the control samples as patient samples but select **Quality Control** on the display of the QuikRead go instrument. The result will be stored as a control measurement.

## 8 Quality control

Regular use of QuikRead go CRP Control (Cat. no. 153764), QuikRead go CRP Control High (Cat. no. 153763) and/or QuikRead go CRP Control Low (Cat. no. 153765) and QuikRead go Hb 10 µl Control (Cat. no. 153656) is recommended. Please consult Section 16. "Troubleshooting" and contact Aidian or your local distributor if you observe any changes in test functionality.

- **QuikRead CRP controls:**

All QuikRead go CRP controls are ready-to-use controls. The control values have been determined for the QuikRead go instrument and the control is applied using the same procedure as for a plasma/serum sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (10 µl). The control is stable for 15 minutes in the cuvette.

- **Other commercially available CRP controls:**

Handle and use the controls in accordance with the instructions for each control substance. Perform the test as outlined in Section 7. "Procedure". Acceptable limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells, because these might not haemolyse normally.

- **QuikRead go Hb control:**

The QuikRead go Hb 10 µl Control is a ready-to-use control. The control values have been determined for the QuikRead go instrument and the control is applied using the same procedure as for a whole blood sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (10 µl).

- **Other commercially available Hb controls:**

Handle and use controls in accordance with the instructions for each control substance. Perform the test as outlined in Section 7. "Procedure". Acceptable limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells because these might not haemolyse normally. Controls in which the Hb is present in a form other than oxyhaemoglobin, or which contain compounds that may modify Hb, might yield results that are not aligned with those obtained by other methods.

## 9 Results' interpretation

### CRP

A rise in the CRP level is non-specific and should not be interpreted without a complete clinical history. Serial CRP measurements may be required to monitor a patient's condition or response to therapy and when measuring CRP levels in newborns. Additionally, establishing the patient-specific baseline CRP level is recommended when assessing slightly elevated CRP levels.

CRP test result	Interpretation of result
< 10 mg/l	<ul style="list-style-type: none"> <li>● Excludes many acute inflammatory diseases but does not specifically exclude inflammatory processes.<sup>21</sup></li> <li>● Possible inflammations in newborn babies are not necessarily excluded.<sup>7-9</sup></li> <li>● Already with elevated initial levels below 10 mg/l the possibility of e.g. a myocardial infarction or recurrent coronary events must be taken into account.<sup>10-12</sup></li> <li>● Possible inflammations or exacerbations in conditions with systemic inflammation are not necessarily excluded.<sup>13-17</sup></li> </ul>
10–50 mg/l	Elevated concentrations are found in acute disease that occur in the presence of a slight to moderate inflammatory process. <sup>21</sup>
> 50 mg/l	Indicates high and extensive inflammatory activity. <sup>21</sup>

If diluted plasma or serum samples have been used, the final results must be calculated in accordance with the instructions in Section 6. "Sample Collection and Preparation". If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

### Hb

Hb test results	Interpretation of result <sup>22</sup>	
Women	< 120 g/l (12.0 g/dl, 7.4 mmol/l)	Below reference value
	> 150 g/l (15.0 g/dl, 9.3 mmol/l)	Above reference value
Men	< 130 g/l (13.0 g/dl, 8.1 mmol/l)	Below reference value
	> 170 g/l (17.0 g/dl, 10.6 mmol/l)	Above reference value
Children 1–12 years*	< 110 g/l (11.0 g/dl, 6.8 mmol/l)	Below reference value
	> 155 g/l (15.5 g/dl, 9.6 mmol/l)	Above reference value

\* Children's Hb values increase progressively reaching gradually adult levels.

For Hb reference values, see Section 11. "Expected values". Hb result levels are affected by patients' diets and the elevation from sea level of their place of residence. It is therefore recommended that each laboratory determines its own reference values for Hb.

When diagnosing patients, the patient's overall clinical status should always be taken into account. If necessary, please consult Section 16. "Troubleshooting" and repeat the test.

## 10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

## 11 Expected values

### CRP

Reference value: <3 mg/l<sup>7,23,24</sup>, median value in apparently healthy adults: 0.8 mg/l<sup>23</sup>. CRP levels in newborns can be affected by e.g. prematurity, duration of delivery and complications related to the delivery. Immediately after birth the CRP levels in term neonates are generally low and will increase to some extent (median below 2 mg/l) over the next few days<sup>7</sup>.

### Hb

Reference group	Reference values <sup>22</sup>
Women	120–150 g/l (12.0–15.0 g/dl, 7.4–9.3 mmol/l)
Men	130–170 g/l (13.0–17.0 g/dl, 8.1–10.6 mmol/l)
Children 1–12 years*	110–155 g/l (11.0–15.5 g/dl, 6.8–9.6 mmol/l)

\* Children's Hb values increase progressively reaching gradually adult levels.

## 12 Performance characteristics

### CRP

#### Precision and reproducibility

A precision study was performed in line with the Clinical and Laboratory Standards Institute (CLSI) guidelines, EP5-A3.

Within-run, between-day and total precision						
Sample material	Sample number	Number of days	Mean CRP (mg/l)	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Whole blood	Sample 1	20	5.8	2.6	2.0	3.4
	Sample 2	20	55	3.8	1.1	4.0
	Sample 3	20	181	5.8	2.8	6.9
Plasma	Sample 1	20	18	4.6	3.4	5.7
	Sample 2	20	63	3.0	1.0	3.5
	Sample 3	20	120	6.1	2.5	6.6
Control 1	Sample 1	20	2.5	2.7	4.2	5.3
Control 2	Sample 1	20	27	2.3	1.5	3.0
Control 3	Sample 1	20	88	5.4	2.9	5.8

### Interference

Interfering substance	No interference found up to concentration
Bilirubin	400 µmol/l
Vitamin C	200 µmol/l
Triglycerides	11.5 mmol/l
Cholesterol	10.0 mmol/l
Rheumatoid factors (RF)	525 IU/ml
Leucocytes	87x10 <sup>9</sup> cells/l
Anticoagulants (Li-heparin or EDTA)	No interference

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test because the assay antibodies lack the Fc part. In rare cases, interference by IgM myeloma protein has been observed.

### Limits of detection

The detection limit of the assay is < 0.4 mg/l.

### Measurement range

For whole blood samples, the measurement range for CRP is 0.5–200 mg/l at the hematocrit level of 40%. If the hematocrit is lower or higher than 40%, the measurement range changes in the manner shown in the table below. The CRP result is not displayed if the level of hematocrit is outside the range 15–75%.

Sample type	Hematocrit %	Sample volume	Measurement range (mg/l CRP)	
Whole blood sample	15–19	10 µl	0.5–140	In cases when the result is outside the measurement range, the result will be displayed as "CRP < 0.5 mg/l" or, for instance, as "CRP > 200 mg/l"
	20–24		0.5–150	
	25–29		0.5–160	
	30–33		0.5–170	
	34–36		0.5–180	
	37–39		0.5–190	
	40–45		0.5–200	
	46–51		0.5–220	
	52–57		0.5–250	
	58–61		0.5–280	
	62–67		0.6–310	
	68–71		0.7–370	
	72–75		0.8–420	
Plasma/serum sample	–	10 µl	0.5–120	

It is possible to turn off the sliding result range, in which case the fixed measurement range for whole blood samples is 0.8–140 mg/l. In order to turn off the sliding measurement range, you have to change the settings of the QuikRead go instrument. The change is made to wrCRP test parameters. Please see the instrument manual.

### Antigen excess

CRP concentrations of less than 1000 mg/l do not give falsely low results.

### Linearity

The recoveries in serial dilutions of CRP plasma samples have ranged from 96% to 111%.

### Hb

#### Precision and reproducibility

A precision study was performed according to the Clinical and Laboratory Standards Institute (CLSI) guidelines, EP5-A3.

Within-run, between-day and total precision						
Sample material	Sample number	Number of days	Mean Hb (g/l)	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Whole blood	Sample 1	20	65	3.1	2.3	3.8
	Sample 2	20	140	2.6	1.7	3.1
	Sample 3	20	196	2.5	2.0	3.6
Control 1	Sample 1	20	128	2.3	0.8	2.4

## Interference

Interfering substance	No interference found up to concentration
Bilirubin	200 µmol/l
Vitamin C	200 µmol/l
Triglycerides	11.5 mmol/l
Rheumatoid factors (RF)	525 IU/ml
B-carotene	3.7 µmol/l
Ibuprofen	500 mg/l
Salicylic acid / Acetylsalicylic acid	650 mg/l
Tetracyclines	200 mg/l
Urea	5000 mg/l
Anticoagulants (Li-heparin or EDTA)	No interference

## Measurement range

The measurement range of haemoglobin in whole blood samples is 50–220 g/l (5.0–22.0 g/dl, 3.1–13.7 mmol/l).

## Linearity

The recoveries in serial dilutions of Hb samples have ranged from 96% to 103%.

## 13 Traceability

### CRP

The calibrators used to calibrate the CRP assay of the QuikRead go wrCRP+Hb test are traceable to the ERM®-DA474 reference material.

### Hb

The Hb assay of the QuikRead go wrCRP+Hb test is traceable to the ICSH (cyanmethaemoglobin) standard 1995 and the CRM BCR-522 reference material.<sup>25,26</sup>

## 14 Method comparison

### CRP

Patient plasma samples were analysed via clinical laboratory method, a point-of-care test and the QuikRead go wrCRP+Hb method. A summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method	Point-of-care test
$y = 0.97x + 0.2$ $r = 0.996$ $n = 97$	$y = 0.99x + 0.7$ $r = 0.995$ $n = 82$

### Whole blood versus plasma comparability

In comparison of 98 clinical samples, whole blood (y) and plasma (x) yielded comparable results:  $y = 1.00x + 0.04$

### Hb

Patient whole blood samples were analysed via clinical laboratory method, a point-of-care test and the QuikRead go wrCRP+Hb method. A summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method	Point-of-care test
$y = 1.03x - 2.1$ $r = 0.996$ $n = 115$	$y = 1.07x - 8.2$ $r = 0.980$ $n = 120$

## 15 Disposal

- Dispose of contents according to national and local law.
- All patient samples, sampling devices, controls, used tubes, caps and cuvettes should be handled and disposed of as potentially infectious material.
- Materials of the components of the QuikRead 101 and QuikRead go kits:
  - Paper: Instructions for use
  - Cardboard: Kit box including its inner parts
  - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, extraction vials and tubes, swabs, plunger and capillary tubes
  - Glass: Capillaries
  - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
- Not to be recycled: Lids of reagent cap tubes (several) and magnetic cards (PVC)
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

## 16 Troubleshooting

### Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the QuikRead go instrument instructions for use.

Error messages	Corrective actions
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner turquoise part of the cap is not pressed down.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25°C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18...25°C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. If this error message appears again, please make sure that the CRP reagent caps have not been exposed to humidity. Also ensure that the reagent cap tube is properly closed.
No result. Hematocrit too low.	If the level of hematocrit in a whole blood sample is below the hematocrit measurement range of the instrument (see Section 12), no result is shown.
No result. Hematocrit too high.	If the level of hematocrit in whole blood samples is above the hematocrit measurement range of the instrument (see Section 12), no result is shown.

### Unexpected low and high results




Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected low result.	Too small sample volume.	Perform a new test. Make sure that the capillary is completely filled. Avoid air bubbles in the capillary.
Unexpected high result.	Too large sample volume.	Perform a new test. Make sure that the capillary is properly filled. Wipe away any excess sample from the outside of the capillary.
	Too low buffer volume.	Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
	The cuvette is dirty.	Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette.
Unexpected low/high result.	Use of components from different kit lots or tests	Perform a new test. Make sure that all the reagents are from the same kit lot.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored in the manner described in the instructions for use.

- Price CP et al. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med.* 1999; 37 (2): 109–13.
- van Leeuwen MA & van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillieres Clin Rheumatol.* 1994; 8 (3): 531–52.
- Olaison L et al. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med.* 1997; 157 (8): 885–92.
- Peltola H et al. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics.* 1997; 99 (6): 846–50.
- Philip AGS & Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics.* 2000; 106 (1): 1–5.
- Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, (eds). *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003; 2: 150–6.
- Chiesa C et al. C reactive protein and procalcitonin: Reference intervals for preterm and term newborns during the early neonatal period. *Clin Chim Acta.* 2011; 412 (11-12): 1053–9.
- Hofer N et al. Non-infectious conditions and gestational age influence C-reactive protein values in newborns during the first 3 days of life. *Clin Chem Lab Med.* 2011; 49 (2): 297–302.
- Wasunna A et al. C-reactive protein and bacterial infection in preterm infants. *Eur J Pediatr.* 1990; 149 (6): 424–7.
- Haverkate F et al. Production of C-reactive protein and risk of coronary events in stable and unstable angina. European Concerted Action on Thrombosis and Disabilities Angina Pectoris Study Group. *Lancet.* 1997; 349 (9050): 462–6.
- Haverkate F et al. C-reactive protein and cardiovascular disease. *Fibrinolysis & Proteolysis.* 1997; 11 (Suppl.1): 133–4.
- Liuzzo G et al. The Prognostic value of C-reactive protein and serum amyloid A protein in severe unstable angina. *N Engl J Med.* 1994; 331 (7): 417–24.
- Pinto-Plata VM et al. C-reactive protein in patients with COPD, control smokers and non-smokers. *Thorax.* 2006; 61 (1): 23–8.
- Mannino DM et al. Obstructive and restrictive lung disease and markers of inflammation: data from the Third National Health and Nutrition Examination. *Am J Med.* 2003; 114 (9): 758–62.
- Dahl M et al. C-reactive protein as a predictor of prognosis in chronic obstructive pulmonary disease. *Am J Respir Crit Care Med.* 2007; 175 (3): 250–5.
- Dessein PH et al. High sensitivity C-reactive protein as a disease activity marker in rheumatoid arthritis. *J Rheumatol.* 2004; 31 (6): 1095–7.
- Poole CD et al. The association between C-reactive protein and the likelihood of progression to joint replacement in people with rheumatoid arthritis: a retrospective observational study. *BMC Musculoskelet Disord.* 2008; 9: 146. /..

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
18. Pearson TA et al. Markers of inflammation and cardiovascular disease: application to clinical and public health practice: A statement for healthcare professionals from the Centers for Disease Control and Prevention and the American Heart Association. *Circulation*. 2003; 107 (3): 499-511
19. Greenland P et al. ACCF/AHA guideline for assessment of cardiovascular risk in asymptomatic adults: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. *J Am Coll Cardiol*. 2010; 56 (25): e50-103.
20. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples (WHO/DIL/LAB/99.1 Rev.2), 2002. Available: [http://whqlibdoc.who.int/hq/2002/WHO\\_DIL\\_LAB\\_99.1\\_Rev.2.pdf](http://whqlibdoc.who.int/hq/2002/WHO_DIL_LAB_99.1_Rev.2.pdf). Accessed 21.8.2014.
21. Thomas L. C-reaktives Protein (CRP) In: Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780–90.
22. Bain BJ et al. Dacie and Lewis Practical Haematology, 11th ed. Edinburgh, Churchill Livingstone, 2012.
23. Shine et al. Solid phase immunoassays for human C-reactive protein. *Clin Chim Acta*. 1981; 117: 13-23.
24. Schlebush H et al. High Sensitive CRP and Creatinine: Reference intervals from Infancy to Childhood. *J Lab Med* 2002; 26 (5/6): 341-6.
25. International Council for Standardisation in Haematology: Expert Panel on Haemoglobinometry. Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1995) and specifications for international haemoglobinocyanide standard (4th edition). *J Clin Pathol* 1996; 49: 271–4.
26. Bull, BS et al. Reference and selected procedures for the quantitative determination of hemoglobin in blood; approved standard – third edition", NCCLS 20 (28), document H15-A3.

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