QuikRead go[®] Hb 10 µl Control, Cat. no. 153656	<i>English</i>
------------------------------------------------------------------	----------------

Intended use

QuikRead go[®] Hb 10 µl Control is intended for quality control of the Hb assays of the QuikRead go[®] wrCRP+Hb tests by the QuikRead go[®] instrument.
For *in vitro* diagnostic use.

Principles of the procedure

Use the control in the same way as unknown **whole blood** specimen would be used, following the instructions for the assay procedure.

Reagents

The control is from bovine origin. It should be handled in the same way as patient sample.

Preparation

The control is ready-to-use. Allow the material to stabilize for 30 minutes at room temperature (18...25°C).
Mix well. Do not shake.

Storage

Stable until the expiry date when stored at 2...8°C.
Opened vial: Stable for 1 month at 2...8°C.
Do not use any vial that has visible evidence of microbial growth.

Expected values

The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.

QuikRead go[®] Hb 10 µl Control, Kat. -Nr. 153656	<i>Deutsch</i>
-------------------------------------------------------------------	----------------

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Warnings and precautions

• Avoid using the last drops of the control because opening of the vial many times may result in evaporation, which may bias the results. If the control is not in the acceptable range, take a new vial.

• Do not use the product after the expiry date marked on the outer package.

• Do not exceed the stability periods for opened reagents.
• QuikRead go Hb 10 µl Control is from bovine origin. As a precaution, handle QuikRead go Hb 10 µl Control as if capable of transmitting infectious agents.

Disposal

• Dispose of contents according to national and local law.
• Materials of the components
Paper: Instructions for use
Cardboard: Outer box
Glass: Control vials
Plastic: Vial caps
• When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

QuikRead go[®] Hb 10 µl Control, Kat. -Nr. 153656	<i>Deutsch</i>
-------------------------------------------------------------------	----------------

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