AIDIAN

UniQ[®] PIIINP RIA

Intact N-terminal propeptide of type III
procollagen, radioimmunoassay

Cat. No. 68570

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Intact N-terminal propeptide of type III collagen, radioimmunoassay

1. INTENDED USE

UniQ PIIINP RIA is a quantitative radioimmunoassay designed for measurement of the concentration of intact aminoterminal propeptide of type III procollagen in human serum. **For *in vitro* diagnostic use.**

2. SUMMARY AND EXPLANATION

Type III collagen is derived from a larger protein – type III procollagen – which has propeptide extensions at both ends of the molecule. Some of the aminoterminal propeptides are set free during the synthesis and deposition of type III collagen and some are retained in the molecules which remain on the surface of the collagen fibrils. This antigen, when found in serum, can thus be derived from the synthesis (including late processing) of new type III collagen or from the degradation of existing type III collagen fibrils.

The further metabolism of PIIINP is not known in detail. Smaller degradation products related to PIIINP are found in serum and excreted in urine. Antigens with a higher molecular weight than the propeptide are also found in serum. Their exact molecular structure is unknown.

Increased concentrations of PIIINP are found in a number of conditions where accumulation and/or degradation of connective tissue takes place. These include fibroproliferative, haematological, endocrinological and malignant diseases. Thus the changes in PIIINP are not specific for a particular disease but reflect the involvement and altered metabolism of type III collagen.

In growth deficient children, the response to growth hormone therapy is reflected in an early increase in serum PIIINP thus providing a means to identify those who respond to the therapy.

Liver fibrosis and cirrhosis of various aetiologies increase serum concentrations of PIIINP. PIIINP provides a non-invasive means of monitoring methotrexate treated psoriatics, thereby reducing the need for liver biopsies. In chronic hepatitis C patients with advanced fibrosis, PIIINP is a useful tool for the follow-up of changes in fibrosis during and after antiviral therapy, and it can provide prognostic information in these patients. In fibrotic conditions of lungs, such as idiopathic lung fibrosis or pulmonary sarcoidosis, serum PIIINP levels may be normal, but a high concentration can be found in bronchoalveolar lavage fluid.

Serum PIIINP is increased in myeloproliferative diseases, particularly during the active phases of myelofibrosis. In polycythaemia vera, a high PIIINP value may predict transformation into myeloid metaplasia.

Increased serum PIIINP levels are found in several malignant diseases, including cancers involving the liver. In breast carcinoma with bone metastases, changes in PIIINP can serve as an indicator of response to therapy. In ovarian carcinoma, changes in serum PIIINP correlate with the development of the disease.

The concentration of PIIINP in serum reflects the healing process after myocardial infarction and predicts left ventricular functional and volume changes. It has also been shown to correlate with scar formation, left ventricular dysfunction and survival rate. The increase in serum PIIINP in the four days following infarction has been shown to be higher in patients with poor outcomes. Furthermore, PIIINP has value as a prognostic predictor of death and hospitalisation.

3. PRINCIPLES OF THE TEST

UniQ PIIINP RIA kit is based on the competitive radioimmunoassay technique. A known amount of labelled PIIINP and an unknown amount of unlabelled PIIINP in the sample compete for a limited number of high-affinity binding sites on the antibody. After separating the free antigen, the amount of labelled PIIINP in the sample tube is inversely proportional to the amount of PIIINP in the sample. The concentrations in unknown samples are obtained from a calibration curve.

4. REAGENTS

UniQ® PIIINP RIA, Cat. No. 68570

Component name and origin	Symbol	∇ 100	Stability at 2...8°C	
			Unopened	Opened
Tracer	Ag ¹²⁵ I	22 ml Radioactivity < 200 kBq	Until the expiration date marked on the label	5 weeks
Antiserum ORIG RAB	Ab	22 ml	Until the expiration date marked on the label	5 weeks
Calibrator ORIG HUM	CAL xxµg/l	7 x 1.5 ml	Until the expiration date marked on the label	5 weeks
Control ORIG HUM	CONTROL	1.5 ml	Until the expiration date marked on the label	5 weeks after reconstitution
Procollagen Separation Reagent	REAG SEPAR	60 ml	Until the expiration date marked on the label	5 weeks
Instructions for use				

Reagents contain sodium azide as a preservative. See section 5 “Warnings and precautions”.

Reagent preparation

The control is lyophilized. All other reagents are ready to use.

Allow the **control** vial to reach room temperature (18...25°C) before opening. Reconstitute the control by adding 1.5 ml distilled or deionised water to the vial. Cap and mix well by gentle swirling or inverting to avoid foaming. Allow to stand for 30 minutes before use.

Reagent deterioration

Alterations in the physical appearance of the reagents or in the slope of the calibrator curve, or values of control sera outside the manufacturer’s acceptable range may be an indication of reagent instability. Do not use reagents if you suspect reagent deterioration.

5. WARNINGS AND PRECAUTIONS

Health and safety information

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- All patient samples and controls should be handled as potentially infectious material.
- The controls contain material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle the controls as if capable of transmitting infectious agents.
- Lyophilised reagents contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Contact with acids liberates very toxic gas (EUH032). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Azides may react with metal plumping, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Reconstituted and liquid reagents contain <0.1% sodium azide, which is not considered a harmful concentration.
- The kit contains ¹²⁵I (half-life $t_{1/2}$ of 60 days) emitting 35.5 keV gamma rays and 27–32 keV X-rays. Radioactive material contained in this kit may be received, acquired, possessed and used only by physicians, clinical laboratories, research institutions or hospitals and used only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom to human beings or animals. All radioactive materials should be stored and handled only in specifically designated areas. The purchase, storage, use or exchange of radioactive products is subject to the laws in force in the user's country.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water and decontaminate, following your institution's guidelines for accidental spills.
- Spills of radioactive materials should be cleaned immediately with a detergent and the contaminated cleaning materials transferred to the radioactive waste container or handled according to the recommendations of the radiation safety office.
- Refer to the radiation regulations and procedures applicable to your institution for handling, storage and disposal of radioactive materials.

Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Do not mix components from different lot numbers or different tests.

6. SAMPLE COLLECTION AND PREPARATION

No special preparation of the patient is required prior to sample taking.

Serum should be used and a normal venous blood sample should be collected using an acceptable medical technique, following the normal precautions for venipuncture. The assay is not to be used with plasma samples. Highly lipemic or hemolyzed patient samples should be avoided.

Sample handling and storage

After the specimen collection, allow the blood to clot and separate serum by centrifugation. Serum samples may be stored for up to 5 days at 2...8°C. For longer periods, store the samples at -20°C or colder. Repeated freezing and thawing should be avoided.

Samples with high PIINP concentrations may be diluted using 0.9 % saline.

7. PROCEDURE

See the section 4. "Reagents", for a list of reagents.

Materials required but not provided

- Pipettes (200 µl, 500 µl and 1.5 ml)
- Test tubes (polypropylene or polystyrene)
- Absorbent paper
- Paraffin film
- Vortex Mixer
- Water bath (capable of maintaining 37 ±1°C)
- Decanting rack
- Centrifuge (capable of at least 2000 g)
- Gamma counter
- Distilled water

Instruments

Any gamma counter capable of measuring ¹²⁵I may be used. The counting efficiency of the instrument should be regularly checked as well as the background counts. It is recommended that the instrument be calibrated systematically. Consult the operations manual supplied by the instrument's manufacturer for details of installation, use, maintenance and calibration.

Details of the procedure

1	Bring all reagents and samples to room temperature (18...25°C) at least 30 minutes before use.
2	Label test tubes in duplicate for non-specific binding (NSB), calibrators, control, patient samples and totals.
3	Pipette 200 µl of calibrator, control or patient sample into the appropriate tubes. Any patient sample may be pipetted into the NSB tubes.
4	Pipette 200 µl of tracer (red) into all tubes.
5	Pipette 200 µl of antiserum (blue) into all tubes except NSB and total. Pipette 200 µl of distilled water into the NSB tubes.
6	Mix all tubes using a vortex mixer, cover tubes with paraffin film and incubate for 2 hours at 37 ±1°C.
7	Mix the separation reagent thoroughly by gentle inversion and add 500 µl to all tubes except totals.
8	Mix the tubes using a vortex mixer and incubate for 30 minutes at room temperature (18...25°C).
9	Centrifuge all tubes except totals for 15 minutes at 2000 g (minimum) at 4...20°C (4°C is recommended).
10	Decant the supernatants from all tubes except the totals by simultaneous inversion. Gently blot the tubes against absorbent material and allow to drain for a few seconds to remove any remaining liquid. Care should be taken not to decant any of the precipitate. Do not reinvert the tubes once they have been turned upright. Failure to decant and blot tubes may result in poor replication and erroneous values.
11	Count each tube using a gamma counter for at least 1 minute or until 10,000 counts per tube have accumulated.
12	Calculate the results as described in the "Results' interpretation" section.

Assay procedure – summary

(all values given in µl)

	Total	NSB	Calibrator	Control & unknown
Pipette sample		200	200	200
Pipette tracer (red)	200	200	200	200
Pipette antiserum (blue)			200	200
Pipette distilled water		200		
Mix		x	x	x
Incubate for 2 hours (37°C)		x	x	x
Pipette separation reagent		500	500	500
Mix		x	x	x
Incubate for 30 min at 18...25°C		x	x	x
Centrifuge for 15 min at 2000 g (4...20°C)		x	x	x
Decant		x	x	x
Count 1 min. or 10 000 counts	x	x	x	x

8. QUALITY CONTROL

An internal quality assurance programme is recommended for all clinical laboratories. For this reason, it is suggested that control sera in both the lower and upper portions of the calibration curve be analysed to enable the performance of the procedure to be monitored. The UniQ PIIIINP RIA assay provides one control material: PIIIINP RIA Control. The expected range of the control is indicated on a separate sheet provided with the assay. It is recommended that each laboratory establish its own acceptable ranges for the control it uses. The mean of the control should fall within the acceptable range stated by Aidian. It is advised to use at least two separate kit lot numbers for the control value assignment. Once the control ranges have been established, controls should be run in duplicate with each assay.

9. RESULTS' INTERPRETATION

For automatic result processing spline function curve fitting is recommended. Results can also be produced manually on semi-log graph paper.

- **Calculate** the mean count for all calibrators, samples and controls.
- **Calculate** the %B/B₀ from:

$$\%B/B_0 = \frac{\text{calibrator or sample count} - \text{NSB}}{\text{calibrator 0 count} - \text{NSB}} \times 100$$

where B is the bound radioactivity of a sample or calibrator, B₀ is the bound radioactivity of the calibrator 0 and NSB is the radioactivity of the non-specific binding.

- **Draw** a calibration curve on semi-log graph paper with %B/B₀ values on the ordinate (Y axis) and PIIIINP concentrations (µg/l) of the calibrators on the abscissa (X axis).
- **Determine** the PIIIINP concentrations of the unknowns from the calibration curve.
- **Any sample reading lower than the lowest calibrator should be reported as such.**
- **Any sample reading higher than the highest calibrator should be diluted and re-assayed.** After dilution, values must be corrected for the dilution factor.

Fig. 1 Typical calibration curve
This calibration curve is for illustration only.


Table 1 Calculation of results using typical data.
This data is for illustration only and must not be used for the calculation of any sample result.

Tube	mean cpm	cpm-NSB	%B/B ₀	PIIINP (µg/l)
Total	48807			
NSB	679			
Cal 0	30976	30297		
Cal 1.0	26705	26026	85.9	
Cal 2.5	22047	21368	70.5	
Cal 5.0	17163	16484	54.4	
Cal 10	12100	11421	37.7	
Cal 25	7241	6562	21.7	
Cal 50	4795	4116	13.6	
Unkn 1	17342	16663	55.0	4.9
Unkn 2	8468	7789	25.7	19.0

10. LIMITATIONS OF PROCEDURE

As a basis for diagnostic and therapeutic decisions the results of any single diagnostic test must be backed up with other clinically relevant data.

11. EXPECTED VALUES

Due to ethnic, dietary and age variations, the reference limits given may not apply to all populations. Therefore each laboratory should establish its own representative reference limits.

PIIINP values of apparently healthy adults (19–65 years) were measured with the UniQ PIIINP RIA kit. The non-parametric method recommended by IFCC was used to calculate the reference intervals for the upper and lower limits.

Table 2 Serum PIIINP reference limits

	Adults
Number of subjects	232
Mean	3.7 µg/l
Reference interval*	2.3–6.4 µg/l
0.90-confidence intervals for lower reference limit	1.9–2.6 µg/l
for upper reference limit	5.5–6.7 µg/l

*Reference interval = 0.025 & 0.975 fractiles

Reference values for children were derived from a population of healthy children aged 0 months to 19 years (modified from Trivedi P et al. J Pediatr 1989; 114:225-230).

Table 3 Serum PIIINP values for children.

Age group	n (F)	Female ($\mu\text{g/l}$)	n (M)	Male ($\mu\text{g/l}$)
1–3 months	4	53.22 \pm 10.87	5	62.05 \pm 15.59
3–6 months	5	31.25 \pm 6.53	13	33.94 \pm 9.35
6–9 months	7	19.71 \pm 3.53	4	18.14 \pm 3.26
9–12 months	3	19.88 \pm 6.16	4	15.25 \pm 4.65
1–2 years	9	13.45 \pm 2.25	7	11.71 \pm 4.30
2–3 years	10	9.85 \pm 2.65	9	8.47 \pm 2.53
3–10 years *)	63	7.75 \pm 2.17	79	7.75 \pm 2.18
10–11 years	7	8.02 \pm 2.56	12	7.78 \pm 1.61
11–12 years	12	11.23 \pm 2.99	9	7.27 \pm 2.15
12–13 years	5	10.85 \pm 3.77	6	8.96 \pm 2.94
13–14 years	2	5.84 \pm 2.86	4	13.25 \pm 5.58
14–15 years	5	8.09 \pm 1.57	2	12.92 \pm 4.10
15–16 years	2	7.34 \pm 0.11	7	15.42 \pm 4.28
16–17 years	3	3.98 \pm 0.99	–	–
17–18 years	4	4.07 \pm 1.08	5	4.88 \pm 1.06

*) no significant change between 3–10 years.

12. PERFORMANCE CHARACTERISTICS

The results presented below were obtained at Aidian and represent typical performance data.

Measurement range

1.0–50 $\mu\text{g/l}$.

Technical detection limit

Approximately 0.4 $\mu\text{g/l}$, defined as twice the standard deviation of the 0-binding value.

Linearity

Linearity was studied by diluting three samples with serum and three samples with 0.9 % NaCl. Mean recovery of the samples was 102 % for 0.9 % NaCl and 100 % for serum.

Recovery

Known amounts of PIIINP were added to four serum sample pools containing 2.9–4.5 $\mu\text{g/l}$ of PIIINP. Recoveries were in the range of 83–100 % with a mean value of 93 %.

Precision and reproducibility

A precision study was performed in line with the Clinical and Laboratory Standards Institute (CLSI) guidelines, EP5-A3.

Table 4 Repeatability, between-day and within laboratory precision (20 days)

Sample	Mean PIIINP (µg/l)	N	Repeatability		Between day		Within laboratory (Total)	
			SD (µg/l)	CV% (≤ 10%)	SD (µg/l)	CV% (≤ 10%)	SD (µg/l)	CV% (≤ 10%)
Serum pool 1	2.81	80	0.11	4.1	0.13	4.7	0.18	6.3
Serum pool 2	6.38	80	0.32	4.9	0.31	4.9	0.44	6.9
Serum pool 3	6.64	80	0.29	4.4	0.22	3.3	0.37	5.5
Serum pool 4	16.4	80	1.44	8.7	0.46	2.8	1.53	9.3

Interfering substances

Serum bilirubin concentration < 400 µmol/l does not interfere. Serum haemoglobin concentration up to 5 g/l does not interfere. Triglycerides up to 30 g/l do not interfere.

Specificity

The intact UniQ PIIINP RIA kit is not sensitive to the smaller degradation products found in blood. It measures propeptide and its higher molecular weight forms. UniQ PIIINP RIA does not cross-react with antigens used in the UniQ PINP RIA kit.

13. TRACEABILITY

The measurand (analyte) in the UniQ PIIINP RIA Calibrators is traceable to the manufacturer's working calibrators, which are prepared from highly purified PIIINP whose value is assigned using amino acid composition analysis.

14. DISPOSAL

- Dispose of contents according to national and local law.
- All patient samples and controls should be handled and disposed of as potentially infectious material.
- Radioactive waste material must be discarded according to the applicable laws.
- Materials of the components:
 - Paper: Instructions for use
 - Cardboard: Kit box
 - Plastic: Separation reagent bottle, vial caps
 - Rubber: Stoppers
 - Glass: Tracer, antiserum, calibrator, control vials
 - Foam rubber
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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LITERATURE

- Alla F et al.** Early changes in serum markers of cardiac extra-cellular matrix turnover in patients with uncomplicated hypertension and type II diabetes. *Eur J Heart Fail* 2006;8:147-153.
- Chalmers RJG et al.** Replacement of routine liver biopsy by procollagen III aminopeptide for monitoring patients with psoriasis receiving long-term methotrexate: a multicentre audit and health economic analysis. *Br J Dermatol* 2005;152:444-450.
- Fontana RJ et al.** Serum fibrosis marker levels decrease after successful antiviral treatment in chronic hepatitis C patients with advanced fibrosis. *Clin Gastroenterol Hepatol* 2009;7:219-226.
- Host NB et al.** The aminoterminal propeptide of type III procollagen provides new information on prognosis after acute myocardial infarction. *Am J Cardiol* 1995;76:869-873.
- Lammi L et al.** Type III and type I procollagen markers in fibrosing alveolitis. *Am J Respir Crit Care Med* 1999;159:818-823.
- Longobardi S et al.** Growth hormone (GH) effects on bone and collagen turnover in healthy adults and its potential as a marker of GH abuse in sports: A double blind, placebo-controlled study. *J Clin Endocrinol Metab* 2000;85:1505-1512.
- Martos R et al.** Diagnosis of heart failure with preserved ejection fraction: improved accuracy with the use of markers of collagen turnover. *Eur J Heart Fail* 2009;11:191-197.
- Maurice PDL et al.** Monitoring patients on methotrexate: hepatic fibrosis not seen in patients with normal serum assays of aminoterminal peptide of type III procollagen. *Br J Dermatol* 2005;152:451-458.
- Poulsen S et al.** Relationship between serum amino-terminal propeptide of type III procollagen and changes of left ventricular function after acute myocardial infarction. *Circulation* 2000;101:1527-1532.
- Risteli J et al.** Rapid equilibrium radioimmunoassay for the aminoterminal propeptide of human type III procollagen. *Clin Chem* 1988;34:715-718.
- Risteli L and Risteli J.** Analysing connective tissue metabolites in human serum. Biochemical, physiological and methodological aspects. *J Hepatol* 1995;22:77-81.
- Sartorio A et al.** Long-term monitoring of rec-GH treatment by serial determination of serum aminoterminal propeptide of type III procollagen in children and adults with GH deficiency. *J Endocrinol Invest* 1999;22:169-175.
- Satta J et al.** Increased turnover of collagen in abdominal aortic-aneurysms, demonstrated by measuring the concentration of the aminoterminal propeptide of type III procollagen in peripheral and aortic blood-samples. *J Vascular Surgery* 1995;22:155-160.
- Soylemezoglu O et al.** Urinary and serum type III collagen: markers of renal fibrosis. *Nephrol Dial Transplant* 1997;12:1883-1889.
- Tapanainen P et al.** Collagen metabolites in the prediction of response to GH therapy in short children. *Eur J Endocrinol* 1997;137:621-625.
- Trivedi P et al.** Growth velocity, growth hormone therapy, and serum concentrations of the amino-terminal propeptide of type III procollagen. *J Pediatr* 1989;114:225-230.
- Uusimaa P et al.** Collagen scar formation after acute myocardial infarction: Relationships to infarct size, left ventricular function, and coronary artery patency. *Circulation* 1997;96:2565-2572.
- Vihervuori E et al.** Collagen formation and degradation increase during growth hormone therapy in children. *Bone* 1997;20:133-138.
- Zachariae H et al.** The value of amino-terminal propeptide of type III procollagen in routine screening for methotrexate-induced liver fibrosis: a 10-year follow-up. *Br J Dermatol* 2001;143:100-103.
- Zannad F et al.** Limitation of excessive extracellular matrix turnover may contribute to survival benefit of spironolactone therapy in patients with congestive heart failure. *Circulation* 2000;102:2700-2706.

UniQ[®] is a registered trademark of Aidian Oy.

AIDIAN

Aidian Oy
Koivu-Mankkaan tie 6 B, FI-02200 Espoo, Finland
P.O. Box 83, FI-02101 Espoo, Finland
+358 10 309 3000, aidian.eu