**AIDIAN**

# UniQ<sup>®</sup> ICTP RIA





C-terminal telopeptide of type 1 collagen,  
radioimmunoassay

Cat. No. 68601

132033-20

**Explanation of symbols • Erläuterung der Symbole • Explication des symboles  
Verklaring van symbolen • Explicación de los símbolos**

	English	German	Français	Nederlands	Español
	<i>In vitro</i> diagnostic medical device	<i>In-vitro</i> Diagnostikum	Dispositif médical de diagnostic <i>in vitro</i>	Alleen voor <i>in vitro</i> diagnostisch gebruik	Producto sanitario para diagnóstico <i>in vitro</i>
	Catalogue number	Artikelnummer	Référence	Catalogus nummer	Número de catálogo
	Batch code	Chargenbezeichnung	Lot	Batch code	Número de lote
	Use by	Verwendbar bis	Utiliser jusqu'au	Gebruik door	Fecha de caducidad
	Temperature limit	Temperaturbegrenzung	Limites de température	Temperatuurlimiet	Limitación de temperatura
	Consult instructions for use	Nachschlagen in der Bedienungsanleitung	Consulter la notice d'utilisation	Raadpleeg gebruiksinstructie	Consulte las instrucciones de uso
	Manufacturer	Hersteller	Fabricant	Producent	Fabricante
	Sufficient for	Ausreichend für	Suffisant pour	Voldoende voor	Válido para
	Tracer	Tracer	Traceur	Tracer	Trazador
	Antiserum	Antiserum	Antiserum	Antiserum	Antisuero
	Separation Reagent	Trennreagenz	Réactif de séparation	Scheidingsreagens	Reactivo de separación
	Control	Kontrolle	Contrôle	Controle	Control
	Calibrator concentration	Kalibrator-konzentration	Concentration de l'étalon	Kalibrator concentratie	Concentración del calibrador
	Origin: rabbit	Ursprung: Kaninchen	Origine: lapin	Oorsprong: konijn	Origen: conejo
	Origin: human	Ursprung: menschlich	Origine: humaine	Oorsprong: humane	Origen: humano
	Contains sodim azide	Enthält Natriumazid	Contient de l'azide de sodim	Bevat natrium azide	Contiene azida sódica

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### C-terminal telopeptide of type I collagen, radioimmunoassay

#### 1. INTENDED USE

UniQ ICTP RIA is a quantitative radioimmunoassay designed for measurement of the concentration of carboxyterminal cross-linked telopeptide of type I collagen in human serum. **For *in vitro* diagnostic use.**

#### 2. SUMMARY AND EXPLANATION

Type I collagen is the most abundant collagen type in the body and the only collagen type found in mineralised bone, where it accounts for more than 90 % of the organic matrix. In addition, type I collagen is found in loose connective tissues together with other collagen types such as types III, V and VI. However, type I collagen is the most prevalent type in these tissues.

Type I collagen in tissues is present in fibres, the structure of which shows some variation according to tissue. In bone, type I collagen molecules are cross-linked via three residues of hydroxylysine, lysine or their derivatives which form a fluorescent, cyclic pyridoline-structure and non-fluorescent unknown structures linking three different collagen polypeptide chains together. In loose connective tissues, such as skin, the major mature cross-link of type I collagen is non-fluorescent and contains histidine as one of the amino acid residues.

ICTP is the carboxyterminal telopeptide region of type I collagen, joined via trivalent cross-links and liberated during the degradation of mature type I collagen. This peptide is found in an immunochemically intact form in blood, where it seems to be derived from bone resorption and degradation of loose connective tissues. It has been shown that the ICTP antigen is produced through the action of matrix metalloproteinases, which are enzymes involved in tissue destruction in various pathological conditions.

Increased serum concentrations of ICTP are hence seen in conditions associated with increased lysis of bone, such as multiple myeloma, osteolytic metastases, rheumatoid arthritis and immobilisation. As ICTP is not produced through physiological Cathepsin-K-mediated bone resorption, its concentration is less affected by conditions such as menopause.

#### 3. PRINCIPLES OF THE TEST

UniQ ICTP RIA kit is based on the competitive radioimmunoassay technique. A known amount of labelled ICTP and an unknown amount of unlabelled ICTP in the sample compete for a limited number of high-affinity binding sites on the antibody. After separating the free antigen, the amount of labelled ICTP in the sample tube is inversely proportional to the amount of ICTP in the sample. The concentrations in unknown samples are obtained from a calibration curve.

## 4. REAGENTS

UniQ® ICTP RIA, Cat. No. 68601

Component name and origin	Symbol	 100	Stability at 2...8°C	
			Unopened	Opened
Tracer	<b>Ag</b> <sup>125</sup> I	22 ml Radioactivity < 200 kBq	Until the expiration date marked on the label	6 weeks
Antiserum <b>ORIG</b> <b>RAB</b>	<b>Ab</b>	22 ml	Until the expiration date marked on the label	6 weeks
Calibrator	<b>CAL</b> xxµg/l	7 x 0.75 ml	Until the expiration date marked on the label	6 weeks
Control <b>ORIG</b> <b>HUM</b>	<b>CONTROL</b> #	2 x 1.0 ml	Until the expiration date marked on the label	1 week after reconstitution.
ICTP RIA Separation Reagent	<b>REAG</b> <b>SEPAR</b>	60 ml	Until the expiration date marked on the label	6 weeks
Instructions for use				

Reagents contain sodium azide as a preservative. See section 5 "Warnings and precautions".

### Reagent preparation

The controls are lyophilized. All other reagents are ready to use.

Allow the **control** vials to reach room temperature (18...25°C) before opening. Reconstitute the controls by adding 1.0 ml distilled or deionised water to each vial. Cap and mix well by gentle swirling or inverting to avoid foaming. Allow to stand for 30 minutes before use.

### Reagent deterioration

Alterations in the physical appearance of the reagents or in the slope of the calibrator curve, or values of control sera outside the manufacturer's acceptable range may be an indication of reagent instability. Do not use reagents if you suspect reagent deterioration.

## 5. WARNINGS AND PRECAUTIONS

### Health and safety information

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- All patient samples and controls should be handled as potentially infectious material.
- The controls contain material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a

precaution, handle the controls as if capable of transmitting infectious agents.

- Lyophilised reagents contain <1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501).
- Reconstituted and liquid reagents contain <0.1% sodium azide, which is not considered a harmful concentration.
- Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- The kit contains  $^{125}\text{I}$  (half-life  $t_{1/2}$  of 60 days) emitting 35.5 keV gamma rays and 27–32 keV X-rays. Radioactive material contained in this kit may be received, acquired, possessed and used only by physicians, clinical laboratories, research institutions or hospitals and used only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom to human beings or animals. All radioactive materials should be stored and handled only in specifically designated areas. The purchase, storage, use or exchange of radioactive products is subject to the laws in force in the user's country.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water and decontaminate, following your institution's guidelines for accidental spills.
- Spills of radioactive materials should be cleaned immediately with a detergent and the contaminated cleaning materials transferred to the radioactive waste container or handled according to the recommendations of the radiation safety office.
- Refer to the radiation regulations and procedures applicable to your institution for handling, storage and disposal of radioactive materials.

### **Analytical precautions**

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Do not mix components from different lot numbers or different tests.

## **6. SAMPLE COLLECTION AND PREPARATION**

No special preparation of the patient is required prior to sample taking.

Serum should be used and a normal venous blood sample should be collected using an acceptable medical technique, following the normal precautions for venipuncture. The assay is not to be used with plasma samples.

Highly lipemic, hemolyzed or microbiologically contaminated patient samples should be avoided.

### **Sample handling and storage**

After the specimen collection, allow the blood to clot and separate serum by centrifugation. Serum samples may be stored for up to 5 days at 2...8°C. For longer periods, store the samples at -20°C or colder. Repeated freezing and thawing should be avoided.

Samples with high ICTP concentrations may be diluted using 0.9 % saline.

## 7. PROCEDURE

See the section 4. "Reagents", for a list of reagents.

### Materials required but not provided

- Pipettes (100  $\mu$ l, 200  $\mu$ l, 500  $\mu$ l and 1000  $\mu$ l)
- Test tubes (polystyrene)
- Absorbent paper
- Paraffin film
- Vortex Mixer
- Water bath (capable of maintaining  $37 \pm 1^\circ\text{C}$ )
- Decanting rack
- Centrifuge (capable of at least 2000 g)
- Gamma counter
- Distilled water

### Instruments

Any gamma counter capable of measuring  $^{125}\text{I}$  may be used. The counting efficiency of the instrument should be regularly checked as well as the background counts. It is recommended that the instrument be calibrated systematically. Consult the operations manual supplied by the instrument manufacturer for details of installation, use, maintenance and calibration.

### Details of the procedure

1	Bring all reagents, controls and samples to room temperature ( $18...25^\circ\text{C}$ ) at least 30 minutes before use.
2	Label test tubes in duplicate for NSB (non-specific binding), calibrators, controls, patient samples and totals.
3	Pipette 100 $\mu$ l of calibrator, control or patient sample into the appropriate tubes. Any patient sample may be pipetted into the NSB tubes.
4	Pipette 200 $\mu$ l of tracer (red) into all tubes.
5	Pipette 200 $\mu$ l of antiserum (blue) into all tubes except NSB and totals. Pipette 200 $\mu$ l of distilled water into the NSB tubes.
6	Mix all tubes using a vortex mixer, cover with paraffin film and incubate for 2 hours at $37 \pm 1^\circ\text{C}$ .
7	Mix the separation reagent thoroughly by gentle inversion and add 500 $\mu$ l to all tubes except totals.
8	Mix the tubes using a vortex mixer and incubate for 30 minutes at $18...25^\circ\text{C}$ .
9	Centrifuge all tubes except totals for 30 minutes at 2000 g (minimum) at $4...20^\circ\text{C}$ ( $4^\circ\text{C}$ is recommended).
10	Decant the supernatants from all tubes except the totals by simultaneous inversion. Gently blot the tubes against absorbent material and allow to drain for a few seconds to remove any remaining liquid. Care should be taken not to decant any of the precipitate. <b>Do not reinvert the tubes once they have been turned upright.</b> Failure to decant and blot tubes may result in poor replication and erroneous values.
11	Count each tube using a gamma counter for at least 1 minute or until 10 000 counts per tube have accumulated.
12	<b>Calculate</b> the results as described in the "Results' interpretation" section.

## Assay procedure – summary

(all values given in  $\mu\text{l}$ )

	Total	NSB	Calibrator	Control & unknown
Pipette sample		100	100	100
Pipette tracer (red)	200	200	200	200
Pipette antiserum (blue)			200	200
Pipette distilled water		200		
Mix		x	x	x
Incubate 2h (37°C)		x	x	x
Pipette separation reagent		500	500	500
Mix		x	x	x
Incubate for 30 min at 18...25°C		x	x	x
Centrifuge for 30 min at 2000 g		x	x	x
Decant		x	x	x
Count 1 min. or 10 000 counts	x	x	x	x

## 8. QUALITY CONTROL

An internal quality assurance programme is recommended for all clinical laboratories. For this reason, it is suggested that control sera in both the lower and upper portions of the calibration curve be analysed to enable the performance of the procedure to be monitored. The UniQ ICTP RIA assay provides two control materials: ICTP RIA Control 1 & 2. The expected ranges of the controls are indicated on a separate sheet provided with the assay. It is recommended that each laboratory establish its own acceptable ranges for the controls it uses. The mean of the control should fall within the acceptable range stated by Aidian. It is advised to use at least two separate kit lot numbers for the control value assignment. Once the control ranges have been established, controls should be run in duplicate with each assay.

## 9. RESULTS' INTERPRETATION

For automatic result processing, spline function curve fitting is recommended. Results can also be produced manually on semi-log graph paper.

- Calculate the mean count for all calibrators, samples and controls.
- Calculate the  $\%B/B_0$  from:


$$\%B/B_0 = \frac{\text{calibrator or sample count} - \text{NSB}}{\text{calibrator 0 count} - \text{NSB}} \times 100$$

where B is the bound radioactivity of a sample or calibrator,  $B_0$  is the bound radioactivity of the calibrator 0 and NSB is the radioactivity of the non-specific binding.

- Draw a calibration curve on semi-log graph paper with  $\%B/B_0$  values on the ordinate (Y axis) and ICTP concentrations ( $\mu\text{g/l}$ ) of the calibrators on the abscissa (X axis).
- Determine the ICTP concentrations of the unknowns from the calibration curve.
- Any sample reading lower than the lowest calibrator should be reported as such.

- Any sample reading higher than the highest calibrator should be diluted and re-assayed. After dilution, values must be corrected for the dilution factor.

**Fig. 1 Typical calibration curve**  
This calibration curve is for illustration only.

**Table 1 Calculation of results using typical data.**  
This data is for illustration only and must not be used for the calculation of any sample result.

Tube	mean cpm	cpm-NSB	%B/B <sub>0</sub>	ICTP (µg/l)
Total	39445			
NSB	1397			
Cal 0	17333	15936		
Cal 1.0	15875	14478	90.9	
Cal 2.5	14078	12681	79.6	
Cal 5.0	11215	9818	61.6	
Cal 10	8489	7092	44.5	
Cal 25	5482	4085	25.6	
Cal 50	3334	1937	12.2	
Unkn	9494	8097	50.8	7.6
Unkn	6046	4649	29.2	20.9

## 10. LIMITATIONS OF THE PROCEDURE

There is a diurnal variation of ICTP the values being higher at night. Renal insufficiency, glomerular filtration rate of 50 ml/min/1.73 m<sup>2</sup> or less leads to elevated ICTP concentrations in blood. The liver does not seem to be involved in the metabolism of the serum ICTP antigen.

As a basis for diagnostic and therapeutic decisions the results of any single diagnostic test must be backed up with other clinically relevant data.

## 11. EXPECTED VALUES

Due to ethnic, dietary and age variations, the reference limits given may not apply to all populations. Therefore, each laboratory should establish its own representative reference limits.

ICTP values of apparently healthy adults (19–74 years) were measured with the ICTP RIA kit. The non-parametric method recommended by IFCC was used to calculate the reference intervals for the upper and lower limits.

**Table 2 Serum ICTP reference limits**

	Female	Male
Number of subjects	146	125
Mean	3.5 µg/l	3.2 µg/l
Reference interval*	2.1–5.6 µg/l	2.1–5.0 µg/l
0.90-confidence intervals for lower reference limit	2.1–2.4 µg/l	1.9–2.2 µg/l
for upper reference limit	4.9–6.2 µg/l	4.4–5.7 µg/l

\*Reference interval = 0.025 & 0.975 fractiles

ICTP values for children were derived from a population of healthy children aged 0 months to 18 years (unpublished data). The values are indicative only.

**Table 3 Typical serum ICTP values for children**

Age groups		Lower limit (µg/l)	Mean (µg/l)	Upper limit (µg/l)
Boys and girls	0–3 months	51	61	71
	3–6 months	23	32	40
	6–12 months	15	22	30
	1–2 years	12	16	20
	2–9 years	7	11	16
Boys	9–13 years	8	11	15
	13–16 years	10	16	23
	16–18 years	8	9	11
Girls	9–11 years	9	12	15
	11–14 years	10	13	16
	14–16 years	6	10	14
	16–18 years	4	5	7

## 12. PERFORMANCE CHARACTERISTICS

The results presented below were obtained at Aidian and represent typical performance data.

### Measurement range

1.0–50 µg/l.

### Limit of detection

Better than 0.6 µg/l.

### Limit of quantitation

Approximately 1.8 µg/l.

### Linearity

Five high concentration serum samples were serially diluted with five low concentration serum samples. Mean recovery of the samples was 102 %.

### Recovery

Known amounts of ICTP were added to three low concentration patient serum samples. Mean recovery was 109 %.

### Precision

Intra- and inter-assay variations were determined using serum pools containing different concentrations of ICTP.

**Table 4 Intra-assay / Inter-assay precision**

Intra-assay precision of 20 replicates			Inter-assay precision of 10 duplicate determinations		
Sample	Mean ( $\mu\text{g/l}$ )	CV (%)	Sample	Mean ( $\mu\text{g/l}$ )	CV (%)
1	3.3	7	4	4.3	10
2	7.8	5	5	12.1	8
3	45.0	5	6	24.4	8

**Interfering substances**

Serum bilirubin concentration < 400  $\mu\text{mol/l}$  does not interfere. Serum haemoglobin concentration up to 5 g/l does not interfere. Triglycerides < 5 g/l do not interfere.

**Specificity**

The antigenic determinant of the antibody is within the two hydrophobic phenylalanine-rich sequences of the tri-valently cross-linked ICTP. A trivalent cross-link is necessary because the divalently cross-linked and monomeric peptides containing only one phenylalanine-rich sequence show poor immunoreaction. The ICTP RIA antiserum does not cross-react with antigens used in other UniQ Collagen assays.

**13. TRACEABILITY**

The measurand (analyte) in the UniQ ICTP RIA Calibrators is traceable to the manufacturer's working calibrators which are prepared from highly purified ICTP whose value is assigned using amino acid composition analysis.

**14. DISPOSAL**

- Dispose of contents according to national and local law.
- All patient samples and controls should be handled and disposed of as potentially infectious material.
- Radioactive waste material must be discarded according to the applicable laws.
- Materials of the components:
  - Paper: Instructions for use
  - Cardboard: Kit box
  - Plastic: Separation reagent bottle, vial caps
  - Rubber: Stoppers
  - Glass: Tracer, antiserum, calibrator, control vials
  - Foam rubber
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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## Literature

- Abildgaard N et al.** Comparison of five biochemical markers of bone resorption in multiple myeloma: elevated pre-treatment levels of s-ICTP and u-Ntx are predictive for early progression of the bone disease during standard chemotherapy. *Br J Haematol* 2003;120:235-242.
- Abildgaard N et al.** Serum markers of bone metabolism in multiple myeloma: prognostic value of the carboxy-terminal telopeptide of type I collagen (ICTP). *Br J Haematol* 1997;96:103-110.
- Aruga A et al.** Usefulness of bone metabolic markers in the diagnosis and follow-up of bone metastasis from lung cancer. *Br J Cancer* 1997;76:760-764.
- Berruti A et al.** Differential patterns of bone turnover in relation to bone pain and disease extent in bone in cancer patients with skeletal metastases. *Clin Chem* 1999;45:1240-1247.
- Engler H et al.** Diagnostic and prognostic value of biochemical markers in malignant bone disease: A prospective study on the effect of bisphosphonate on pain intensity and progression of malignant bone disease. *Clin Chem Lab Med* 1998;36:879-885.
- Fonseca R et al.** Prognostic value of serum markers of bone metabolism in untreated multiple myeloma patients. *Br J Haematol* 2000;109: 24-29
- Hakala M et al.** Application of markers of collagen metabolism in serum and synovial fluid for assessment of disease process in patients with rheumatoid arthritis. *Ann Rheum Dis* 1995;54:886-890.
- Hakala M et al.** Increased type I collagen degradation correlates with disease severity in rheumatoid arthritis. *Ann Rheum Dis* 1993;52:866-869.
- Horiguchi T et al.** Usefulness of serum carboxy-terminal telopeptide of type I collagen (ICTP) as a marker of bone metastasis from lung cancer. *Jpn J Clin Oncol* 2000;30:174-179.
- Jakob C et al.** Serum levels of carboxy-terminal telopeptide of type-I collagen are elevated in patients with multiple myeloma showing skeletal manifestations in magnetic resonance imaging but lacking lytic bone lesions in conventional radiography. *Clin Cancer Res* 2003;9:3047-3051.
- Keskikuru R et al.** Preoperative high type I collagen degradation marker ICTP reflects advanced breast cancer. *Anticancer Res* 1999;19:4481-4484.
- Kobayashi T et al.** Type I collagen metabolites as tumor markers in patients with lung carcinoma. *Cancer* 1999;85:1951-1957.
- Koga H et al.** Use of bone turnover marker, pyridinoline cross-linked carboxyterminal telopeptide of type I collagen (ICTP), in the assessment and monitoring of bone metastasis in prostate cancer. *The Prostate* 1999;39:1-7.
- Koizumi M et al.** Bone metabolic markers in bone metastasis of breast cancer. *Int J Clin Oncol* 1999;4:331-337.
- Koizumi M et al.** Bone metabolic markers in the evaluation of bone scan flare phenomenon in bone metastases of breast cancer. *Clin Nucl Med* 1999;24:15-20.
- Kymälä T et al.** Type I collagen degradation product (ICTP) gives information about the nature of bone metastases and has prognostic value in prostate cancer. *Br J Cancer* 1995;71:1061-1064.
- Maeda H et al.** Correlation between bone metabolic markers and bone-scan in prostatic cancer. *J Urology* 1997;157:539-543.
- Maemura M et al.** Serum concentration of pyridinoline cross-linked carboxyterminal telopeptide of type I collagen in patients with metastatic breast cancer. *Oncology Reports* 2000;7:1333-1338.
- Noguchi M et al.** Serum levels of bone turnover markers parallel the results of bone scintigraphy in monitoring bone activity of prostate cancer. *Urology* 2003;61:993-998.

**Okamura T et al.** The efficacy of a serum carboxyterminal pyridinoline cross-linked telopeptide of type I collagen as a quantitative screening marker for bone metastases in patients with urological malignancies. *Cancer Det Prev* 2003;27:14-18.

**Rauchenzauner M et al.** Sex- and age-specific reference curves for serum markers of bone turnover in healthy children from 2 months to 18 years. *J Clin Endocrinol Metab* 2007;92: 443-449.

**Risteli J et al.** Radioimmunoassay for the pyridinoline cross-linked carboxyterminal telopeptide of type I collagen: A new serum marker of bone resorption. *Clin Chem* 1993;39:635-640.

**Santala M et al.** Type I and type III collagen metabolites as predictors of clinical outcome in epithelial ovarian cancer. *Clin Cancer Res* 1999;5:4091-4096.

**Sassi M-L et al.** Assay for cross-linked carboxyterminal telopeptide of type I collagen (ICTP) unlike CrossLaps assay reflects increased pathological degradation of type I collagen in rheumatoid arthritis. *Clin Chem Lab Med* 2003;41:1038-1044.

**Sassi M-L et al.** Immunochemical characterisation of assay for carboxyterminal telopeptide of human type I collagen: Loss of antigenicity by treatment with cathepsin-K. *Bone* 2000;26:367-373.

**Shimozuma K et al.** Biochemical markers of bone turnover in breast cancer patients with bone metastases: A preliminary report. *Jpn J Clin Oncol* 1999;29:16-22.

**St Clair et al.** A cross sectional analysis of 5 different markers of collagen degradation in rheumatoid arthritis. *J Rheumatol* 1998;25:1472-1479.

**Voorzanger-Rousselot N et al.** Association of 12 serum biochemical markers of angiogenesis, tumour invasion and bone turnover with bone metastases from breast cancer: a cross-sectional and longitudinal evaluation. *Br J Cancer* 2006;95:506-514.

**Ylissirniö S et al.** Serum type I collagen degradation markers, ICTP and CrossLaps, are factors for poor survival in lung cancer. *Anticancer Res* 1999;19:5577-5582.

**Åman S et al.** Prediction of disease progression in early rheumatoid arthritis by ICTP, RF and CRP. A comparative 3-year follow-up study. *Rheumatology* 2000;39:1009-1013.

UniQ® is a registered trademark of Aidian Oy.

**AIDIAN**

Aidian Oy  
Koivu-Mankkaan tie 6 B  
P.O. Box 83, FI-02101 Espoo, Finland  
[www.aidian.eu](http://www.aidian.eu)