AIDIAN

UniQ[®] PINP RIA







Intact N-terminal propeptide of type I
procollagen, radioimmunoassay

Cat. No. 67034

132059-20

**Explanation of symbols • Erläuterung der Symbole • Explication des symboles
Verklaring van symbolen • Explicación de los símbolos**

	English	German	Français	Nederlands	Español
	<i>In vitro</i> diagnostic medical device	<i>In-vitro</i> Diagnostikum	Dispositif médical de diagnostic <i>in vitro</i>	Alleen voor <i>in vitro</i> diagnostisch gebruik	Producto sanitario para diagnóstico <i>in vitro</i>
	Catalogue number	Artikelnummer	Référence	Catalogus nummer	Número de catálogo
	Batch code	Chargen-bezeichnung	Lot	Batch code	Número de lote
	Use by	Verwendbar bis	Utiliser jusqu'au	Gebruik door	Fecha de caducidad
	Temperature limit	Temperaturbegrenzung	Limites de température	Temperatuurlimiet	Limitación de temperatura
	Consult instructions for use	Nachschlagen in der Bedienungsanleitung	Consulter la notice d'utilisation	Raadpleeg gebruiksinstructie	Consulte las instrucciones de uso
	Manufacturer	Hersteller	Fabricant	Producent	Fabricante
	Sufficient for	Ausreichend für	Suffisant pour	Volgende voor	Válido para
	Tracer	Tracer	Traceur	Tracer	Trazador
	Antiserum	Antiserum	Antiserum	Antiserum	Antisuero
	Separation Reagent	Trennreagenz	Réactif de séparation	Scheidingsreagens	Reactivo de separación
	Control	Kontrolle	Contrôle	Controle	Control
	Calibrator concentration	Kalibrator-konzentration	Concentration de l'étalon	Kalibrator concentratie	Concentración del calibrador
	Origin: rabbit	Ursprung: Kaninchen	Origine: lapin	Oorsprong: konijn	Origen: conejo
	Origin: human	Ursprung: menschlich	Origine: humaine	Oorsprong: humane	Origen: humano
	Contains sodim azide	Enthält Natriumazid	Contient de l'azide de sodim	Bevat natrium azide	Contiene azida sódica
	Authorised representative in Switzerland	Autorisierter Vertreter in der Schweiz	Représentant autorisé en Suisse	Gemachtigde vertegenwoordiger in Zwitserland	Representante autorizado en Suiza

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Intact N-terminal propeptide of type I procollagen, radioimmunoassay

1. INTENDED USE

UniQ PINP RIA is a quantitative radioimmunoassay designed for the measurement of the concentration of intact aminoterminal propeptide of type I procollagen, an indicator of osteoblastic activity, in human serum. The test is intended to be used as an aid in the management of postmenopausal osteoporosis.

For *in vitro* diagnostic use.

2. SUMMARY AND EXPLANATION

Type I collagen, which is synthesised by fibroblasts and osteoblasts, is the most abundant collagen type in the body and the only collagen type found in mineralised bone, where it accounts for more than 90 % of the organic matrix. Since bone is the major collagenous organ and metabolically highly active throughout life, the majority of the synthesised type I collagen stems from bone osteoblasts.

Bone collagen is derived from a larger protein – type I procollagen –, which consists of three amino acid chains that are intertwined to form a rod-like triple helix. Type I procollagen has propeptide extensions at both ends of the molecule, which are removed by specific proteinases before the collagen molecules thus formed are assembled into collagen fibres. Both propeptides can be found in the circulation, where their concentration reflects the synthesis rate of type I collagen. The UniQ PINP RIA assay measures the serum concentration of the aminoterminal propeptide of type I procollagen (PINP).¹ As the concentration of this extension propeptide is directly proportional to the amount of new collagen laid down in bone, it can be used to assess bone formation.

During bone formation, the bone matrix is produced before mineralisation occurs; hence PINP is an early marker of bone formation. Bone alkaline phosphatase, a cell membrane enzyme present in preosteoblasts and osteoblasts, is another early marker of osteoblast activity.² Osteocalcin is a later marker of osteoblast differentiation that appears in the bone mineralisation phase.³

Bone tissue is remarkably metabolically active and throughout life undergoes constant remodelling. Bone remodelling is achieved by two counteracting processes: bone formation and bone resorption, which under normal conditions are tightly coupled to each other. Metabolic bone diseases are characterised by imbalances in bone turnover and often lead to an uncoupling between bone formation and resorption.⁴ An indicator of type I collagen turnover, such as PINP, is very useful for investigating skeletal remodelling under normal and abnormal conditions.

The UniQ PINP assay provides a sensitive tool for assessing increased bone turnover in postmenopausal women⁵⁻⁹ and may be particularly useful for monitoring the effect of antiresorptive and anabolic therapy on bone metabolism.¹⁰⁻²¹ The determination of PINP concentrations has also been used to detect increases in type I collagen turnover in disease states such as Paget's disease of bone²²⁻²⁴ and bone metastasis of cancer. PINP determination may be useful in assessing bone metastatic activity and in predicting survival.²⁵⁻²⁷

3. PRINCIPLES OF THE TEST

The UniQ PINP RIA kit is based on the competitive radioimmunoassay technique. A known amount of labelled PINP and an unknown amount of unlabelled PINP in the sample compete for a limited number of high affinity binding sites on the polyclonal rabbit anti-PINP antibody. A second antibody, directed against rabbit IgG and coated to kaolin particles, is used to separate the antibody-bound PINP from free PINP. The radioactivity of the bound tracer antigen

is measured in a gamma counter. The amount of labelled PINP in the sample tube is inversely proportional to the amount of PINP in the sample. The concentrations in unknown samples are obtained from a calibration curve, which is based on the concurrent testing of the UniQ PINP RIA calibrators that contain 0 to 250 µg/l PINP.

4. REAGENTS

UniQ® PINP RIA, Cat. No. 67034

Component name and origin	Symbol	Σ 100	Stability at 2...8°C	
			Unopened	Opened
Tracer	Ag ¹²⁵ I	22 ml Radioactivity < 200 kBq	Until the expiration date marked on the label	6 weeks
Antiserum ORIG RAB	Ab	22 ml	Until the expiration date marked on the label	6 weeks
Calibrator	CAL xxµg/l	7 x 0.5 ml	Until the expiration date marked on the label	6 weeks
Control ORIG HUM	CONTROL #	2 x 1.0 ml	Until the expiration date marked on the label	6 weeks after reconstitution.
Procollagen Separation Reagent	REAG SEPAR	60 ml	Until the expiration date marked on the label	6 weeks
Instructions for use				

Reagents contain sodium azide as a preservative. See section 5 "Warnings and precautions".

Reagent preparation

The controls are lyophilized. All other reagents are ready to use.

Allow the **control** vials to reach room temperature (18...25°C) before opening. Reconstitute the controls by adding 1.0 ml distilled or deionised water to each vial. Cap and mix well by gentle swirling or inverting to avoid foaming. Allow to stand for 30 minutes before use.

Reagent deterioration

Alterations in the physical appearance of the reagents or in the slope of the calibrator curve, or values of control sera outside the manufacturer's acceptable range may be an indication of reagent instability. Do not use reagents if you suspect reagent deterioration.

5. WARNINGS AND PRECAUTIONS

Health and safety information

- For *in vitro* diagnostic use only.
- **Rx only:**
 - Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner.
 - For prescription use only.
 - For professional use only.
- Do not pipette by mouth.
- All patient samples and controls should be handled as potentially infectious material.
- The controls contain material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle the controls as if capable of transmitting infectious agents.
- Lyophilised reagents contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Contact with acids liberates very toxic gas (EUH032). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Azides may react with metal plumping, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration.
- The kit contains ¹²⁵I (half-life $t_{1/2}$ of 60 days) emitting 35.5 keV gamma rays and 27–32 keV X-rays. Radioactive material contained in this kit may be received, acquired, possessed and used only by physicians, clinical laboratories, research institutions or hospitals and used only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom to human beings or animals. All radioactive materials should be stored and handled only in specifically designated areas. The purchase, storage, use or exchange of radioactive products is subject to the laws in force in the user's country.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of water and decontaminate, following your institution's guidelines for accidental spills.
- Spills of radioactive materials should be cleaned immediately with a detergent and the contaminated cleaning materials transferred to the radioactive waste container or handled according to the recommendations of the radiation safety office.
- Refer to the radiation regulations and procedures applicable to your institution for handling, storage and disposal of radioactive materials.

Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Do not mix components from different lot numbers or different tests.

6. SAMPLE COLLECTION AND PREPARATION

No special preparation of the patient is required prior to sample taking.

Serum should be used and a normal venous blood sample should be collected using an acceptable medical technique, following the normal precautions for venipuncture. The assay is not to be used with plasma samples.

Highly lipemic, hemolyzed or microbiologically contaminated patient samples should be avoided.

Sample handling and storage

After the specimen collection, allow the blood to clot and separate serum by centrifugation. Serum samples may be stored for up to 5 days at 2...8°C. For longer periods, store the samples at -20°C or colder. Repeated freezing and

thawing should be avoided.

Samples with high PINP concentrations may be diluted using 0.9 % saline.

7. PROCEDURE

See the section 4, “Reagents”, for a list of reagents.

Materials required but not provided

- Pipettes (for accurately delivering 50 µl, 200 µl, 500 µl and 1000 µl)
- Test tubes (polystyrene, e.g. 12 x 75 mm)
- Absorbent paper
- Paraffin film
- Vortex mixer
- Water bath (capable of maintaining 37 ±1°C)
- Decanting rack
- Centrifuge (capable of at least 2000 g, preferably refrigerated)
- Gamma counter
- Distilled or deionised water

Instruments

Any gamma counter capable of measuring ¹²⁵I may be used. The counting efficiency of the instrument should be regularly checked as well as the background counts. It is recommended that the instrument be calibrated systematically. Consult the operations manual supplied by the instrument manufacturer for details of installation, use, maintenance and calibration.

Details of the procedure

1	Bring all reagents and samples to room temperature (18...25°C) prior to use. Mix all reagents well by gentle inversion prior to use.
2	Label polystyrene test tubes in duplicate for NSB (non-specific binding), calibrators, controls, patient samples and total counts.
3	Pipette 50 µl of calibrator, control or patient sample into the appropriate tubes. Any patient sample or the 0 µg/l calibrator may be pipetted into the NSB tubes.
4	Pipette 200 µl of Tracer (red) into all tubes.
5	Pipette 200 µl of Antiserum (blue) into all tubes except NSB and total counts. Pipette 200 µl of distilled or deionised water into the NSB tubes.
6	Mix all tubes on a vortex mixer, cover with paraffin film and incubate for 2 hours at 37 ±1°C.
7	Mix the Procollagen Separation Reagent thoroughly by gentle inversion. Add 500 µl to all tubes except total counts.
8	Mix the tubes using a vortex mixer and incubate for 30 minutes at room temperature (18...25°C).
9	Centrifuge all tubes except total counts for 15 minutes at 2000 g at 4...20°C (4°C is recommended).
10	Decant the supernatants from all tubes except total counts by simultaneous inversion. Gently blot tubes against absorbent material and allow to drain for a few seconds to remove any remaining liquid. Care should be taken not to decant any of the precipitate. Do not reinvert the tubes once they have been turned upright. Failure to decant and blot tubes may result in poor replication and erroneous values.
11	Count each tube in a gamma counter for 1 minute.
12	Calculate the results as described in the “Calculation of Results” section.

Assay procedure – summary

(All volumes given in µl)

	Total	NSB	Calibrator	Control and unknown
Pipette sample		50	50	50
Pipette tracer (red)	200	200	200	200
Pipette antiserum (blue)			200	200
Pipette distilled / deionised water		200		
Mix		x	x	x
Incubate for 2 hours at 37°C		x	x	x
Pipette separation reagent		500	500	500
Mix		x	x	x
Incubate for 30 min at room temperature (18...25°C)		x	x	x
Centrifuge for 15 min at 2000 g		x	x	x
Decant		x	x	x
Count 1 min	x	x	x	x

8. QUALITY CONTROL

An internal quality assurance programme is recommended for all clinical laboratories. For this reason, it is suggested that control sera in both the lower and upper portions of the calibration curve be analysed to enable the performance of the procedure to be monitored. The UniQ PINP RIA assay provides two control materials: PINP RIA Control 1 & 2. The expected ranges of the controls are indicated on a separate sheet provided with the assay. It is recommended that each laboratory establish its own acceptable ranges for the controls it uses. The mean of the control should fall within the acceptable range stated by Aidian. It is advised to use at least two separate kit lot numbers for the control value assignment. Once the control ranges have been established, controls should be run in duplicate with each assay.

9. RESULTS' INTERPRETATION

For automatic result processing spline function curve fitting is recommended. Results can also be produced manually on semi-log graph paper.

- **Calculate** the mean count for all calibrators, samples, controls and non-specific binding.
- **Calculate** the %B/B₀ from:

$$\%B/B_0 = \frac{\text{calibrator or sample count} - \text{NSB}}{\text{calibrator 0 count} - \text{NSB}} \times 100$$

where B is the bound radioactivity of a sample or calibrator, B₀ is the bound radioactivity of the calibrator 0 and NSB is the radioactivity of the non-specific binding.

- **Draw** a calibration curve on semi-log graph paper with %B/B₀ values on the ordinate (Y axis) and PINP concentrations (µg/l) of the calibrators on the abscissa (X axis).

- **Determine** the PINP concentrations of the unknowns from the calibration curve.
- **Any sample reading lower than the lowest calibrator should be reported as such.**
- **Any sample reading higher than the highest calibrator should be diluted and re-assayed.** After dilution, values must be corrected for the dilution factor.

Fig. 1 Typical calibration curve.
This calibration curve is for illustration only.


Table 1 Calculation of results using typical data.
This data is for illustration only and must not be used for the calculation of any sample result.

Tube	mean cpm	cpm-NSB	%B/B ₀	PINP (µg/l)
Total	50970			
NSB	845			
Cal 0	28330	27485		
Cal 5.0	25343	24498	89.1	
Cal 10	22143	21289	77.5	
Cal 25	15344	14499	52.8	
Cal 50	9897	9052	32.9	
Cal 100	6095	5250	19.1	
Cal 250	3379	2534	9.2	
Unkn 1	10725	9880	35.9	44.6
Unkn 2	4572	3727	13.6	157

10. LIMITATIONS OF THE PROCEDURE

Results of this assay should be used in conjunction with other pertinent clinical information for making diagnostic and therapeutic decisions. Therefore, the UniQ PINP assay is not recommended for use as a screening procedure to detect the presence of osteoporosis in the general population.

When evaluating subsequent samples, collect the samples at the same time of the day, as there is a diurnal variation of PINP with the values being higher at night.

Type I collagen is also synthesised by fibroblasts and therefore can be found in loose connective tissues together with other collagen types. However, since bone is the major collagenous organ and also metabolically active throughout life, the majority of the circulating PINP is of osteoblast origin.²⁸

PINP is metabolised in the liver.²⁹ Severe liver disease may affect clearance from the circulation and give rise to elevated serum PINP levels.³⁰

As with most assays utilising secondary antibodies, samples containing antibodies (such as heterophilic antibodies) that react with rabbit/goat IgG may cause erroneous results.

11. EXPECTED VALUES

Due to ethnic, dietary and age variations, the reference limits given may not apply to all populations. The ranges of observations reported below are representative only, and do not necessarily reflect the ranges that will be observed in a particular clinical setting.

PINP values of apparently healthy white Finnish adults were measured with the UniQ PINP RIA kit. The non-parametric method recommended by IFCC was used to calculate the reference intervals for the upper and lower limits.³¹ Each laboratory should establish its own representative expected values for the local population.


Table 2 Serum PINP expected values

	Female	Female Postmenopausal *	Male
Age, years	25–44	52–62	22–65
Number of subjects	210	383	323
Mean	42 µg/l	48 µg/l	45 µg/l
Median	39 µg/l	45 µg/l	43 µg/l
Reference interval**	19–83 µg/l	16–96 µg/l	22–87 µg/l
0.95-confidence intervals			
for lower reference limit	14–23 µg/l	15–19 µg/l	20–25 µg/l
for upper reference limit	71–96 µg/l	87–103 µg/l	77–94 µg/l

*Min. 3 years postmenopausal

**Reference interval = 0.025 & 0.975 fractiles

Fig. 2 Box-Whisker plots of the distribution of PINP concentrations in apparently healthy adults. The horizontal line in each box represents the median value while the lower and upper limits of the boxes indicate the interquartile ranges, i.e., values 25 % above and below the median. The tips of the whiskers show the reference range limits (0.025 & 0.975 fractiles).

The results presented above demonstrate that PINP concentrations in a population of postmenopausal women are elevated compared to concentrations in women aged 25–44 years ($p = 0.0005$, Mann Whitney). This increase in PINP reflects the increase in bone remodelling associated with oestrogen deficiency in a postmenopausal population compared to a younger population. However, there is considerable overlap in PINP concentrations in the two populations.

CLINICAL STUDIES

Bisphosphonate treatment

The usefulness of UniQ PINP RIA to monitor changes in bone remodelling during bisphosphonate treatment (10 mg daily alendronate) was studied in a multi-site, prospective, double blind placebo-controlled study in postmenopausal women of the US population with defined osteoporosis, at least 5 years since menopause.³² Data were available from 124 patients randomised into two groups: 62 receiving the treatment (TRM) and 62 receiving placebo (CTR). All study subjects received 500 mg/d calcium supplementation. Patients were analysed at baseline and at 3, 6 and 12 months of treatment. Bone mineral density (BMD) was measured with dual-energy X-ray absorptiometry (DXA) and intact PINP was measured by UniQ PINP RIA.

The data are presented in Table 3 and Figure 3. At baseline, the BMD was 0.740 ± 0.083 g/cm² (mean \pm SD) in the CTR and 0.727 ± 0.104 g/cm² in the TRM group. In the TRM group the mean BMD increased by 6.0 % from baseline at 12 months. The change was significant at all time points compared to baseline ($p < 0.0001$). The increase in the CTR group was smaller and transient returning back to baseline by 12 months.


At baseline the concentration of PINP was 62.5 ± 18.5 μ g/l in the CTR group and 61.1 ± 30.1 μ g/l in the TRM group. The mean concentration of PINP in the TRM group decreased by 70.8 % at 12 months of treatment. At all three time points the decrease was significant ($p < 0.0001$) compared to baseline. The decrease in the CTR group, also receiving the calcium supplementation, was smaller, at 9.8 % at 12 months.

The changes observed in BMD and PINP during bisphosphonate treatment reflect the effects this treatment has on bone remodelling. The results were as expected and similar to earlier reports, which have shown significant decreases in PINP levels during bisphosphonate treatment.^{11, 33-34}

Table 3 Response expressed as a mean percentage change in BMD and PINP in patients receiving bisphosphonate or placebo

Test	Treatment group	Mean percentage change (SD) from baseline value		
		3 months	6 months	12 months
BMD	Placebo	+0.3 (3.1)	+1.2 (3.3)	-0.1 (3.2)
	Bisphosphonate	+3.5 (3.2)	+5.1 (4.4)	+6.0 (3.3)
PINP	Placebo	-8.6 (25.9)	-14.0 (25.5)	-9.8 (27.4)
	Bisphosphonate	-55.4 (34.0)	-71.6 (15.6)	-70.8 (15.5)

Fig. 3 Response as mean percent change (\pm SE bars) in (A) BMD and (B) PINP concentration in patients receiving bisphosphonate (●) or placebo (○).

In order to have clinical significance, the percentage change between two measurements must exceed both the analytical (CV_a) and intra-individual (CV_i) variation. This threshold is referred to as the minimum percentage change and it can be calculated as follows:

$$\text{Minimum percentage change} = \sqrt{2} \times 1.96 \times \sqrt{(CV_i^2 + CV_a^2)}$$

Several values for PINP have been published in literature. For example, Hannon and colleagues¹² have reported a minimum percentage change of 21 % for PINP.

Hormone replacement therapy

The ability of the UniQ PINP RIA to reflect changes in bone remodelling during hormone replacement therapy (HRT) in postmenopausal women was studied in a randomised, double blind dose-finding study.³⁵⁻³⁶ Study subjects were apparently healthy Finnish women at least 3 years since menopause. Data were available from the treatment arm of 44 subjects receiving the lowest dose of treatment: 1 mg oral oestradiol valerate and 2.5 mg oral medroxyprogesterone acetate daily. Patients were analysed at baseline and at 6 and 12 months of treatment. BMD was measured with dual-energy X-ray absorptiometry and intact PINP was measured using UniQ PINP RIA.


The data are presented in Table 4 and Figure 4. At baseline the BMD was 1.083 ± 0.182 g/cm² (mean \pm SD) and the concentration of PINP was 42.2 ± 15.3 μ g/l. At 6 months BMD had increased by mean 1.4 % and at 12 months the increase was 3.1 % from the baseline value. PINP decreased from the baseline value by 12.3 % at 6 months and by 19.8 % at 12 months. The changes in BMD and PINP were significant ($p < 0.002$) at both time points compared to the baseline value.

The changes observed in BMD and intact PINP during HRT reflect the effect this treatment has on bone remodelling even with the low dose treatment used in this study. The results were as expected and similar to earlier reports which all showed significant decreases in PINP levels during different types of hormone replacement therapies.^{12-13, 15}

Table 4 Response expressed as a mean percentage change in BMD and PINP concentration in postmenopausal patients receiving low dose combined hormone replacement therapy.

Test	Mean percent change (SD) from the baseline value	
	6 months	12 months
BMD	+1.4 (2.3)	+3.1 (3.8)
PINP	-12.3 (39.4)	-19.8 (46.5)

Fig. 4 Response expressed as a mean percentage change (\pm SE bars) in BMD (o) and PINP concentration (●) in postmenopausal women receiving low dose combined hormone replacement therapy.

Teriparatide treatment

Teriparatide, human parathyroid hormone (PTH) 1–34, of recombinant DNA origin, stimulates bone formation. The ability of the UniQ PINP RIA to reflect changes in bone remodelling following teriparatide therapy in the management of osteoporosis was studied at 19 clinical sites using a randomised double blind design.³⁷ Study subjects were ambulatory post-menopausal women with defined osteoporosis, at least 5 years since menopause. A subset of subjects received once-daily doses of 20 µg teriparatide with daily supplementation of calcium (1000 mg) and vitamin D (400–800 IU) throughout the study. Retrospectively, data were collected from 62 patients age 64.9 ± 7.8 years (mean \pm SD) receiving teriparatide and having full sets of BMD and PINP results available at the scheduled sampling points. Lumbar spine (L1–L4) BMD was measured in all patients by DXA using Lunar or Hologic densitometers at baseline and 3, 6, 12 and 18 months of teriparatide treatment. Intact PINP was measured at baseline and 1, 3, 6 and 12 months of treatment using the UniQ PINP RIA assay.

Teriparatide therapy resulted in a mean BMD increase throughout the study period with the increases being statistically significant at all time points ($p < 0.0001$) as compared with baseline. Mean BMD increased by 0.025 g/cm² at 3 months, 0.039 g/cm² at 6 months, 0.061 g/cm² at 12 months and 0.079 g/cm² at 18 months.


The mean PINP concentration increased by 59 µg/l from baseline at 1 month, 82 µg/l at 3 months, 117 µg/l at 6 months and 103 µg/l at 12 months. The increase in the PINP concentration compared to baseline was statistically significant at all time points ($p < 0.0001$).

The percentage change from baseline for BMD and PINP was calculated at each time point. The mean percentage changes are shown in Table 5 and graphically depicted in Figure 5. The mean values for PINP and BMD increased in response to the teriparatide treatment. PINP peaked at 6 months followed by a slight decline at 12 months as compared to the 6-month time point. Mean BMD increased in response to treatment with the maximum value being observed at the end of the study period.

Table 5 Response expressed as a mean percentage change (\pm SD) in BMD and PINP concentration in 62 postmenopausal women with osteoporosis receiving once-daily doses of 20 μ g teriparatide.

Test	Change (mean \pm SD)	Change (mean \pm SD)	Change (mean \pm SD)	Change (mean \pm SD)	Change (mean \pm SD)
	1 month	3 months	6 months	12 months	18 months
BMD	NA	3.3 \pm 3.4 %	5.2 \pm 3.8 %	8.1 \pm 5.1 %	10.6 \pm 5.6 %
PINP	132.6 \pm 118.2 %	187.8 \pm 164.2 %	278.5 \pm 257.1 %	249.5 \pm 221.7 %	NA

Fig. 5 Response expressed as a mean percentage change (\pm SE bars) in BMD and PINP concentration in 62 postmenopausal women with osteoporosis receiving once-daily doses of 20 μ g teriparatide

Changes observed in BMD and intact PINP levels reflect the anabolic effect of teriparatide therapy on bone remodelling. These results were as expected and similar to other reports, which have shown significant increases in PINP levels during teriparatide treatment.¹⁷⁻²⁰ Partly based on the same clinical material as the above study²⁷, an algorithm for using PINP to monitor patients treated with teriparatide has been developed, which underlines the clinical utility of PINP measurements in the management of postmenopausal osteoporosis.²¹

12. PERFORMANCE CHARACTERISTICS

The results presented below were obtained at Aidian and represent typical performance data.

Linearity

The linearity of the UniQ PINP RIA assay was studied by diluting samples with the 0 µg/l calibrator of the kit. The recovery percentage was calculated by dividing the measured concentration by the expected concentration and multiplying by 100.

Table 6 Linearity of the PINP assay

Sample	Dilution	Expected (µg/l)	Measured mean (µg/l, n=4)	Calculated for undiluted (µg/l)	Recovery (%)	
1	Neat	–	211.5	–	–	
	4:5	169.2	153.3	191.6	91	
	3:5	126.9	103.7	172.8	82	
	2:5	84.6	76.6	191.5	91	
	1:5	42.3	40.3	201.4	95	
				Mean	189.3	
				SD	12	
			CV %	6.3		
2	Neat	–	80.9	–	–	
	4:5	64.7	64.5	80.6	100	
	3:5	48.5	48.0	79.9	99	
	2:5	32.4	32.0	82.4	102	
	1:5	16.2	15.3	76.7	95	
				Mean	79.9	
				SD	2.4	
			CV %	3		

Precision

Intra- and inter-assay precision of the assay was determined using serum pools and control samples containing different concentrations of PINP.

Table 7 Intra-assay / Inter-assay precision

Intra-assay precision of 18 or 20 replicates			Inter-assay precision of 10 duplicate determinations		
Sample	Mean (µg/l)	CV (%)	Sample	Mean (µg/l)	CV (%)
1 (n=20)	26.6	5.4	5	26.0	5.5
2 (n=18)	33.1	3.2	6	32.9	8.3
3 (n=20)	53.1	5.4	7	52.7	9.5
4 (n=20)	149.2	9.6	8	142.9	8.9

Recovery

The recovery was studied by adding known amounts of purified PINP antigen to a sample containing 19.1 µg/l of endogenous PINP. The sample was measured in duplicate with the UniQ PINP RIA assay and the recovery percentage was calculated. The recovery percentage equals the increase in concentration divided by the spike concentration multiplied by 100.

Table 8 The recovery of PINP assay

Spike (µg/l)	Measured (µg/l)	Increase (µg/l)	Recovery (%)
–	19.1	–	–
109.0	120.0	100.9	93
68.8	87.8	68.7	100
41.3	59.1	40.1	97

Specificity

The intact UniQ PINP RIA assay is highly specific for intact trimeric forms of N-terminal propeptide of type I procollagen. The small molecular weight monomeric degradation products of the PINP propeptide (Col1) are detected at low levels (cross-reactivity of 2.7 %). UniQ PINP RIA does not cross-react with the N-terminal propeptide of type III procollagen (PIIINP).

Interfering substances

Serum bilirubin concentration < 400 µmol/l (< 23 mg/dl) does not interfere with the assay. Serum haemoglobin concentration up to 5 g/l (500 mg/dl) or triglycerides up to 30 g/l (3000 mg/dl) do not interfere with the assay.

Detection limit

The detection limit of the UniQ PINP RIA assay was determined by measuring 20 replicates of the 0 µg/l calibrator and the calibration curve (4 replicates) three times. The minimum detectable concentration of PINP in the assay was less than 2 µg/l, defined as twice the standard deviation of the 0-binding value.

Measurement range

The measurement range of the UniQ PINP RIA is 5–250 µg/l.

13. TRACEABILITY

The measurand (analyte) in the UniQ PINP RIA Calibrators is traceable to the manufacturer's working calibrators which are prepared from highly purified PINP whose value is assigned using amino acid composition analysis.

14. DISPOSAL

- Dispose of contents according to national and local law.
- All patient samples and controls should be handled and disposed of as potentially infectious material.
- Radioactive waste material must be discarded according to the applicable laws.
- Materials of the components:
 - Paper: Instructions for use
 - Cardboard: Kit box
 - Plastic: Separation reagent bottle, vial caps
 - Rubber: Stoppers
 - Glass: tracer, antiserum, calibrator, control vials
 - Foam rubber
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

UniQ® PINP RIA Intaktes N-terminales Propeptid von Typ I Prokollagen, Radioimmunoassay Kit

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Literature

1. **Melkko J et al.** Immunoassay for intact amino-terminal propeptide of human type I procollagen. *Clin Chem* 1996;42:947-954.
2. **Kress B et al.** Use of bone alkaline phosphatase to monitor alendronate therapy in individual postmenopausal osteoporotic women. *Clin Chem* 1999;45:1009-1017
3. **Rosenquist C et al.** Measurement of a more stable region of osteocalcin in serum by ELISA with two monoclonal antibodies. *Clin Chem* 1995;41:1439-1445.
4. **Seibel M J.** Bone metabolism, mineral homeostasis and its pharmacological modulation. *Clin Lab* 2004;50:255-264.
5. **Cabrera DC et al.** Biochemical markers of bone formation in the study of postmenopausal osteoporosis. *Osteoporosis Int* 1998;8:147-151.
6. **Saarto T et al.** Aminoterminal propeptide of type I procollagen (PINP) correlates to bone loss and predicts the efficacy of antiresorptive therapy in pre- and post-menopausal non-metastatic breast cancer patients. *Br J Cancer* 1998;78:240-245.
7. **Garnero P et al.** Markers of bone turnover predict postmenopausal forearm bone loss over 4 years: The OFELY study. *J Bone Miner Res* 1999;14:1614-1621.
8. **Scarioano JK et al.** Diagnostic efficacy of serum cross-linked N-telopeptide (NTx) and aminoterminal procollagen extension peptide (PINP) measurements for identifying elderly women with decreased bone mineral density. *Scand J Clin Lab Invest* 2002;62:237-244.
9. **Scarioano JK et al.** Serum levels of cross-linked N-telopeptides and aminoterminal propeptides of type I collagen indicate low bone mineral density in elderly women. *Bone* 1998;23:471-477.
10. **Garnero P et al.** Short-term effects of new synthetic conjugated estrogens on biochemical markers of bone turnover. *J Clin Pharmacol* 2002;42:290-296.
11. **Fink E et al.** Differences in the capacity of several biochemical bone markers to assess high bone turnover in early menopause and response to alendronate therapy. *Osteoporosis Int* 2000;11:295-303.
12. **Hannon R et al.** Response of biochemical markers of bone turnover to hormone replacement therapy: impact of biological variability. *J Bone Miner Res* 1998;13:1124-1133.
13. **Peris P et al.** Biochemical markers of bone turnover after surgical menopause and hormone replacement therapy. *Bone* 1999;25:349-353.
14. **Reginster J-Y et al.** Reduction in PINP, a marker of bone metabolism, with raloxifene treatment and its relationship with vertebral fracture risk. *Bone* 2004;34:344-351.
15. **Sharp CA et al.** Effects of low- and conventional-dose transcutaneous HRT over 2 years on bone metabolism in younger and older postmenopausal women. *Eur J Clin Invest* 1996;26:763-771.
16. **Suvanto-Luukkonen E et al.** Comparison of three serum assays for bone collagen formation during postmenopausal estrogen-progestin therapy. *Clin Chim Acta* 1997;266:105-116.
17. **Black DM et al.** The effects of parathyroid hormone and alendronate alone or in combination in postmenopausal osteoporosis. *N Engl J Med* 2003;349:1207-1215.
18. **Chen P et al.** Early changes in biochemical markers of bone formation predict BMD response to teriparatide in postmenopausal women with osteoporosis. *J Bone Miner Res* 2005;20:962-97.
19. **Deal C et al.** Combination teriparatide and raloxifene therapy for postmenopausal osteoporosis: results from a 6-month double-blind placebo-controlled trial. *J Bone Miner Res* 2005;20:1905-1911.

20. **Ettinger B et al.** Differential effect of teriparatide on BMD after treatment with raloxifene and alendronate. *J Bone Miner Res* 2004;19:745-751.
21. **Eastell R et al.** Development of an algorithm for using PINP to monitor treatment of patients with teriparatide. *Curr Med Res Opin* 2006;22:61-66.
22. **Alvarez L et al.** Relationship between biochemical markers of bone turnover and bone scintigraphic indices in assessment of Paget's disease activity. *Arthritis Rheum* 1997;40:461-468.
23. **Alvarez L et al.** Components of biological variation of biochemical markers of bone turnover in Paget's bone disease. *Bone* 2000;26:571-576.
24. **Pons F et al.** Quantitative evaluation of bone scintigraphy in the assessment of Paget's disease activity. *Nucl Med Comm* 1999;20:525-528.
25. **Diaz-Martin MA et al.** Aminoterminal propeptide of type I collagen and bone alkaline phosphatase in the study of bone metastases associated with prostatic carcinoma. *Scand J Clin Lab Invest* 1999;59:125-132.
26. **Koizumi M et al.** The serum level of the amino-terminal propeptide of type I procollagen is a sensitive marker for prostate cancer metastasis to bone. *BJU International* 2001;87:348-351.
27. **Jukkola A et al.** Postoperative PINP in serum reflects metastatic potential and poor survival in node-positive breast cancer. *Anticancer Res* 2001;21:2873-2876.
28. **Risteli J and Risteli L.** Products of bone collagen metabolism. In: Seibel MJ, Robins SP, Bilezikian JP (eds). *Dynamics of bone and cartilage metabolism*. Academic Press, San Diego, 1999. p 275-287.
29. **Melkko J et al.** Clearance of the aminoterminal propeptides of types I and III procollagen is a physiological function of the scavenger receptor in liver endothelial cells. *J Exp Med* 1994;79:405-412.
30. **Naylor KE and Eastell R.** Measurement of biochemical markers of bone formation. In: Seibel MJ, Robins SP, Bilezikian JP (eds). *Dynamics of bone and cartilage metabolism*. Academic Press, San Diego, 1999. p 401-410.
31. **IFCC, EPTRV:** The theory of reference values. Part 5. Statistical treatment of collected reference values. Determination of reference limits. *J Clin Chem Clin Biochem* 1983;21:749-760.
32. **Liberman U et al.** Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. *N Engl J Med* 1995;333:1437-1443.
33. **Stepan JJ and Vokrouhlicka J.** Comparison of biochemical markers of bone remodelling in the assessment of the effects of alendronate on bone in postmenopausal osteoporosis. *Clin Chim Acta* 1999;288:121-135.
34. **Bauer DC et al.** Change in bone turnover and hip, non-spine, and vertebral fracture in alendronate-treated women: The Fracture Intervention Trial. *J Bone Miner Res* 2004;19:1250-1258.
35. **Heikkinen J et al.** Long-term continuous combined hormone replacement therapy in the prevention of postmenopausal bone loss: A comparison of high- and low-dose estrogen-progestin regimens. *Osteoporos Int* 2000;11:929-937.
36. **Heikkinen J et al.** Optimizing continuous-combined hormone replacement therapy for postmenopausal women: A comparison of six different treatment regimens. *Am J Obstet Gynecol* 2000;182:560-567.
37. **McClung M et al.** Opposite bone remodelling effects of teriparatide and alendronate in increasing bone mass. *Arch Inter Med* 2005;165:1762-1768.
38. **Crofton PM et al.** Procollagen type I amino-terminal propeptide: Pediatric reference data and relationship with procollagen type I carboxy-terminal propeptide. *Clin Chem* 2004;11:2173-2176.

UniQ[®] is a registered trademark of Aidian Oy.

AIDIAN

Aidian Oy
Koivu-Mankkaan tie 6 B, FI-02200 Espoo, Finland
P.O. Box 83, FI-02101 Espoo, Finland
+358 10 309 3000, aidian.eu