# QuikRead go® CRP Verification Set

Cat. No. 145216

Rx only

## Intended use

The QuikRead go CRP Verification Set is designed to be used for calibration verification and for method validation of the QuikRead go CRP system. This assayed verification material is intended for use with the QuikRead go CRP test and the QuikRead go instrument. **For *in vitro* diagnostic use.**

## Summary and explanation of the test

The QuikRead go CRP Verification Set contains CRP material to verify QuikRead go CRP system calibration. Local or state regulations may require that calibration verification be performed at set intervals. Use of calibration verification material can be an important quality-improvement tool.

## Reagents

The QuikRead go CRP Verification Set contains three ready-to-use vials, with different concentrations of CRP: QuikRead go CRP Verification Set Vial 1, Vial 2, and Vial 3. The material is serum-based and of human origin, and it contains <0.1% sodium azide as a preservative and bovine aprotinin as a stabilizer.

Component name and origin	Contents	Stability	
		Opened at 2...8°C (36...46°F)	Unopened at 2...8°C (36...46°F)
Vial 1 ORIG HUM	1 mL	21 days	Until the expiration date marked on the label
Vial 2 ORIG HUM	1 mL	21 days	Until the expiration date marked on the label
Vial 3 ORIG HUM	1 mL	21 days	Until the expiration date marked on the label
Instructions for use			

Do not use any vial that has visible evidence of microbial growth.

## Warnings and precautions

- Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.
- For *in vitro* diagnostic use only.
- For prescription only.
- The reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react

with metal plumbing, forming explosive compounds. Azide buildup can be avoided via flushing with a large volume of water when one is disposing of the reagents.

- The QuikRead go CRP Verification Set should be used in accordance with federal, state, and/or local regulations or accreditation requirements, as applicable.
- Do not use the product after the expiration date marked on the individual vial label.
- Do not exceed the stability periods for opened reagents.
- The QuikRead go CRP Verification Set contains liquid material of human origin that has been tested and found negative for human immunodeficiency virus antibodies (HIV), hepatitis C virus (HCV), and hepatitis B surface antigen (HBsAg). As a precaution, however, handle the calibrators as if they were capable of transmitting infectious agents.

## Preparation

Allow the material to stabilize for 30 minutes at room temperature, 18...25 °C (64...77 °F). Mix well. Do not shake.

## Procedure

The materials provided with the QuikRead go CRP Verification Set are listed in the "Reagents" section. The following materials are required to perform the test but are not provided:

## Materials required but not provided

Material	Cat. No.
QuikRead go® CRP kit	145215
QuikRead go® instrument	145218

## Assay procedure

The verification set is to be used in accordance with the instructions accompanying the QuikRead go CRP test. Use the verification set material in the same way as unknown whole blood specimen would be used, following the instructions for the QuikRead go CRP test procedure.

The QuikRead go CRP Verification Set values are uniquely specified for a sample volume of 20 µL, with the same instructions for use as for a whole blood sample.

Instructions for use  
146380-2  
English

## Limitations of the procedure

The QuikRead go CRP Verification Set is to be used at room temperature: 18...25 °C (64...77 °F).

## Expected values

The verification set materials represent the low, middle, and high level of the reportable range. The expected values are printed on the vial labels. Each individual result should fall within the range indicated on the respective vial label. The target values have been derived from multiple analyses of representative samples of the product and are specific to each verification set lot.

## Disposal

Disposal of the verification set materials and their containers should be in accordance with federal, state, and local requirements. Discard all biohazardous waste in container designed as for hazardous biowaste. Dispose of capillary tubes in a sharps container labeled with the biohazard symbol.

## Trademarks

QuikRead go is a registered trademark of Aidian.

## Explanation of symbols


*In vitro* diagnostic medical device

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Aidian Oy