# QuikRead go® CRP Control Set

Cat. No. 145217

Rx only

Instructions for use  
146381-3  
English

## Intended use



The QuikRead go CRP Control Set is intended for use as assayed quality-control material for monitoring the performance of the quantitative QuikRead go CRP assay with the QuikRead go Instrument. **For *in vitro* diagnostic use.**

## Summary and explanation of the test

The QuikRead go CRP Control Set contains CRP controls that enable the users to monitor the QuikRead go CRP test performance. The controls should be handled in the same way as patient samples.

## Reagents

The QuikRead go CRP Control Set features two ready-to-use controls, with different concentrations of CRP: QuikRead go CRP Control Low and QuikRead go CRP Control High. They both are human-based material with added colour, and contain < 0.1% sodium azide as a preservative.

Component name and origin	Con-tents	Stability	
		In-use stability at 2...8°C (36...46°F)	Shelf-life at 2...8°C (36...46°F)
Control Low 	1 mL	1 month	Until the expiration date marked on the label
Control High 	1 mL	1 month	Until the expiration date marked on the label
Instructions for use			

Do not use any vial that has visible evidence of microbial growth.

## Warnings and precautions

- Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For *in vitro* diagnostic use only.
- For prescription only.
- The reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide buildup can be avoided by flushing with a large volume of water when disposing of the reagents.
- The QuikRead go CRP Control Set should be used in accordance with federal, state, and/or local regulations or accreditation requirements, as applicable.

- Do not use the product after the expiration date marked on the label of the individual vial.
- Do not exceed the stability periods for opened reagents.
- QuikRead go CRP controls contain material of human origin that has been tested and found negative for human immunodeficiency virus (HIV 1 and 2) antibodies, hepatitis C virus (HCV), and hepatitis B surface antigen (HBsAg). As a precaution, handle this product as if capable of transmitting infectious agents.

## Control preparation

Allow the control material to stabilize for 30 minutes at room temperature, 18...25 °C (64...77 °F). Mix well. Do not shake.

## Procedure

The materials provided with the QuikRead go CRP Control Set are listed in the section "Reagents." The following materials are required to perform the test but are not provided.

## Materials required but not provided

Material	Cat. no.
QuikRead go® CRP kit	145215
QuikRead go® Instrument	145218

## Assay procedure

The QuikRead go CRP Control Set is to be used in accordance with the instructions accompanying the QuikRead go CRP test. Use the control set's controls in the way in which an unknown whole blood specimen would be used, in accordance with the instructions for the assay procedure. The QuikRead go CRP control values are uniquely specified for a control sample volume of 20 µL, with the same instructions for use as for a whole blood sample.

## Limitations of the procedure

The QuikRead go CRP Control Set is to be used at room temperature: 18...25 °C (64...77 °F).

## Expected values

Each individual result should fall within the range indicated on the respective control vial label. Each laboratory should establish its own mean and precision data. The target concentrations and ranges specified on the control vials have been derived from multiple analyses with the QuikRead go Instrument and QuikRead go CRP tests.


## Disposal

Disposal of the controls and their containers should be in accordance with federal, state, and local requirements. Discard all biohazardous waste into a container for hazardous biowaste. Discard used capillaries in a sharps container labeled with the biohazard symbol.

## Trademarks

QuikRead go is a registered trademark of Aidian Oy.

## Explanation of symbols

	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Batch code
	Use-by date
	Temperature limit
	Manufacturer
	Contents
	Control
	Origin: human
	Contains sodium azide
Rx only	For prescription only

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QuikRead go es una marca registrada de Aidian Oy.