






146382-3

- English
- Español

# QuikRead go<sup>®</sup> CRP

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**8**

Measure	RESULT	
CRP:	20 mg/l	
Patient ID: XXXXXXXX	Measurement time: 2013-06-11 12:19	
Test: CRP	Result info	
<b>i</b> Choose <b>Result info</b> to view result information. Remove cuvette to perform a new measurement.		
Exit	Print	New measurement

## 1 Intended use

The QuikRead go<sup>®</sup> CRP test is an immunoturbidimetric assay for the in vitro quantitative determination of C-reactive protein (CRP) in K<sub>2</sub>-EDTA and lithium heparin whole blood, K<sub>2</sub>-EDTA and lithium heparin plasma, and in serum samples. The test is carried out by means of the QuikRead go instrument.

Measurement of C-reactive protein aids in the evaluation of injury to body tissues, and infection and inflammatory disorders. The instrument and assay are for use by trained professionals in the clinical laboratory. *For in vitro diagnostic use only.* Not for point-of-care use.

## 2 Summary and explanation of the test

Infections, tissue injuries, and inflammatory disorders may be accompanied by a non specific elevation of the CRP level in the patient's blood. The rise in CRP levels is rapid, and increased levels can be detected within 6–12 hours of the onset of the inflammatory process.<sup>1</sup> In acute inflammatory processes the CRP levels can reach approximately 20–500 mg/L.<sup>3,4</sup>

Quantitative measurement of CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections, as well as an effective tool in controlling and monitoring postoperative infections.<sup>1, 4-7</sup>

The advantages of measuring CRP lie in it being a sensitive, early, and reliable reflector of the acute-phase response. It can also be used for accurate monitoring of the patient's condition and response to therapy. However, the rise in CRP level is non-specific.

## 3 Principles of the procedure

QuikRead go CRP is an immunoturbidimetric test based on the agglutination reaction. Microparticles are coated with anti-human CRP F(ab)<sub>2</sub> fragments, and the CRP present in the sample reacts with the microparticles. The resultant change in the turbidity of the solution is measured by the QuikRead go instrument.

The sample is added to the buffer in prefilled cuvettes, and the blood cells in a whole blood sample are hemolyzed. The calibration data of the test are found in the barcode label of the cuvette, which the QuikRead go unit reads automatically before the start of the test. For whole blood sample, the CRP value is automatically corrected on the basis of the sample's hematocrit level. Hematocrit results are not displayed but are used in the calculations. The CRP measurement range of the assay is 5–150 mg/L with whole blood samples and 5–120 mg/L with serum or plasma samples.

## 4 Reagents

### QuikRead go CRP Cat. no. 145215

Component name and origin	Symbol	$\nabla$ <sub>50</sub>	Stability			
			Opened at 2...8 °C (36...46 °F)	Opened at 18...25 °C (64...77 °F)	Unopened at 2...8 °C (36...46 °F)	Unopened at 18...25 °C (64...77 °F)
CRP Reagent Caps <b>ORIG</b> <b>SHP</b>	<b>REAG</b> <b>CPS</b>	50	Until the expiration date marked on the label	24 h/day for 1 month, 7.5 h day for 3 months	Until the expiration date marked on the label	24 h/day for 1 month, 7.5 h/day for 3 months
Buffer	<b>BUF</b>	50 x 1 mL	Foil pouch: 6 months	Foil pouch: 3 months	Until the expiration date marked on the label	Until the expiration date marked on the label
			Single prefilled cuvette: 2 hours	Single prefilled cuvette: 2 hours	–	–
Capillaries (20 µL)	<b>CAPIL</b>	50	–	–	–	–
Plungers	<b>PLUN</b>	50	–	–	–	–

The reagents contain sodium azide as a preservative; please see Section 5, "Warnings and precautions."

### Reagent preparation and storage conditions

All reagents are ready to use.

Keep the CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

### Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

## 5 Warnings and precautions

### Health and safety information

- Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For *in vitro* diagnostic use only.
- For prescription only.
- For professional use only.
- Do not smoke, eat, or drink on premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash your hands thoroughly when finished performing the test.
- Avoid contact with the skin and eyes. After skin contact, wash immediately with plenty of water.
- All patient samples and controls should be handled as potentially infectious material.
- The lyophilized reagent contains <1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long-lasting effects (H412). Avoid release to the environment (P273). Dispose of contents in accordance with national and local law (P501). The reconstituted and liquid reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide buildup can be avoided via flushing with a large volume of water when one is disposing of the reagents.
- For disposal considerations, see Section 15.

### Analytical precautions


- Do not use the product after the expiration date marked on the outer package.
- Do not exceed the stated stability periods for opened reagents.
- The QuikRead go CRP kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers. The components are disposable; never reuse components already used for performing a test.
- When opening a kit for the first time, ensure that the foil pouches (2 pcs.) protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces of the lower part of the cuvette (the optical part). Discard any cuvettes with fingerprints.
- Keep the CRP Reagent Caps away from moisture. Close the aluminum tube immediately after taking out the required number of reagent caps.
- Do not splash liquid into the measuring well of the instrument.

## 6 Sample collection and preparation

### Sample material, collection and volume

#### Use of capillaries and plungers

The QuikRead go CRP kit contains pre-measured 20  $\mu\text{L}$  capillaries and plungers. The plunger is inserted on the blue stripe end. The capillary is then filled up to the white stopper. Take care to ensure that there is no air in the capillary and that the capillary is filled all the way to the white stopper. Use a soft tissue to wipe away any residual sample from the outside of the capillary. Make sure the soft tissue does not draw any sample from inside the capillary.

Sample material	Sample collection
Anticoagulated whole blood	Use a venous blood sample collected in a tube containing Li-heparin or $\text{K}_2\text{-EDTA}$ . Mix whole blood by inverting the tube several times, and take 20 $\mu\text{L}$ with a capillary. Only the capillaries included in the kit are recommended for adding of the sample.
Plasma	Use $\text{K}_2\text{-EDTA}$ /Li-heparin plasma. Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Only the capillaries included in the kit are recommended for adding of the sample.
Serum	Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Only the capillaries included in the kit are recommended for adding of the sample.

### Sample dilution

Sample type	Instructions
Whole blood	Do not dilute whole blood samples.
Plasma/serum	Plasma or serum samples may be diluted with a 0.9% NaCl solution before the sample is added to the cuvette. The recommended dilution ratio is 1+3 (by volume, 1 part sample + 3 parts 0.9% NaCl). Take 20 $\mu\text{L}$ of the diluted sample into a cuvette, analyze the sample, and multiply the obtained result by 4.

### Sample storage

Sample material	Storage
Anticoagulated whole blood ( $\text{K}_2\text{-EDTA}$ and Li-heparin)	2...8 $^{\circ}\text{C}$ (36...46 $^{\circ}\text{F}$ ) for 3 days
Plasma ( $\text{K}_2\text{-EDTA}$ and Li-heparin)	2...8 $^{\circ}\text{C}$ (36...46 $^{\circ}\text{F}$ ) for 7 days
Serum	2...8 $^{\circ}\text{C}$ (36...46 $^{\circ}\text{F}$ ) for 7 days
Sample ( $\text{K}_2\text{-EDTA}$ and Li-heparin whole blood, $\text{K}_2\text{-EDTA}$ and Li-heparin plasma, serum) in buffer	18...25 $^{\circ}\text{C}$ (64...77 $^{\circ}\text{F}$ ) for 2 hours

Allow the samples to reach room temperature (18...25  $^{\circ}\text{C}$ , 64...77  $^{\circ}\text{F}$ ) before testing.

## 7 Procedure

### Materials required but not provided

The following materials are required to perform the test but are not provided. The materials provided are listed in Section 4. "Reagents".

Material	Cat. No	Remarks
QuikRead go <sup>®</sup> instrument with software version 6.3.6 or newer	145218	–
QuikRead go <sup>®</sup> CRP Control Set	145217	Recommended for quality control
QuikRead go <sup>®</sup> CRP Verification Set	145216	Recommended for calibration verification

### Assay procedure

Before starting an assay, read and follow the QuikRead go instrument and QuikRead go CRP test instructions for use. Note that marketing materials that may be enclosed (including multimedia presentations) do not replace the written instructions for use; they merely provide an additional tool contributing to the understanding of the QuikRead concept.

Open the foil pouch of the cuvette rack, and mark the opening date on the cuvette rack.

### Sample collection (see figures 1–5)

Only the capillaries (20 µL, cat. no. 132153) and plungers (cat. no. 67966) included in the QuikRead go CRP kit (cat. no. 145215) are recommended for adding the sample to a prefilled cuvette.

1. A prefilled cuvette must reach room temperature (18...25 °C, 64...77 °F) before use. This takes 15 minutes for an individual refrigerated (2...8 °C, 36...46 °F) prefilled cuvette. Do not touch the clear flat surfaces of the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on the results. The test must be performed within two (2) hours from opening of the cuvette.
2. Insert the plunger into the capillary at the end with the blue stripe.
3. Fill the capillary with the sample (20 µL), up to the white stopper. Ensure that there are no air bubbles in the capillary. Use a soft cloth or tissue to wipe any residual sample from the outside of the capillary.
4. Place the sample capillary into the buffer solution in the cuvette, and dispense by pressing down the plunger. Make sure that the capillary is completely empty.
5. Close the cuvette tightly with a CRP Reagent Cap. Do not press down the blue inner part of the cap. After the sample has been added to the buffer, the solution is stable for two (2) hours. Keep the cuvette upright, and do not shake.

### Analyzing the sample (see figures 6–8)

Begin the test by turning on the instrument. For full instructions, see the QuikRead go instrument's instructions for use. The display will guide you through the test.

6. Select *Measure* on the display of the QuikRead go instrument.
7. Place the cuvette in the instrument, with the barcode facing you, as shown in Figure 7. The instrument detects the sample type. The display shows how the measurement is progressing.
8. The result appears on the display when the measurement is completed, and the cuvette automatically rises from the measurement well.

### Analyzing the control

Please see the instructions for use provided with the QuikRead go CRP Control Set (cat. no. 145217). Analyze the control samples as patient samples, but select Quality Control on the display of the QuikRead go instrument. The result will be stored as a control measurement.

## 8 Quality control

External controls may be used to verify that the reagents and procedures are performing properly. Use of QC materials is recommended for each day of testing, in accordance with federal, state, and/or local regulations or accreditation requirements. The QuikRead go CRP Control set is recommended to be used with the QuikRead go CRP test.

### QuikRead go CRP Control Set (cat. no. 145217):

The QuikRead go CRP Control Set contains two ready-to-use controls, with different levels of CRP. Perform the test as outlined in Section 7. "Procedure". The sample volume is 20 µL. The 20 µL capillaries are recommended for adding of the sample.

The QuikRead go CRP tests are factory-calibrated. The QuikRead go CRP Verification Set (cat. no. 145216) should be used for verification of the calibration in line with CLIA regulations. Please consult Section 16. "Troubleshooting" and contact Orion Diagnostica or your local distributor if you observe any changes in test functionality.

### Other commercially available controls:

Handle and use the controls in accordance with the instructions for each control substance. Acceptable limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells, because these might not hemolyze normally.

## 9 Results' interpretation

After the sample has been successfully analyzed, the CRP results are presented on the screen of the QuikRead go instrument. If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

The rise in CRP level is non-specific and should not be interpreted without a complete clinical history. Serial CRP measurements may be required in monitoring of patient condition or response to therapy.

Test values are considered to be clinically significant at levels above 10 mg/L (cutoff value)<sup>2, 9, 10</sup>.

If diluted plasma or serum samples have been used, the final results need to be calculated in accordance with the instructions in Section 6. "Sample Collection and Preparation".

## 10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

CRP results should never be used alone, without a complete clinical evaluation. Intra-individual variations in CRP are significant and should be taken into account – e.g., by means of serial measurements – in interpretation of the values.

The QuikRead go CRP test is not intended for use as high sensitivity CRP or cardiac CRP.

## 11 Expected values

The assay reference range limit was determined by using 143 apparently healthy adults (59 males and 84 females) with age 19-65 according to CLSI EP28-A3C guideline. Based on the results the 95 % reference limit was concluded to be  $\leq 5$  mg/L for each sample type. Results are based on a non parametric approach. However, each laboratory is recommended to establish a range of normal values for their population in their region.

## 12 Performance characteristics

### Precision

A precision study was performed in line with Clinical and Laboratory Standards Institute (CLSI) guidelines, EP05-A3.

Within-run, between-day, and total precision						
Sample material	Sample number	Number of days	Mean CRP (mg/L)	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Li-heparin whole blood	Sample 1	20	11	6.7	2.2	8.2
	Sample 2	20	21	4.0	2.3	4.7
	Sample 3	20	47	2.5	2.9	3.9
	Sample 4	20	146	2.5	1.8	3.3
K <sub>2</sub> -EDTA whole blood	Sample 1	20	8	7.7	2.2	9.2
	Sample 2	20	23	3.4	2.3	6.1
	Sample 3	20	53	1.3	2.9	2.6
	Sample 4	20	131	1.7	1.9	2.6
Li-heparin plasma	Sample 1	20	9	4.9	2.7	5.6
	Sample 2	20	20	2.8	2.3	3.9
	Sample 3	20	43	2.2	1.7	2.8
	Sample 4	20	104	1.8	0.0	2.2
K <sub>2</sub> -EDTA plasma	Sample 1	20	10	2.8	1.2	3.4
	Sample 2	20	19	1.7	1.3	2.6
	Sample 3	20	51	0.8	0.5	1.4
	Sample 4	20	107	1.5	1.2	2.0
Serum	Sample 1	20	9	4.0	3.0	5.0
	Sample 2	20	21	2.9	1.4	3.2
	Sample 3	20	51	1.0	0.1	1.5
	Sample 4	20	95	2.0	2.3	3.2
Control 1	Sample 1	20	25	2.7	2.1	3.4
Control 2	Sample 1	20	76	1.0	0.9	1.3

### Interference

Interferent	No interference found up to concentration
Bilirubin	230 mg/L
Vitamin C	35 mg/L
Triglycerides	10 g/L
Rheumatoid factors (RF)	525 IU/mL
Anticoagulants (Li-heparin or EDTA)	No interference
Hemoglobin (plasma and serum)	20 g/L

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test, because the assay antibodies lack the Fc part. In rare cases, interference by IgM myeloma protein has been observed.

### Limits of detection and quantitation

The assay limit of detection is 1.9 mg/L.

The assay limit of quantitation is 3.6 mg/L.

The lower limit of measuring range is 5 mg/L.

### Measurement range

#### Plasma and serum samples

The QuikRead go CRP test is linear with plasma and serum samples within the measurement range of 5–120 mg/L. When the CRP result is outside the measurement range, the result will be displayed as "CRP < 5 mg/L" or "CRP > 120 mg/L".

#### Whole blood samples

The QuikRead go CRP test is linear with whole blood samples (with hematocrit values between 15–75%) within the measurement range of 5–150 mg/L. When the CRP result is outside this measurement range, the result will be displayed as "CRP < 5 mg/L" or "CRP > 150 mg/L". The CRP result is not displayed if the level of hematocrit is outside the range 15–75%.

### Prozone

CRP concentrations of less than 1000 mg/L do not give falsely low results.

## 13 Traceability

The calibrators used to calibrate the QuikRead go CRP test are traceable to the ERM®-DA474 reference material.

## 14 Method comparison

Method comparison studies were performed with all sample matrices. The QuikRead go CRP test was compared against two clinical laboratory methods. The performance characteristics are in the table below:

Sample type	Clinical laboratory method	Number of samples	Intercept	Slope	R
K <sub>2</sub> -EDTA whole blood	ABX CRP REA	68	-2.98	0.90	0.992
K <sub>2</sub> -EDTA plasma	Tina-Quant C-reactive protein Gen 3.	62	-0.24	1.01	0.998
Li-heparin plasma	Tina-Quant C-reactive protein Gen 3.	62	-0.12	0.95	0.998
Serum	Tina-Quant C-reactive protein Gen 3.	62	-0.34	0.95	0.997

### Whole blood versus plasma comparability

In a comparison of 62 patient samples, K<sub>2</sub>-EDTA whole blood (y) and respective K<sub>2</sub>-EDTA plasma (x) yielded comparable results:  $y = 0.96x - 1.24$

## 15 Disposal

- Dispose of contents in accordance with federal, state, and local requirements.
- All patient samples and used cuvettes should be handled and disposed of as biohazardous waste and discarded in a container designated as for hazardous biowaste. Discard lancets and capillaries in a sharps container labeled with the biohazard symbol.
- Materials used in the components
  - Paper: Instructions for use
  - Cardboard: Kit box
  - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, plunger and capillary tubes
  - Glass: Mylar coated capillaries
  - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
  - Several (not to be recycled): Lids of reagent cap tubes
- When used in accordance with Good Laboratory Practice, good occupational hygiene, and the instructions for use, the reagents supplied should not present a health hazard.

## 16 Troubleshooting

### Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the QuikRead go instrument instructions for use.

Error messages	Corrective actions
Cuvette position not correct. Remove cuvette.	Remove the cuvette and reinsert it, in the right position.
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner blue part of the cap is not pressed down.
Measurement prohibited. Kit lot expired.	Discard the expired kit lot. Use a new one.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25 °C, 64...77 °F). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool to room temperature (18...25 °C, 64...77 °F). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed, or the sample may contain interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been a problem during reagent addition. Make sure that the cap is properly closed.
Test cancelled. Instrument failure.	Perform a new test. If this message appears often, contact your local distributor or call Customer Service.
No result. Unable to do hematocrit correction.	Test the same cuvette again. If the error message appears again, perform a new test. If the test cannot be completed at all, the sample may contain interfering substances or be otherwise compromised.
No result. Hematocrit too low.	If the level of hematocrit in a whole blood sample is below the hematocrit measurement range of the instrument (see Section 12.), no result is shown.
No result. Hematocrit too high.	If the level of hematocrit in whole blood samples is above the hematocrit measurement range of the instrument (see Section 12.), no result is shown.
Error code xx. Please restart the instrument.	Restart the instrument. If the error message appears again, contact your local distributor or call Customer Service and give the error code.
Error code xx. Please call customer service.	Contact your local distributor or call Customer Service, reporting the error code.

## References • Referencias

1. van Leeuwen MA & van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillieres Clin Rheumatol.* 1994; 8 (3): 531–52.
2. Weitkamp JH and Aschner JL. Diagnostic use of C-reactive protein (CRP) in Assessment of Neonatal Sepsis. *Neoreviews* 2005;6(11): 508-15.
3. Hind CRH, Pepys MB. The role of C-reactive protein (CRP) measurement in clinical practice. *Int. Med* 1984;5: 112-151.
4. Olaison L et al. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med.* 1997; 157 (8): 885–92.
5. Peltola H et al. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics.* 1997; 99 (6): 846–50.
6. Philip AGS & Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics.* 2000; 106 (1): 1–5.
7. Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, (eds). *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003; 2: 150–6.
8. Guidance for Industry and FDA staff: Review Criteria for Assessment of C-reactive Protein (CRP), High Sensitivity C-Reactive Protein (hsCRP) and Cardiac C-Reactive Protein (cCRP) assays. Document issued on September 22, 2005.
9. Myers et al. CDC/AHA Workshop on Markers of Inflammation and Cardiovascular Disease: Application to Clinical and Public Health Practice: laboratory tests available to assess inflammation—performance and standardization: a background paper. *Circulation* 2004;110(25): 545-549.

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QuikRead go® is a registered trademark of Aidian Oy.