CULTURA[®] M

D deutsch

GB english

F français

I italiano

S svenska

NL nederlands

CE

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Bedienungselemente

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OPERATING INSTRUCTIONS

Operating controls/parts

- Main switch with control lamp
- Temperature regulator

Positioning

The CULTURA can be used in the vertical or horizontal position (switch upside).

Picture legend

(see pictures on back cover)

- The CULTURA can be used in either vertical or horizontal position (switch upside)
- Main switch with control lamp
- Temperature regulator

Setting up

- remove from packaging and plug into the mains supply.
- place the shelf in the desired position and stand the thermometer on the shelf.
- switch on and check that the control lamp is lit.
- after 30 mins. the set temperature is reached (approx 35°C).

Temperature

The working temperature can be changed by rotating the regulating screw with a screw driver. Turning to the left decreases and to the right increases the temperature (15° turn = approx. 1°C).

Tips

Make sure that the door is properly closed, and do not open it more than 90 degrees. Avoid direct sun light shining on the door. Do not place plastic tubes directly on the floor of the CULTURA – always use the shelf provided.

Cleaning

The CULTURA can be easily cleaned with standard cleaning and disinfectant materials. Important: Fluids containing methanol, acetone, etc. should never be used for cleaning purposes.

Failure

The CULTURA is fitted with an overheating fuse for safety in the event of accidental overheating. This fuse cannot be reset or replaced. In case of a failure, please contact your CULTURA supplier.

Maintenance/Regular checks

The Cultura is maintenance free. Regular visual inspections should be carried out. Especially the power cable and the housing should be checked for damage. In case of damages unplug the Cultura and do not use any longer. The operating temperature should be constantly controlled and readjusted when necessary.

TECHNICAL SPECIFICATIONS

- Registered name:** CULTURA
- Article number:** 70700 = 220–240 V
70701, 70713 = 100–120 V
70702, 70712 = 12 V
- Temperature range:** 25–45°C
- Temperature accuracy:** +/- 1°C
- Available voltages:** (see label on equipment)
220–240 V = 30 W
100–120 V = 30 W
12 V = 20 W
- Frequency:** 50Hz / 60Hz
- Safetyclass:** II
- Construction:** Flame resistant PBDE-free PC & ABS
- Door:** PC, transparent
- Thermometer:** up to 60°C, in special container
- Base:** ABS
- Guarantee:** Two years from date of purchase
- External Dimensions:** L x H x D:
310 x 155 x 168 mm
- Internal Dimensions:** L x H x D:
220 x 120 x 150 mm
- Weight:** 1,1 kg

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GB CULTURA Guarantee

This guarantee covers the above product and lasts for two years from the date of purchase. We guarantee that in the event of a fault in the CULTURA, resulting from material or manufacturing defect, the CULTURA will be replaced free of charge within the period of the guarantee. The CULTURA should be returned in the original packaging. This guarantee does not cover faults caused by improper use. Please keep this guarantee in a safe place and return it with the CULTURA in the event of a claim. The CULTURA is constructed to reliably operate over many years. Please check, in the event of a breakdown, that the power supply is properly connected and that the fuses are in order.

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
2

3

ALMEDICA
MONITORING & SAFETY IN HYGIENE

Almedica AG
Hauptstrasse 76
3285 Galmiz
026 672 90 90
almedica.ch
office@almedica.ch