Hygicult® TPC

Intended use

Hygicult TPC slides are intended for rapid monitoring of micro-biological hygiene in different types of materials, both solid and liquid. The test can be performed on-site, or the slides can be used as convenient transport media for samples.

The slide is covered on both sides with Total Plate Count Agar which supports rapid growth of most common bacteria and fungi. It is impossible to give exact cut-off values for microbial counts indicating excessive contamination of the monitored material because "normal" levels vary in a wide range depending on the application. The main significance of the test is that elevation of total microbial counts can be detected. Normal levels must therefore be first established. It should also be emphasised that total microbial count is not a direct indicator of the safety of food. It is usually just a sign of shorter storage life or mishandling of the product.

Note: The limit values for microbial count in normal drinking water are too low to be detected by the Hygicult method.

Contents of the kit

Hygicult TPC	Cat. No. 68010
Test slides	10 pcs
Labels	10 pcs
Instructions for use	1 pc

Typical formulation

Total Plate Count Agar	
Tryptose	Tween 80
Yeast extract	Agar agar
Glucose D	Water
Lecithin	

Warnings and precautions

Do not use product beyond the expiry date marked on the kit.

Do not use the kit if you notice

- discoloration or dehydration of the growth medium
- detachment of the growth media from the plastic slide
- evidence of bacterial or fungal growth

Do not touch the growth because any colony growing on the slide may be pathogenic.

Storage

Store the kit at room temperature (18...25°C / 64...77°F) protected from draught, temperature fluctuations and light sources.Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

Sampling

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested. After sampling screw the slide tightly back into the tube.

Contact inoculation (Fig. 1a, 1b)

Solid surfaces can be tested by pressing each side of the slide firmly against the surface for three or four seconds. The slide should be held still during pressing. The hinged design offers ease of use.

Dipping (Fig. 2)

Fluid samples are tested by dipping the slide in the liquid for three or four seconds. Blot the last drops on absorbent paper.

Instructions for use • English

Swabbing (Fig. 3)

Semisolid materials or objects that are difficult to reach can be tested by carefully rolling a sterile swab over an area delimited using e.g. a frame. If the object is dry, the swab should first be moistened with sterile water. The moistened swab can also be used for obtaining samples from powders (e.g. spices) or viscous fluids. After swabbing the sample area, roll the swab gently over the agar surfaces of the slide from left to right and from bottom to top.

Incubation (Fig. 4)

Incubate the slide tightly enclosed in its tube

- at 35...37° C for one day or
- at 27...30° C for two days or
- at 18...25° C for up to five days.

Some of the most common yeasts and moulds do not grow at 35...37° C, which should be taken into account when choosing the incubation temperature. When the incubation time exceeds one day, it is advisable to also read the results at day 1, as swarming strains of *Proteus* and *Bacillus* are often easier to read after one day's incubation. Some slow-growing organisms may not yet be visible after one day's incubation.

Interpretation of results (Fig. 5)

Remove the slide from its tube after incubation and determine the microbial count (number of colony forming units, CFU) by comparing the density of growth on the slide with the model density chart. The following levels¹ can be considered as a rough basis for evaluating the degree of total microbial contamination of foodstuffs. The limit values for total microbial count are lower for cooked food than unprepared food.

	Good	< 10 ⁵ CFU/g
	Tolerable	10 ⁵ – 10 ⁷ CFU/g
	Poor	> 10 ⁷ CFU/g

Limit values for total microbial count on worktops and other surfaces are equally difficult to give since the criteria depend on the use. Nevertheless, if a surface after cleaning yields 5 colonies/cm², it can be considered improperly cleaned². This figure corresponds to approximately 45 colonies/side on Hygicult TPC. The total microbial count gives no indication as to the nature of the growth, and any investigation for the presence of pathogens should be carried out separately.

Limitations of the method

When used as a contact slide, Hygicult TPC equals the contact plate method in sensitivity³, whereas the dip and swab procedures have a detection limit of 1000 CFU/ml. The allowed total microbial concentration of normal drinking water is too low to be reliably detected using Hygicult TPC.

Results obtained with different inoculation systems should not be compared. Valid comparisons can only be made among results obtained using the same technique on the same type of material.

Disposal

- Dispose of contents according to national and local law.
- All used components should be handled and disposed of as potentially pathogenic material.
- Materials of the components:
 - Paper: Instructions for use, patient labels
 - Cardboard: Kit box
 - Plastic: Tubes, caps and dipslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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