
136262-10








- English
- German
- French
- Italian

QuikRead go[®] Strep A

Assay procedure • Testverfahren

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1 Intended use

QuikRead go® Strep A is a qualitative test for the detection of *Streptococcus pyogenes* (Strep A) in throat samples, when Strep A pharyngitis (tonsillitis) is suspected. The test is carried out using the QuikRead go® instrument. For *in vitro* diagnostic use only.

2 Summary and explanation of the test

A sore throat is usually due to an inflammatory response to viral or bacterial infection and a common reason for visiting a doctor^{1,2,3,4}. The most likely bacterial cause is *Streptococcus pyogenes*, which causes a sore throat and mostly affects school-age children and young adults^{2,3,5}. A Strep A sore throat usually lasts less than 10 days, and the patients are infectious during the period they have symptoms and for approximately one week thereafter^{1,3}. Antibiotic therapy reduces the duration of symptoms^{1,6}.

The symptoms of a Strep A sore throat do not differ significantly from those of a sore throat caused by other agents^{2,4}. Therefore, symptoms alone cannot be used for an accurate Strep A sore throat diagnosis, so clinical and epidemiological data should be supported by laboratory tests^{2,3,4}.

3 Principles of the procedure

QuikRead go Strep A is an immunoturbidimetric test based on microparticles coated with rabbit anti Strep A antibodies. Strep A antigens present in the sample react with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead go Instrument.

The throat sample is collected with a HydraFlock® sterile flocked collection device (later referred as HydraFlock flocked swab) and processed in a separate extraction tube. During extraction the bacteria contained in the sample break down and the antigen being tested is released into the swab. When the extraction is finished, the swab is placed in a prefilled cuvette and the antigen is released into the solution by swirling the swab vigorously in the solution. Colour changes that are in accordance with the instructions indicate successful extraction and neutralisation. Test calibration data is in the barcode label of the cuvette and read automatically by the QuikRead go instrument before the test is started.

4 Reagents

Content of the kit

Component name and origin	Symbol	QuikRead go® Strep A Cat. no. 135883 50 tests
Strep A Reagent Caps ORIG RAB	REAG CPS	2 x 25
Buffer in prefilled cuvettes	BUF	2 x 25 x 0.8 ml
Extraction Reagent 1	REAG EXTR 1	1 x 6 ml
Extraction Reagent 2	REAG EXTR 2	1 x 6 ml
Extraction Tubes	TUBE EXTR	50
Positive Control	CONTROL +	1 x 1 ml
Negative Control	CONTROL -	1 x 1 ml
HydraFlock sterile flocked collection devices (Puritan Medical Products Company, REF 25-3406-H)	SWAB	50
Instructions for use		

Reagents contain sodium azide.

Extraction Reagent 1 contains acetic acid and Extraction Reagent 2 contains sodium nitrite. See section 5 "Warnings and precautions".

Storage

Reagent	Storage temperature 2...8°C	Storage temperature 18...25°C
Reagent Caps, Extraction Reagents and Controls in unopened tubes/vials	Until the expiry date on the Reagent tube/vial	Until the expiry date on the Reagent tube/vial
Reagent Caps, Extraction Reagents and Controls after first opening	12 months	12 months
Prefilled cuvettes in an unopened foil pouch	Until the expiry date on the cuvette pouch	Until the expiry date on the cuvette pouch
Prefilled cuvettes after first opening of the foil pouch	6 months	3 months
Opened prefilled cuvette	2 hours	2 hours

Mark the date of opening the foil pouch on the cuvette rack.

Reagent preparation and storage conditions

All reagents are ready to use. Keep the Strep A Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

Reagent deterioration

Do not use a cuvette with visible dirt in the buffer.

5 Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- Do not breath spray.
- All patient samples and controls should be handled as potentially infectious material.
- Lyophilised reagent inside a reagent cap contains < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.

- Extraction Reagent 1 contains < 10% acetic acid. Safety data sheet available on request (EUH210).
- Extraction Reagent 2 contains < 25% sodium nitrite (Acute Tox. 4, Eye Irrit. 2). Harmful if swallowed (H302). Causes serious eye irritation (H319). Wash hands thoroughly after handling (P264). Do not eat, drink or smoke when using this product (P270). Wear eye protection/protective gloves/protective clothing (P280). If eye irritation persists: Get medical advice/attention (P337+P313). Dispose of contents according to national and local law (P501).
- Disposal, see Section 14.

Warning

Analytical precautions

- Do not mix reagents from different lot numbers.
- Use only HydraFlock flocked swabs provided with a QuikRead go Strep A reagent kit. Other swabs may affect test performance.
- When opening a reagent kit for the first time, make sure that the protective foil pouch of the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes in it. Before using a cuvette, always check that its protective foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- Protect the QuikRead go Strep A Reagent Caps from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

6 Sample collection and preparation

Sample material, sample collection and storage

Sample material	Secretion collected from the tonsils with a swab provided with the kit.
Sample collection	The sample is collected from the tonsils from both sides of the throat and/or from the back of the throat. Avoid touching other parts of the mouth.
Sample storage	The swab sample can be stored for 48 hours at 2...25°C before testing.

7 Procedure

Materials required but not supplied

	Cat. no.
QuikRead go [®] Instrument with software 4.1 or newer	133893, 149915

Assay procedure

Open the foil pouch protecting the cuvette rack and mark the rack with the date when the foil pouch was opened.

Sample collection and preparation (figures 1 to 8)

- 1 A prefilled cuvette must reach room temperature (18...25°C) before use. It will take 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on results. The test must be performed within two (2) hours from opening the cuvette.
- 2 Take a sterilised HydraFlock flocked swab from its cover. Ask the patient to open his/her mouth, press down the tongue with a wooden spatula and start sample collection.
- 3 Remember to take the sample from tonsils on both sides of the throat. Rolling the swab while collecting the sample ensures that the sample is representative.
- 4 Place the swab in the extraction tube provided with the kit.
- 5 Add two (2) drops of colourless Extraction Reagent 1 into the tube and then two (2) drops of red Extraction Reagent 2. The solution turns yellow-orange.
- 6 Swirl the swab in the solution for 30 seconds and then leave it in the solution for at least 90 seconds, but no longer than 15 minutes.
- 7 Lift the swab from the extraction tube into the prefilled cuvette, transferring as much liquid as possible, and stir vigorously. The solution turns red. This colour change indicates neutralisation of the extraction solution and the sample being transferred into the solution. Note: If a large volume of liquid remains in the extraction tube, you should pour it into the cuvette. See step 8.
- 8 Lift the swab out of the liquid, and press it slightly against the inner wall of the cuvette to release all of the liquid for measurement. Take care not to splash any solution on the outer surface of the cuvette. Discard the swab. Pour any extra liquid from the extraction tube into the cuvette. Do not proceed if the solution remains orange-yellow (indicating acidity; see the section "Warnings and precautions – sodium azide").

Sample analysis (figures 9 to 12)

For full instructions, see QuikRead go Instrument instructions for use. The display of the instrument will guide you through the test.

- 9 Close the cuvette tightly with a Strep A Reagent Cap. Do not press down the pink coloured inner part of the reagent cap. After the sample has been added to the buffer, the solution is stable for four (4) hours.
- 10 On the QuikRead go instrument display, select **Measure**.
- 11 Place the cuvette in the instrument. The barcode should be facing you, as shown in figure 11. The display shows how the measuring progresses. The instrument first measures the sample blank and then the Strep A concentration for 1–3 minutes.
- 12 The result appears on the display when the measurement is completed and the cuvette will raise up automatically from the measurement well.

Control analysis (figures 5 to 12)

Start control analysis by placing the swab in the extraction tube. Add two (2) drops of control onto the bottom of the extraction tube. Continue the process from step 5 of section "Sample collection and preparation".

Analyse the control samples in the same way as the other samples, but select **Quality control** on the instrument display. The result will be saved as a control measurement.

8 Quality control

Regular use of ready-to-use controls provided with the QuikRead go Strep A kit is recommended.

Commercial controls: Handle and use controls according to instructions of use. When liquid controls are used, please note that the varying absorbance capacity affects the amount of bacteria transferred to the swab. This may affect the result especially in samples close to cut-off value. The varying absorbance capacity does not affect patient samples which are not liquid. Swabs pretreated with controls may give inaccurate results due to properties of the swab used.

9 Results' interpretation

Strep A result	Interpretation of test result
Strep A Positive	Strep A antigen concentration in the sample corresponds to bacterial concentration of at least 7×10^4 CFU/swab.
Strep A Negative	The sample does not contain a measurable concentration of Strep A antigen.

When Strep A test results are interpreted, the patient's clinical symptoms and background information should always be taken into account.

10 Limitations of the procedure

Proper sample collection and careful following of assay procedures as specified in the instructions for use are necessary for reliable test results.

The test does not separate acutely infected individuals from asymptomatic carriers. Low bacterial concentrations in the sample may show a negative result.

11 Expected values

A reasonable estimate of the prevalence of Strep A infection in an adult population with acute pharyngitis is 5% to 10% and in a paediatric population with acute pharyngitis 20% to 25%³. In a healthy population, *Streptococcus pyogenes* should not be present. The possibility exists, however, that some people are asymptomatic carriers of this bacterium⁵.

12 Performance characteristics

Sensitivity and specificity

The clinical performance of QuikRead go Strep A was evaluated in a multi-centre study of samples from 279 patients. The Strep A throat sample was collected with the swab provided with the QuikRead go Strep A kit at six physician offices. All patients showed symptoms of pharyngitis.

The comparison throat sample was taken with another swab and used to inoculate the Streptocult culture medium. After inoculation, a bacitracin disc was applied to the Streptocult medium as described in the instructions for use. The Streptocult slide was then sent to a microbiology laboratory for further culturing and identification of bacteria. β -haemolytic streptococcus colonies growing on the slide were confirmed by means of commercially available typing reagents. The density of growth and the presence of β -haemolytic streptococci other than of group A were recorded also. The QuikRead go Strep A test was carried out on the swab in accordance with the instructions for use contained in the kit.

The test detects bacteria in amounts corresponding to 7×10^4 CFU/swab. Therefore, Streptocult cultures containing less than 10 group A streptococcus colonies were excluded when the sensitivity and the specificity of the test were calculated.

Streptocult cultures having 10 or more group A streptococcus colonies on the slide were considered positive.

Table 1.	Culture positive	Culture negative
QuikRead go Strep A positive	74	5
QuikRead go Strep A negative	15	177
Sensitivity	83% (73.7-90.2%)*	
Specificity	97% (93.7-99.1%)*	
Agreement	93%	
PPV	94%	
NPV	92%	

* 95% confidence interval

When the sensitivity of the QuikRead go Strep A test was reported on the basis of the density of group A streptococcus colonies on the slides, the following figures were obtained:

Table 2.	Sensitivity / Number of CFU			
Sensitivity according to the number of colonies	< 10	≥ 10	≥ 51	≥ 100
	12.5%	83.1%	91.0%	95.5%

Reproducibility

Reproducibility was tested through analysis of three samples, with different Strep A concentrations, daily over five days. One of the samples did not contain Strep A antigen (negative), one had low amounts of Strep A antigen (low positive), and one contained large amounts of Strep A antigen (high positive).

Table 3.	Reproducibility of QuikRead go Strep A										
	Day 1		Day 2		Day 3		Day 4		Day 5		
	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	Test 1	Test 2	
Negative	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.	neg.
Low positive	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.
High positive	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.	pos.

The results show high reproducibility.

Interfering substances

Cross-reactivity

Tested organism	Bacterial strain	Cross-reactivity
<i>Streptococcus B</i>	ATCC 12386	On basis of the test results it can be estimated that the organisms in question, possibly detected in the mouth, do not cause any cross-reactivity in the QuikRead go Strep A test.
<i>Streptococcus C</i>	ATCC 12388	
<i>Streptococcus F</i>	ATCC 12392	
<i>Streptococcus G</i>	ATCC 12394	
<i>Staphylococcus aureus (Cowan)</i>	ATCC 12598	
<i>Candida albicans</i>	ATCC 14053	
<i>Neisseria sicca</i>	ATCC 29259	
<i>Pseudomonas aeruginosa</i>	ATCC 27853	
<i>Haemophilus influenzae type B</i>	ATCC 9795	
<i>Streptococcus pneumonia</i>	ATCC 6303	
<i>Moraxella catarrhalis</i>	ATCC 25238	

Antigen excess

No antigen excesses have been detected during the testing of samples in which *Streptococcus pyogenes* levels are exceptionally high (10^9 CFU/ml).

Measurement range

The test detects bacteria in amounts corresponding to 7×10^4 CFU/swab.

13 Traceability

QuikRead go Strep A has been tested against reference material produced from *Streptococcus pyogenes* strain ATCC 19615.

14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, sampling devices, controls, used tubes, caps and cuvettes should be handled and disposed of as potentially infectious material.
- Materials of the components of the QuikRead go kits:
 - Paper: Instructions for use
 - Cardboard: Kit box including its inner parts
 - Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, extraction vials and tubes, plunger and capillary tubes, swabs
 - Glass: Capillaries
 - Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps
 - Not to be recycled: Lids of reagent cap tubes (several)
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

15 Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For further details of error messages, consult the QuikRead go Instrument instructions for use.




Error message	Corrective action
Cuvette position not correct. Remove cuvette.	Remove the cuvette and re-insert it into right position.
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner pink part of the cap is not pressed down.
Measurement prohibited. Kit lot expired.	Discard expired kit lot. Use a new one.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25°C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18...25°C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. Make sure that the cap is properly closed.
Test cancelled. Instrument failure.	Perform a new test. If this message appears often, contact your local distributor or call customer service.
Error code xx. Please restart the instrument.	Restart the instrument. If the error message appears again, contact your local distributor or call customer service and give the error code.
Error code xx. Please call customer service.	Contact your local distributor or call customer service and give the error code.

Unexpected negative or positive result

Possible reasons for unexpected results are listed in the table below.

Possible cause	Corrective action
Samples that have not been collected according to instructions may contain little or no Strep A bacteria and thus give an incorrect negative result.	Take a new sample. Make sure the sample is collected from the tonsils and/or back of the throat.
If the swab touches the cheeks, gums or tongue, interfering materials may be included in the sample giving an incorrect result.	Take a new sample. Make sure the sample is collected from the tonsils and/or back of the throat.
Using swabs other than those provided with the QuikRead go Strep A kit may give an incorrect result.	Take a new sample using a HydraFlock flocked swab provided with the QuikRead go Strep A kit.
Using other than QuikRead go Strep A extraction tubes provided with the kit may give an incorrect negative result.	Retest. Make sure that you use an extraction tube provided with the kit.
Mixing reagents coming from different lot numbers and/or swabs intended for different tests may cause an incorrect result.	Retest. Make sure that all reagents originate from the same reagent package.
Storing reagents at an unsuitable temperature may cause an incorrect result.	Retest. Make sure that the reagents have been stored according to instructions.
Dirt on the optical surface of the cuvette can cause an incorrect positive result.	Retest. Do not touch the clear flat surfaces on the lower part of the cuvette.
Keeping the sample for too long in the swab or in the solution may cause an incorrect result.	Retest. Follow the time limits given in the instructions.

Sample analysis • Probenanalyse

12	Measure	RESULT
	Strep A:	POSITIVE
Patient ID:	XXXXXXXX	Measurement time: 2011-11-11 12:19
Test:	Strep A	Result info
ⓘ Choose Result info to view result information. Remove cuvette to perform a new measurement.		
<input type="button" value="Exit"/> <input type="button" value="Print"/> <input type="button" value="New measurement"/>		

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Swab / Abstrichtupfer / Ecouvillons / Tamponi:

Puritan Medical Products Company
31 School Street, Guilford, Maine 04443-0149, USA

EC REP Emergo Europe
Arnhem, The Netherlands

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


12 Measure		RESULT
Strep A:		POSITIVE
Patient ID:	XXXXXXXX	Measurement time: 2011-11-11 12:19
Test:	Strep A	Result info
Choose Result info to view result information. Remove cuvette to perform a new measurement.		
Exit		Print
		New measurement

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1. Del Mar CB, Glasziou PP, Spinks AB. Antibiotics for sore throat (Cochrane review). In: The Cochrane library, Issue2, 2004. Chichester, UK: John Wiley & Sons Ltd.
2. Bisno AL, Peter GS, Kaplan EL. Diagnosis of Strep throat in adults: Are Criteria Really Good Enough? Clin Infect Dis 2002; 35:126-9.
3. Ebel MH, Smith MA, Barry HC, Ives K, Carey M. The rational clinical examination. Does This Patient Have Strep A Throat? JAMA 2000; 284:2912-18.
4. Bisno AL, Gerber MA, Gwaltney JM Jr, Kaplan EL, Schwartz RH. Practise Guidelines for the Diagnosis and Management of Group A Streptococcal Pharyngitis. CID 2002; 35:113-25.
5. Cunningham MW. Pathogenesis of group A Streptococcal Infections. Clin microbial Rev 2000; 13(3):470-511.
6. Snellman LW, Stang HJ, Stang JM, Johnson DR, Kaplan EL. Duration of Positive Throat Cultures for Group A Streptococci After initiation of Antibiotic Therapy. Pediatrics, 1993; 91:1166-70.
7. Data on file, Aidian Oy.

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