Revision History/ Revisionshistorie/ Historique des révisions/ Cronologia delle revisioni

REF 140068 QuikRead go CRP+Hb

Material number and release version/ Materialnummer und Version und Veröffentlichungsdatum/ Numéro et version du document et date de publication/ Numero materiale e versione e data di rilascio	Change classification/ Klassifizierung ändern/ Classification de la modification/ Variazione della classificazione	Section of change/ Abschnitt der Änderung/ Section de la modification/ Sezione della variazione
140119-11, 06/2026	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	1. Intended purpose/ Verwendungszweck/ Utilisation prévue/ Finalità d'uso
		2. Summary and explanation of the analyte/ Zusammenfassung und Beschreibung des Analyten/ Résumé et explication de l'analyte/ Riepilogo e spiegazione dell'analita
		3. Principles of the procedure/ Verfahrensprinzipien/ Principe de la procédure/ Principi della procedura
		4. Reagents/ Reagenzien/ Réactifs/ Reagenti
		5. Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni
		6. Sample collection and preparation/ Probennahme und -vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni
		7. Procedure/ Verfahren/ Procédure/ Procedura
		8. Quality control/ Qualitätskontrolle/ Contrôle de qualité/ Controllo di qualità
		9. Interpretation of results/ Interpretation der Testergebnisse/ Interprétation des résultats/ Interpretazione dei risultati
		12. Performance characteristics/ Leistungsparameter/ Performances/ Caratteristiche delle prestazioni
		13. Traceability/ Nachweisbarkeit/ Traçabilité/ Tracciabilità
		14. Method comparison/ Methodenvergleich/ Comparaison de méthodes/ Comparaison de méthodes
		Explanation of symbols/ Zeichenerklärung/ Explication des symboles/ Spiegazione dei simboli
		Revision History/ Revisionshistorie/ Historique des révisions/ Cronologia delle revisioni
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Assay procedure/ Testverfahren/ Réalisation du test/ Procedura del test
		1. Intended purpose/ Verwendungszweck/ Utilisation prévue/ Finalità d'uso

Material number and release version/ Materialnummer und Version und Veröffentlichungsdatum/ Numéro et version du document et date de publication/ Numero materiale e versione e data di rilascio	Change classification/ Klassifizierung ändern/ Classification de la modification/ Variazione della classificazione	Section of change/ Abschnitt der Änderung/ Section de la modification/ Sezione della variazione
		2. Summary and explanation of the analyte/ Zusammenfassung und Beschreibung des Analyten/ Résumé et explication de l'analyte/ Riepilogo e spiegazione dell'analita
		3. Principles of the procedure/ Verfahrensprinzipien/ Principe de la procédure/ Principi della procedura
		4. Reagents/ Reagenzien/ Réactifs/ Reagenti
		5. Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni
		6. Sample collection and preparation/ Probennahme und -vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni
		7. Procedure/ Verfahren/ Procédure/ Procedura
		8. Quality control/ Qualitätskontrolle/ Contrôle de qualité/ Controllo di qualità
		9. Interpretation of results/ Interpretation der Testergebnisse/ Interprétation des résultats/ Interpretazione dei risultati
		11. Expected values/ Erwartete Ergebnisse/ Valeurs attendues/ Valori attesi
		12. Performance characteristics/ Leistungsparameter/ Performances/ Caratteristiche delle prestazioni
		13. Traceability/ Nachweisbarkeit/ Traçabilité/ Tracciabilità
		14. Method comparison/ Methodenvergleich/ Comparaison de méthodes/ Comparaison de méthodes
		16. Troubleshooting/ Fehlerdiagnose/ Résolution des problèmes/ Risoluzione dei Problemi
	General change/ Allgemeine Änderungen/ Modification générale	Assay procedure/ Testverfahren/ Réalisation du test/ Procedura del test
		3. Principles of the procedure/ Verfahrensprinzipien/ Principe de la procédure/ Principi della procedura
		4. Reagents/ Reagenzien/ Réactifs/ Reagenti
		5. Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni
		6. Sample collection and preparation/ Probennahme und -vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni
		7. Procedure/ Verfahren/ Procédure/ Procedura


Material number and release version/ Materialnummer und Version und Veröffentlichungsdatum/ Numéro et version du document et date de publication/ Numero materiale e versione e data di rilascio	Change classification/ Klassifizierung ändern/ Classification de la modification/ Variazione della classificazione	Section of change/ Abschnitt der Änderung/ Section de la modification/ Sezione della variazione
		8. Quality control/ Qualitätskontrolle/ Contrôle de qualité/ Controllo di qualità 9. Interpretation of results/ Interpretation der Testergebnisse/ Interprétation des résultats/ Interpretazione dei risultati 10. Limitations of the procedure/ Einschränkungen des Verfahrens/ Limites de la procédure/ Limitazioni della procedura 11. Expected values/ Erwartete Ergebnisse/ Valeurs attendues/ Valori attesi 12. Performance characteristics/ Leistungsparameter/ Performances/ Caratteristiche delle prestazioni 15. Disposal/ Entsorgung/ Mise au rebut/ Smaltimento 16. Troubleshooting/ Fehlerdiagnose/ Résolution des problèmes/ Risoluzione dei Problemi References/ Referenzen/ Références/ Riferimenti bibliografici
	Correction/ Korrektur/ Correction/ Correzione	5. Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni 7. Procedure/ Verfahren/ Procédure/ Procedura 12. Performance characteristics/ Leistungsparameter/ Performances/ Caratteristiche delle prestazioni 14. Method comparison/ Methodenvergleich/ Comparaison de méthodes/ Comparaison de méthodes 15. Disposal/ Entsorgung/ Mise au rebut/ Smaltimento 16. Troubleshooting/ Fehlerdiagnose/ Résolution des problèmes/ Risoluzione dei Problemi

Classification/ Klassifizierung/ Classification/ Classificazione	Description/ Beschreibung/ Description/ Descrizione
New instructions for use/ Neue Anwendungshinweise/ Nouveau mode d'emploi/ Nuove istruzioni per l'uso	First version of a new instructions for use/ Erste Version der neuen Anwendungshinweise/ Première version d'un nouveau mode d'emploi/ Prima versione delle nuove istruzioni per l'uso/
Correction/ Korrektur/ Correction/ Correzione	Correction of documentation anomalies (minor typographical, grammar, and formatting changes are not included)/ Korrektur von Anomalien in der Dokumentierung (kleinere typografische, grammatikalische und Formatierungsänderungen sind nicht eingeschlossen)/ Correction d'anomalies

	dans la documentation (les modifications typographiques, grammaticales et de formatage mineures ne sont pas incluses)/ Correzione delle anomalie della documentazione (i piccoli refusi, la grammatica e le variazioni di formattazione non sono inclusi)
Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Addition, revision and/or removal of information related to the product/ Hinzufügung, Änderung und/oder Entfernung von Informationen über das Produkt/ Ajout, révision et/ou suppression d'informations liées au produit/ Aggiunta, revisione e/o rimozione delle informazioni relative al prodotto
Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Changes due to regulations, directives, laws and/or registration/ Änderungen aufgrund von Verordnungen, Richtlinien, Gesetzen und/oder Registrierung/ Modifications dues aux réglementations, directives, lois et/ou à l'enregistrement/ Variazioni dovute a normative, direttive, leggi e/o registrazione
General change/ Allgemeine Änderungen/ Modification générale/ Variazione generale	Implementation of non-technical changes/ Implementierung von nichttechnischen Änderungen/ Mise en œuvre de modifications non technique/ Implementazione di variazioni non tecniche

140119-11

- English
- Deutsch
- Français
- Italiano

QuikRead go
Instrument


QuikRead go Plus
Instrument

QuikRead go


CRP+Hb

Assay procedure • Testverfahren


1


2


3


4


5


6

7

8

CRP 20 mg/l
Hb 125 g/l

QuikRead go CRP+Hb

Cat. no. 140068

English

To be used with the QuikRead go Instrument and the QuikRead go Plus Instrument. Hereafter referred to as the QuikRead go instrument.

1

Intended purpose

The QuikRead go CRP+Hb test is an immunoturbidimetric assay for the quantitative determination of C-reactive protein (CRP) values from capillary blood, anticoagulated (EDTA or heparin) venous whole blood or plasma, or serum sample, and a spectrophotometric assay for the quantitative determination of hemoglobin (Hb) from capillary blood or anticoagulated (EDTA or heparin) venous whole blood sample. The test is carried out using an automated instrument belonging to the QuikRead go system. The assay is intended to be used by healthcare professionals in the clinical laboratory and near-patient testing (NPT) setting (also known as point of care setting). CRP is an acute phase protein present in low concentrations in the blood of healthy individuals. Infections, tissue injuries and inflammatory disorders may cause a non-specific elevation of the CRP level in the patient's blood. Quantitative CRP measurement is a sensitive and reliable indicator of the acute phase response. The measurement of CRP aids in the diagnosis, evaluation, and monitoring of infection, tissue injury, inflammatory disorders, and associated diseases.^{1,2} Hb is the iron-containing protein in red blood cells that carries oxygen. Determining its concentration is a common procedure in both primary healthcare and acute care. Hemoglobin measurement is used as an aid in the diagnosis and screening of anemia. For *in vitro* diagnostic use only.

2

Summary and explanation of the analyte

C-reactive protein (CRP)

Infections, tissue injuries and inflammatory disorders may be accompanied by a non-specific elevation of the CRP level in the patient's blood. The rise of CRP concentration is rapid, and increased levels can be detected within 6 hours of the onset of the inflammatory process.^{1,2} Quantitative measurement of CRP NPT decreases antibiotic prescribing for approximately 25% in different acute respiratory infections and therefore can be used to guide antibiotic treatment.³ CRP has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections⁴⁻⁷. CRP is also suitable for the evaluation of neonatal infections^{8,9} and it is an effective tool in ruling out postoperative infections^{10,11}. In exacerbations of chronic obstructive pulmonary disease (COPD), CRP is used to guide antibiotic prescribing¹². In systemic inflammation, such as rheumatoid arthritis and inflammatory bowel disease, CRP measurement is routinely used as part of the initial assessment of the disease and for monitoring of disease activity, progress, and treatment response^{13,14}.

Hemoglobin (Hb)

Hemoglobin is the iron-containing protein in red blood cells that carries oxygen. Determining its concentration is a common procedure in both primary healthcare and acute care. A decrease in the hemoglobin concentration may follow for example an iron deficiency or loss of blood.¹⁵ Quantitative measurement of hemoglobin is used as an aid in diagnosis of anemia.

3

Principles of the procedure

CRP

The CRP measurement is immunoturbidimetric and it is based on agglutination reaction. Microparticles are coated with anti-human CRP F(ab)₂ fragments, and the CRP present in the sample reacts with microparticles. The resultant change in the turbidity of the solution is automatically measured by the QuikRead go instrument. The CRP assay correlates well with results obtained from several other immunoturbidimetric methods.

The sample is added to the buffer in a prefilled cuvette, and the blood cells in a whole blood sample are hemolyzed. The calibration data of the test are found on the barcode label of the cuvette, which the instrument reads automatically before the start of the test. The CRP value is automatically calculated by the instrument. When whole blood sample is used, the CRP value is automatically corrected according to the sample's hematocrit level. Hematocrit results are not displayed but are used in calculations. The measuring range of the CRP assay is 5–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 5–120 mg/l with serum or plasma samples.

Hb

The Hb measurement is based on measuring oxyhemoglobin photometrically at two wavelengths. The sample is added to the buffer in prefilled cuvette, and the blood cells in a whole blood sample are hemolyzed, allowing the absorption of released Hb to be measured. The calibration data of the test are found on the barcode label of the cuvette, which the instrument reads automatically before the start of the test. The Hb value is automatically calculated by the instrument. The instrument corrects for the effect of lipids or other turbidity-inducing factors on the results. The measuring range of the Hb assay is 50–245 g/l (5.0–24.5 g/dl, 3.1–15.2 mmol/l).

4

Reagents

Content of the kit

Component name and origin	Symbol	QuikRead go CRP+Hb Cat. no. 140068, 50 tests
CRP Reagent Caps ORIG SHP	REAG CPS	2 x 25
QuikRead go CRP+Hb Buffer (in prefilled cuvettes)	BUF	2 x 25 x 1 ml
QuikRead Capillaries (20 µl)	CAPIL HEP	50
QuikRead Plungers	PLUN	50
Instructions for use		

The reagents contain preservatives, please see Section 5. "Warnings and precautions".

Storage

The shelf-life of the unopened QuikRead go CRP+Hb kit at 2–8°C is until the expiry date marked on the kit box.

Component	Shelf-life at 2–8°C	Shelf-life at 18–25°C	In-use stability at 2–8°C	In-use stability at 18–25°C
Reagent caps	Until the expiry date of the kit in an unopened reagent cap tube	24 h per day for 1 month or 7.5 h per day for 3 months in an unopened reagent cap tube	Until the expiry date of the kit in a closed reagent cap tube	24 h per day for 1 month or 7.5 h per day for 3 months in a closed reagent cap tube
Prefilled cuvettes	Until the expiry date of the kit in an intact foil pouch		6 months	3 months
			Prefilled cuvette without foil cover: 2 hours	
Capillaries	Until the expiry date of the capillary tube			
Plungers	Until the expiry date of the plunger tube			

Reagent preparation and storage conditions

All reagents are ready to be used. Let the reagents reach room temperature. Keep the reagent caps away from moisture. Close the tube immediately after taking out the required number of reagent caps.

Reagent deterioration

Do not use a cuvette with visible dirt in the buffer.

5

Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- QuikRead Capillaries are intended for non-invasive sampling without skin contact.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- All patient samples and controls should be handled as potentially infectious material.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/ protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagent inside a reagent cap contains < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032). Contains synthetic polymer microparticles. Do not release the cap reagent outside the cuvette. Microparticles are contained within the closed cuvette when used according to instructions.
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration.
- Disposal: see Section 15.
- In case of a serious incident, please report it to the manufacturer or its representative and national authority. Summary of safety and performance is available on Eudamed.

Warning

Analytical precautions

- The QuikRead go CRP+Hb kit is intended solely to be used with the QuikRead go instrument.
- Do not use the product after the expiry date marked on the outer package.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- Do not exceed the stability periods for opened reagents.
- Hold the cuvette from the sides. Do not touch the data matrix or the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- CRP Reagent Caps are colour-coded blue to distinguish them from other QuikRead analytes.
- Keep the reagent caps away from moisture. Close the tube immediately after taking out the required number of reagent caps. Do not use reagent caps if they have not been stored according to the instructions.
- Do not attach labels to the cuvette. If needed, the identification information can be written on the opposite side of the label using a permanent pen.
- Do not splash liquid into the measuring well of the instrument.

6

Sample collection and preparation

Sample material, sample volume and collection

CRP results can be obtained from whole blood, plasma and serum samples, whereas Hb results can only be obtained from whole blood samples. QuikRead Capillaries (Cat. no. 67962) included in the QuikRead go CRP+Hb kit (Cat. no. 140068) are recommended for adding the sample.

Sample material	Sample volume	Obtainable result	Sample collection
Capillary blood	20 µl	CRP and Hb	Prick the clean and dry finger/heel with a lancet. Wipe off the first drop and collect 20 µl of blood from the second drop into a heparinized capillary.
Anticoagulated whole blood	20 µl	CRP and Hb	Use a venous blood sample collected in a tube containing heparin or EDTA. Mix whole blood by inverting the tube several times and take 20 µl with a capillary.
Plasma	20 µl or 12 µl *	CRP	Use EDTA/heparin plasma. Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Take 20 µl with a capillary.
Serum	20 µl or 12 µl *	CRP	Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Take 20 µl with a capillary.

*It is possible to pipette a sample volume of 12 µl in order to have a wider measuring range. In that case, you have to change the test parameter settings of the instrument. Please see the instrument instructions for use.

Sample dilution

QuikRead go CRP+Hb test does not require sample dilution. However, if the displayed CRP result is marked with a greater than (>) symbol, you can dilute the sample to obtain a result within the measuring range.

Sample material	Instructions
Whole blood	Do not dilute whole blood samples.
Plasma/ serum	Plasma or serum samples may be diluted with a 0.9% NaCl solution before the sample is added to the cuvette. The recommended dilution ratio is 1+3 (by volume, 1 part sample + 3 parts 0.9% NaCl). Take 20 µl of the diluted sample into a cuvette, analyse the sample and multiply the obtained result by 4.

Sample storage

Sample material	Short-term storage	Long-term storage
Capillary blood	Max. 15 minutes in the heparinized capillary	Do not store
Anticoagulated whole blood	2–8 °C for 3 days	Separate plasma and store as indicated below
Plasma	2–8 °C for 7 days	Below –20°C for 3 years ¹⁶
Serum	2–8 °C for 7 days	Below –20°C for 3 years ¹⁶
Diluted plasma/serum in 0.9% NaCl	18–25 °C for 30 minutes	Do not store
Patient sample in a cuvette	18–25 °C for 2 hours	Do not store

Allow the samples to reach room temperature (18–25 °C) before testing. Frozen samples must be thawed completely, mixed thoroughly, and brought to room temperature prior to testing. If the thawed plasma or serum samples contain clots, the samples should be centrifuged. Samples should not be repeatedly frozen and thawed.

7

Procedure

Materials required but not provided

The following materials are required to perform the test but are not provided or they are available also separately. The materials provided are listed in Section 4. "Reagents".

Material	Cat. no.	Additional info
QuikRead go Instrument with software version 4.1 or newer	133893 149915	QuikRead go instruments are equivalent instruments
QuikRead go Plus Instrument with software version 2.0.0 or newer	155378 155377	
QuikRead Capillaries (20 µl) 50 pcs	67962	Supplied with 140068
QuikRead Plungers 50 pcs	67966	
Pipette	–	Capillaries are recommended for adding the sample
QuikRead go CRP Control	153764	Recommended for quality control
QuikRead go CRP Control High	153763	
QuikRead go Hb Control	141154	
Fingertip lancets	–	–

Assay procedure

Before starting an assay, read and follow the instrument and QuikRead go CRP+Hb test instructions for use. Open the foil pouch of the cuvette rack. Remember to mark the opening date on the cuvette rack and the reagent cap tube.

Sample collection (see figures 1–5)

Glass capillaries (20 µl) and plungers are recommended for adding the sample into a prefilled cuvette. QuikRead Capillaries (Cat. no. 67962) and QuikRead Plungers (Cat. no. 67966) can be purchased separately, whereas the QuikRead go CRP+Hb kit (Cat. no. 140068) contains them already.

1. A prefilled cuvette must reach room temperature (18–25 °C) before use. This will take 15 minutes for an individual refrigerated (2–8 °C) prefilled cuvette. Do not touch the data matrix or the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on results. The test must be performed within two (2) hours from opening the cuvette.
2. Drop the plunger into the capillary at the end with the blue stripe.
3. Fill the capillary with the sample up to the white stopper (20 µl). Ensure that there are no air bubbles in the capillary. Use a soft cloth or tissue to wipe away any residual sample from the outside of the capillary.
4. Place the sample capillary into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.
5. Close the cuvette tightly with a reagent cap. Do not press down the blue inner part of the reagent cap. Close the reagent cap tube immediately after taking out the required number of reagent caps. Once the sample has been added to the buffer, the assay should be run within two (2) hours. Keep the cuvette in an upright position and do not shake.

Note! If a venous sample is used, it is recommended to take the sample directly from the tube. Tilt the tube slightly and collect the sample.

Analysing the sample (see figures 6–8)

Begin the test by turning on the QuikRead go instrument. For full instructions, see instrument instructions for use. The screen will guide you through the test.

6. Choose **Measure / Measure sample** on the screen of the QuikRead go instrument.
7. Place the cuvette into the instrument with the barcode facing you, as shown in Figure 7. The instrument detects the sample type. The screen shows how the measurement is progressing. The instrument first measures the sample blank and then the CRP and Hb concentrations within one minute.
8. The result is displayed on the screen, and the cuvette rises automatically when the measurement is completed.

Analysing the control (Fig. 7–8)

Please see the instructions for use provided with each control vial. Analyse the control samples as patient samples, but choose **Quality Control** on the screen of the QuikRead go instrument. Select the control that matches the test you are using (CRP or Hb). The result will be stored in a separate result file for quality control results.

8

Quality control

Quality control testing must be performed in accordance with applicable national or local regulations or accreditation requirements and laboratory standard Quality Control Procedures. Regular use of QuikRead go CRP Control (Cat. no. 153764) and/or QuikRead go CRP Control High (Cat. no. 153763), and QuikRead go Hb Control (Cat. no. 141154) is recommended.

It is recommended to perform the control at least with each new shipment and new reagent lot and every time an unexpected result is obtained. The control should be used each time new QuikRead go system users are trained. Please consult Section 16. "Troubleshooting" and contact Aidian or your local distributor if you observe any changes in test functionality.

QuikRead go CRP controls: All QuikRead go CRP controls are ready-to-use. The control values have been determined for the QuikRead go instrument, and the control is applied using the same procedure as for a plasma/serum sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (20 µl). If you wish to use sample volume 12 µl, change the test parameter settings of the instrument. Please see the instrument instructions for use. The control is stable for 15 minutes in the cuvette.

QuikRead go Hb control: QuikRead go Hb Control is ready-to-use. The control values have been determined for the QuikRead go instrument and the control is applied using the same procedure as for a whole blood sample. Perform the test as outlined in Section 7. "Procedure", use the same control sample volume as for unknown blood specimens (20 µl). The control is stable for 5 minutes in the capillary and 2 hours in the cuvette.

Other commercially available CRP and Hb controls: Handle and use the controls in accordance with the instructions for each control substance. Perform the test as outlined in Section 7. "Procedure". Acceptance limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells, because these might not hemolyze normally. Controls in which the hemoglobin is present in a form other than oxyhemoglobin, or which contain compounds that may modify Hb, might yield results that are not aligned with those obtained by other methods.

9

Interpretation of results

CRP

Elevation of CRP value is unspecific, and the results should be interpreted in the light of other clinical findings. Serial CRP measurements may be required to monitor a patient's condition or response to therapy and when measuring CRP levels in newborns. Additionally, establishing the patient-specific baseline CRP level is recommended when assessing slightly elevated CRP levels.

CRP test result	Interpretation of result
< 10 mg/l	• Excludes many acute inflammatory diseases but does not specifically exclude inflammatory processes.¹⁷ • Possible inflammations in newborn babies are not necessarily excluded.^{9,18,19} • Possible inflammations or exacerbations in conditions with systemic inflammation are not necessarily excluded.^{20–23}
10–50 mg/l	Elevated concentrations are found in acute disease that occur in the presence of a slight to moderate inflammatory process. ¹⁷
> 50 mg/l	Indicates high and extensive inflammatory activity. ¹⁷

If diluted plasma or serum samples have been used, the final results must be calculated in accordance with the instructions in Section 6. "Sample collection and preparation". If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

Hb

Hb test result		Interpretation of result ²⁴
Women	< 120 g/l (12.0 g/dl, 7.4 mmol/l)	Below reference value
	> 150 g/l (15.0 g/dl, 9.3 mmol/l)	Above reference value
Men	< 130 g/l (13.0 g/dl, 8.1 mmol/l)	Below reference value
	> 170 g/l (17.0 g/dl, 10.6 mmol/l)	Above reference value
Children 1–12 years*	< 111 g/l (11.1 g/dl, 6.9 mmol/l)	Below reference value
	> 155 g/l (15.5 g/dl, 9.6 mmol/l)	Above reference value

* Children's Hb values increase progressively gradually reaching adult levels. Hb values for children < 1 year vary.

For Hb reference values, see Section 11. "Expected values". Hb result levels are affected by patients' diets and the elevation from sea level of their place of residence²⁵. It is therefore recommended that each laboratory determines its own reference values for Hb. If necessary, please consult Section 16. "Troubleshooting", and repeat the test.

10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance characteristics".

CRP and/or Hb test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

11 Expected values

CRP

Reference value: < 3 mg/l^{23,26,27}, median value in apparently healthy adults: 0.8 mg/l²³. CRP levels in newborns can be affected by e.g. prematurity, duration of delivery and complications related to the delivery.

Assay reference range limit: < 5 mg/l. The assay reference range limit was determined by using 143 apparently healthy adults (59 males and 84 females) with age 19–65 according to the Clinical and Laboratory Standards Institute (CLSI) EP28-A3C guideline. Results are based on a non parametric approach. However, each laboratory is recommended to establish a range of normal values for their population in their region.

Hb

Reference group	Reference values ²⁴
Women	120–150 g/l (12.0–15.0 g/dl, 7.4–9.3 mmol/l)
Men	130–170 g/l (13.0–17.0 g/dl, 8.1–10.6 mmol/l)
Children 1–12 years*	111–155 g/l (11.1–15.5 g/dl, 6.9–9.6 mmol/l)

*Children's Hb values increase progressively reaching gradually adult levels.

However, each laboratory is recommended to establish a range of CRP and Hb normal values for their population in their region.

12 Performance characteristics

The performance characteristics of the QuikRead go CRP+Hb test have been demonstrated to be equivalent with both the QuikRead go Instrument and the QuikRead go Plus Instrument.

CRP

Measuring range

The measuring range of the CRP assay is 5–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 5–120 mg/l with serum or plasma samples. If the hematocrit is lower or higher than 40%, the measuring range changes in the manner shown in the table below (i.e. sliding result range). The CRP result is not displayed if the level of hematocrit is outside the range 15–75%.

Sample material	Hematocrit (%)	Sample volume	Measuring range (mg/l CRP)	
Whole blood	15–19	20 µl	5–150	In cases when the result is outside the measuring range, the result will be displayed as "CRP < 5 mg/l" or, for instance, as "CRP > 200 mg/l".
	20–28		5–160	
	29–35		5–180	
	36–41		5–200	
	42–46		5–220	
	47–50		5–240	
	51–53		5–260	
	54–56		5–280	
	57–60		5–300	
	61–62		5–330	
	63–64		5–340	
	65–66		5–360	
	67–68		5–390	
	69–70		5–410	
	71–72		5–440	
	73–74		5–470	
	75		5–510	
Plasma/serum	–	20 µl 12 µl*	5–120 5–200	

*It is possible to pipette a sample volume of 12 µl in order to have a wider measuring range. In that case, you have to change the test parameter settings of the instrument. Please see the instrument instructions for use.

Analytical sensitivity (limit of detection)

The detection limit of the CRP assay is < 5 mg/l. Limit of detection was performed applying the CLSI guideline EP17-A:2004. 6 blank and 6 low level samples were measured in 10 replicates.

Linearity

The recoveries in serial dilutions of CRP plasma samples have ranged from 96% to 101%. A linearity study was performed applying the CLSI guideline EP06-A:2003. 6 dilutions were tested as in 4 replicates.

Precision

A precision study was performed applying the CLSI guideline EP5-A2:2004. QuikRead go CRP Control, QuikRead go CRP Control High and 3 whole blood (20 µl) and plasma (12 µl) samples were analyzed for 20 days, 2 runs per day and 2 replicates per run. In total, 80 replicates per sample were measured.

Sample material	Specimen	Number of days	Mean CRP (mg/l)	Repeatability		Between-day precision		Reproducibility	
				SD	CV (%)	SD	CV (%)	SD	CV (%)
Whole blood	Low	20	9	0.4	4.5	0.4	4.3	0.5	6.2
	Medium		52	0.7	1.4	1.3	2.4	1.6	3.2
	High		177	3.3	1.9	4.9	2.8	6.5	3.6
Plasma	Low	20	10	0.4	3.7	0.0	0.1	0.4	4.2
	Medium		55	0.8	1.5	0.2	0.4	1.0	1.9
	High		148	3.0	2.0	1.5	1.0	3.5	2.4
Control	153764	20	24	0.6	2.3	0.2	0.8	0.7	2.8
Control High	153763	20	77	0.8	1.0	0.3	0.4	1.0	1.2

Antigen excess

CRP concentrations less than 1000 mg/l do not give falsely low results.

Interference

Interference study was performed applying the CLSI guideline EP07-A:2002. The interferents for CRP assay were compared to reference sample or results.

Interferent	No interference found up to concentration
Bilirubin	400 µmol/l
Vitamin C	200 µmol/l
Triglycerides	11.5 mmol/l
Cholesterol	9.0 mmol/l
Rheumatoid factors (RF)	525 IU/ml
Leucocytes	145x10 ⁹ /l
Anticoagulants (Li-heparin or EDTA)	No interference

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test because the assay antibodies lack the Fc part. In rare cases, interference by IgM myeloma protein has been observed.

Hb

Measuring range

The measuring range of the hemoglobin assay in whole blood samples is 50–245 g/l (5.0–24.5 g/dl, 3.1–15.2 mmol/l).

Linearity

The recoveries in serial dilutions of Hb samples have ranged from 100% to 105%. A linearity study was performed applying the CLSI guideline EP06-A:2003. 4 dilutions were tested as in 4 replicates.

Precision

A precision study was performed applying the CLSI guideline EP5-A2:2004. 3 whole blood samples were analyzed for 20 days, 2 runs per day and 2 replicates per run. In total, 80 replicates per sample were measured. QuikRead go Hb Control was analyzed for 5 days, 2 runs per day and 2 replicates per run. In total, 20 replicates per sample were measured.

Sample material	Specimen	Number of days	Mean Hb (g/l)	Repeatability		Between-day precision		Repro-ducibility	
				SD	CV (%)	SD	CV (%)	SD	CV (%)
Whole blood	Low	20	96	1.7	1.7	1.2	1.2	2.2	2.2
	Medium		134	2.3	1.7	1.7	1.3	2.9	2.1
	High		164	2.5	1.5	1.5	0.9	3.6	2.2
Control	141154	5	117	1.8	1.5	1.5	1.3	2.4	2.1

Interference

Interference study was performed applying the CLSI guideline EP07-A:2002. The interferents for Hb assay were compared to reference sample or results.

Interferent	No interference found up to concentration
Bilirubin	200 µmol/l
Vitamin C	200 µmol/l
Triglycerides	11.5 mmol/l
Rheumatoid factors (RF)	525 IU/ml
Creatinine	400 mg/l
Leucocytes	65 x 10 ⁹ /l
Thrombocytes	1000 x 10 ⁹ /l
β-carotene	3.7 µmol/l
Ibuprofen	500 mg/l
Acetylsalicylic acid	650 mg/l
Tetracyclines	200 mg/l
Urea	5000 mg/l
Anticoagulants (Li-heparin or EDTA)	No interference

13 Traceability

CRP

The calibrators used to calibrate the CRP assay of the QuikRead go CRP+Hb test are traceable to the ERM®-DA474/IFCC reference material. The values of the manufacturer’s calibrators are directly transferred to the test measurement. The instrument automatically reads the lot-specific calibration data from the cuvette. Hence, calibration by the user is not needed. The difference between lots was analyzed as average of 3 replicates of ERM®-DA474/IFCC reference material per lot.

Specimen	Certified value (uncertainty) mg/l	QuikRead go CRP+Hb results (mg/l CRP)		
		Lot 1	Lot 2	Lot 3
ERM®-DA474/IFCC	41.2 (2.5)	43.0	42.3	40.7

Hb

The Hb assay of the QuikRead go CRP+Hb test is traceable to the ICSH (International Council for Standardization of Haematology) standard 1995 for cyanmethemoglobin and the BCR-522 bovine blood lysate certified reference material.^{28,29} ICSH comparison study was performed applying the CLSI guideline EP9-A2:2002. In a comparison of 61 patient whole blood samples, QuikRead go CRP+Hb (y) and ICSH reference method (x) were found comparable. Regression equation for the comparability was y = 1.07x – 7.8 and r = 0.996. Samples were measured as 2 replicates for both methods.

The difference between lots was analyzed as average of 5 replicates of hemoglobin control per lot.

Specimen	QuikRead go CRP+Hb results (g/l Hb)		
	Lot 1	Lot 2	Lot 3
Control	110	111	111

14 Method comparison

CRP

Method comparison

Patient plasma samples were measured using two clinical laboratory methods, and the QuikRead go CRP+Hb method. Summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method 1	Clinical laboratory method 2
y = 1.00x + 1.0 r = 0.997 N = 116	y = 0.93x + 1.4 r = 0.993 N = 113

Whole blood versus plasma comparability

In a comparison of 104 patient samples, whole blood (y) and plasma (x) were found comparable. Regression equation for the comparability was y = 1.03x – 0.3 and r = 0.995.

Hb

Method comparison

Patient whole blood samples were measured using a clinical laboratory method, a point of care test and the QuikRead go CRP+Hb method. Summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method	Point of care test
y = 1.07x – 7.8 r = 0.996 N = 61	y = 1.01x – 1.1 r = 0.994 N = 130

15 Disposal

Dispose of contents according to national and local law. Treat and dispose of all patient samples, controls, used capillaries, plungers, reagent caps, and cuvettes as potentially infectious materials. Discard unused reagent caps as hazardous waste containing synthetic polymer microparticles. Prevent release to the environment. For more information, please refer to the component-specific Safety Data Sheets (SDS).

Materials of the components of the QuikRead go CRP+Hb kit:

- Paper: Instructions for use
- Cardboard: Kit box
- Plastic: Cuvette racks, plungers, plunger and capillary tubes, and filling inside of capillary tube
- Glass: Capillaries
- Metal: Reagent cap tubes, plunger and capillary tube caps, and cuvette foil covers
- Several: Reagent caps, lids of reagent cap tubes, foil pouches covering cuvette racks, and unused cuvettes

When used in accordance with Good Laboratory Practice, good occupational hygiene, and the instructions for use, the reagents supplied should not present a hazard to health.

16 Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, consult the instrument instructions for use.

Error messages	Corrective action
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner blue part of the cap is not pressed down.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18–25 °C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18–25 °C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed, the cuvette may contain dirt or the sample interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. If this error message appears again, please make sure that the reagent caps have not been exposed to humidity. Also ensure that the reagent cap tube is properly closed.
No result. Hematocrit too low.	Whole blood sample's hematocrit is below < 15% and the test cannot be completed. In case of plasma/serum sample, the sample is clearly hemolyzed. Please collect a new venous sample.
No result. Hematocrit too high.	Whole blood sample's hematocrit is above > 75% and the test cannot be completed.

Unexpected low and high results

Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected low result.	Too small sample volume.	Perform a new test. Make sure that the capillary is completely filled. Avoid air bubbles.
Unexpected high result.	Too large sample volume.	Perform a new test. Make sure that the capillary is properly filled. Wipe away any excess sample from the outside of the capillary.
	Too low buffer volume.	Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
	The cuvette is dirty.	Perform a new test. Do not touch the lower part of the cuvette.
Unexpected low/high result.	Use of components from different kit lots or tests.	Perform a new test. Make sure that the reagents are from the same kit lot and the sample volume is correct.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored in the manner described in the instructions for use.
	Incorrect sample volume: plasma/serum/ control sample.	Check, that the used sample volume (20 or 12 µl) matches the set plasma/serum sample volume on the instrument.

1.

Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. A rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. Clinical chemistry and laboratory medicine 1999; 37(2): 109-113.

2.

Pepys MB, Hirschfield GM. C-reactive protein: a critical update. The Journal of clinical investigation 2003; 111(12): 1805-1812.

3.

Tonkin-Crane SKG, Tan PS, van Hecke O, et al. Clinician-targeted interventions to influence antibiotic prescribing behaviour for acute respiratory infections in primary care: an overview of systematic reviews (Review). Cochrane Database of Systematic Reviews 2017; 7(9): CD012252.

4.

Bruns AH, Oosterheert JJ, Hak E, Hoepelman AI. Usefulness of consecutive C-reactive protein measurements in follow-up of severe community-acquired pneumonia. The European respiratory journal 2008; 32(3): 726-732.

5.

Coelho L, Povoia P, Almeida E, et al. Usefulness of C-reactive protein in monitoring the severe community-acquired pneumonia clinical course. Critical care (London, England) 2007; 11(4): R92.

6.

Coelho LM, Salluh JI, Soares M, et al. Patterns of c-reactive protein RATIO response in severe community-acquired pneumonia: a cohort study. Critical care (London, England) 2012; 16(2): R53.

7.

Povoia P, Teixeira-Pinto AM, Carneiro AH. C-reactive protein, an early marker of community-acquired sepsis resolution: a multi-center prospective observational study. Critical care (London, England) 2011; 15(4): R169.

8.

Hengst JM. The role of C-reactive protein in the evaluation and management of infants with suspected sepsis. Adv Neonatal Care 2003; 3(1): 3-13.

9.

Benitz WE, Han MY, Madan A, Ramachandra P. Serial serum C-reactive protein levels in the diagnosis of neonatal infection. Pediatrics 1998; 102(4): E41.

10.

Gans SL, Ateama JJ, van Dieren S, Groot Koerkamp B, Boermeester MA. Diagnostic value of C-reactive protein to rule out infectious complications after major abdominal surgery: a systematic review and meta-analysis. International journal of colorectal disease 2015; 30(7): 861-873.

11.

Straatman J, Harmsen AMK, Cuesta MA, Berkhof J, Jansma EP, van der Peet DL. Predictive Value of C-Reactive Protein for Major Complications after Major Abdominal Surgery: A Systematic Review and Pooled-Analysis. PloS one 2015; 10(7): e0132995.

12.

Butler CC, Gillespie D, White P, et al. C-Reactive Protein Testing to Guide Antibiotic Prescribing for COPD Exacerbations. New England Journal of Medicine 2019; 381(2): 111-120.

13.

Smolen JS, Aletaha D, Barton A, et al. Rheumatoid arthritis. Nature reviews Disease primers 2018; 4: 18001.

14.

Chang S, Malter L, Hudesman D. Disease monitoring in inflammatory bowel disease. World journal of gastroenterology 2015; 21(40): 11246-11259.

15.

Freeman AM, Rai M, Morando DW. Anemia Screening. Treasure Island (FL): StatPearls Publishing; 2021.

16.

World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 Rev.2, 2002.

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18. Ehl S et al. C-reactive protein is a useful marker for guiding duration of antibiotic therapy in suspected neonatal bacterial infection. *Pediatrics* 1997; 99:216-221.

19. Meem M et al. Biomarkers for diagnosis of neonatal infections: A systematic analysis of their potential as a point-of-care diagnostics. *J Glob Health* 2011; 1:201-209.

20. Zhang Y et al. Association between C-reactive protein concentration and chronic obstructive pulmonary disease: a systematic review and meta-analysis. *J Int Med Res* 2012; 40:1629-1635.

21. van Durme YM et al. C-reactive protein levels, haplotypes, and the risk of incident chronic obstructive pulmonary disease. *Am J Respir Crit Care Med* 2009; 179:375-382.

22. Dessein PH et al. High sensitivity C-reactive protein as a disease activity marker in rheumatoid arthritis. *J Rheumatol.* 2004; 31 (6):1095-7.

23. Shine et al. Solid phase immunoassays for human C-reactive protein. *Clin Chim Acta.* 1981; 117: 13-23.

24. Bates I. Reference Ranges and Normal Values. *Dacie and Lewis Practical Haematology.* 12th ed: Elsevier; 2017: 8-17.

25. World Health Organization. Haemoglobin concentrations for the diagnosis of anaemia and assessment of severity. *Vitamin and Mineral Nutrition Information System.* Geneva, WHO, 2011.

26. Chiesa C et al. C reactive protein and procalcitonin: Reference intervals for preterm and term newborns during the early neonatal period. *Clin Chim Acta.* 2011; 412 (11-12): 1053–9.

27. Schlebusch H et al. High Sensitive CRP and Creatinine: Reference intervals from Infancy to Childhood. *J Lab Med* 2002; 26 (5/6): 341-6.

28. International Council for Standardisation in Haematology: Expert Panel on Hemoglobinometry: Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1995) and specifications for international haemoglobinocyanide standard (4th edition). *J Clin Pathol* 1996; 49: 271–4.

29. Bull, BS et al. Reference and selected procedures for the quantitative determination of hemoglobin in blood; approved standard – third edition*, *NCCLS* 20 (28), document H15-A3.

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Aidian, QuikRead and QuikRead go are trademarks or registered trademarks of Aidian Oy.