Revision History/ Revisionshistorie/ Historique des révisions/ Cronologia delle revisioni

REF 153287 QuikRead go easyCRP

Material number and release version/ Materialnummer und Version und Veröffentlichungsdatum/ Numéro et version du document et date de publication/ Numero materiale e versione e data di rilascio	Change classification/ Klassifizierung ändern/ Classification de la modification/ Variazione della classificazione	Section of change/ Abschnitt der Änderung/ Section de la modification/ Sezione della variazione
153289-8, 05/2026	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Intended purpose/ Verwendungszweck/ Utilisation prévüe/ Finalità d'uso
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Summary and explanation of the analyte/ Zusammenfassung und Beschreibung des Analyten/ Résumé et explication de l'analyte/ Riepilogo e spiegazione dell'analita
	General change/ Allgemeine Änderungen/ Modification générale/	Principles of the procedure/ Verfahrensprinzipien/ Principe de la procédure/ Principi della procedura
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Reagents/ Reagenzien/ Réactifs/ Reagenti
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Reagents/ Reagenzien/ Réactifs/ Reagenti
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Warnings and precautions/ Warn- und Entsorgungshinweis/ Précautions d'emploi/ Avvertimenti e precauzioni
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Sample collection and preparation/ Probennahme und - vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni

	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Sample collection and preparation/ Probennahme und -vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni
	General change/ Allgemeine Änderungen/ Modification générale/	Sample collection and preparation/ Probennahme und -vorbereitung/ Prélèvement et préparation des échantillons/ Raccolta e preparazione dei campioni
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Procedure/ Verfahren/ Procédure/ Procedura
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Procedure/ Verfahren/ Procédure/ Procedura
	General change/ Allgemeine Änderungen/ Modification générale/	Procedure/ Verfahren/ Procédure/ Procedura
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Quality control/ Qualitätskontrolle/ Contrôle de qualité/ Controllo di qualità
	General change/ Allgemeine Änderungen/ Modification générale/	Quality control/ Qualitätskontrolle/ Contrôle de qualité/ Controllo di qualità
	General change/ Allgemeine Änderungen/ Modification générale/	Interpretation of results/ Interpretation der Testergebnisse/ Interprétation des résultats/ Interpretazione dei risultati
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Performance characteristics/ Leistungsparameter/ Performances/ Caratteristiche delle prestazioni
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Traceability/ Nachweisbarkeit/ Traçabilité/ Tracciabilità
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Method comparison/ Methodenvergleich/ Comparaison de méthodes/ Comparaison de méthodes
	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Disposal/ Entsorgung/ Mise au rebut/ Smaltimento
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Disposal/ Entsorgung/ Mise au rebut/ Smaltimento
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Troubleshooting/ Fehlerdiagnose/ Résolution des problèmes/ Risoluzione dei Problemi
	General change/ Allgemeine Änderungen/ Modification générale/	Troubleshooting/ Fehlerdiagnose/ Résolution des problèmes/ Risoluzione dei Problemi
	Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	References/ Referenzen/ Références/ Riferimenti bibliografici

	Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Explanation of symbols/ Zeichenerklärung/ Explication des symboles/ Spiegazione dei simboli
--	---	--

Classification/ Klassifizierung/ Classification/ Classificazione	Description/ Beschreibung/ Description/ Descrizione
New instructions for use/ Neue Anwendungshinweise/ Nouveau mode d'emploi/ Nuove istruzioni per l'uso	First version of a new instructions for use/ Erste Version der neuen Anwendungshinweise/ Première version d'un nouveau mode d'emploi/ Prima versione delle nuove istruzioni per l'uso/
Correction/ Korrektur/ Correction/ Correzione	Correction of documentation anomalies (minor typographical, grammar, and formatting changes are not included)/ Korrektur von Anomalien in der Dokumentierung (kleinere typografische, grammatikalische und Formatierungsänderungen sind nicht eingeschlossen)/ Correction d'anomalies dans la documentation (les modifications typographiques, grammaticales et de formatage mineures ne sont pas incluses)/ Correzione delle anomalie della documentazione (i piccoli refusi, la grammatica e le variazioni di formattazione non sono inclusi)
Technical change/ Technische Änderungen/ Modification technique/ Variazione tecnica	Addition, revision and/or removal of information related to the product/ Hinzufügung, Änderung und/oder Entfernung von Informationen über das Produkt/ Ajout, révision et/ou suppression d'informations liées au produit/ Aggiunta, revisione e/o rimozione delle informazioni relative al prodotto
Regulatory change/ Regulatorische Änderungen/ Modification réglementaire/ Variazione normativa	Changes due to regulations, directives, laws and/or registration/ Änderungen aufgrund von Verordnungen, Richtlinien, Gesetzen und/oder Registrierung/ Modifications dues aux réglementations, directives, lois et/ou à l'enregistrement/ Variazioni dovute a normative, direttive, leggi e/o registrazione
General change/ Allgemeine Änderungen/ Modification générale/ Variazione generale	Implementation of non-technical changes/ Implementierung von nichttechnischen Änderungen/ Mise en œuvre de modifications non technique/ Implementazione di variazioni non tecniche

QuikRead go easy CRP

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To be used with the QuikRead go Instrument and the QuikRead go Plus Instrument. Hereafter referred to as the QuikRead go instrument.

1. Intended purpose

The QuikRead go easy CRP test is an immunoturbidimetric assay for the quantitative determination of C-reactive protein (CRP) values from a capillary blood, anticoagulated (EDTA or heparin) venous whole blood or plasma, or serum sample. The test is carried out using an automated instrument belonging to the QuikRead go system. The assay is intended to be used by healthcare professionals in the clinical laboratory and near-patient testing (NPT) setting (also known as point of care setting). CRP is an acute phase protein present in low concentrations in the blood of healthy individuals. Infections, tissue injuries and inflammatory disorders may cause a non-specific elevation of the CRP level in the patient's blood. Quantitative CRP measurement is a sensitive and reliable indicator of the acute phase response. The measurement of CRP aids in the diagnosis, evaluation, and monitoring of infection, tissue injury, inflammatory disorders, and associated diseases^{1,2}. For *in vitro* diagnostic use only.

2. Summary and explanation of the analyte

CRP is an acute phase protein present in low concentrations in the blood of healthy individuals. Infections, tissue injuries, and inflammatory disorders may cause a non-specific elevation of the CRP level in the patient's blood. CRP measurement is a sensitive and reliable indicator of the acute phase response. The rise of CRP concentration is rapid, and increased levels can be detected within 6 hours of the onset of the inflammatory process.^{1,2} Quantitative measurement of CRP NPT decreases antibiotic prescribing for approximately 25% in different acute respiratory infections³ and exacerbations of chronic obstructive pulmonary disease (COPD)⁴, and therefore can be used to guide antibiotic treatment. CRP has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections⁵⁻⁸. CRP is also suitable for the evaluation of neonatal infections^{9,10} and it is an effective tool in ruling out postoperative infections^{11,12}. In systemic inflammation, such as rheumatoid arthritis and inflammatory bowel disease, CRP measurement is routinely used as part of the initial assessment of the disease and for monitoring of disease activity, progress, and treatment response^{13,14}.

3. Principles of the procedure

The CRP measurement is immunoturbidimetric and it is based on agglutination reaction. Microparticles are coated with anti-human CRP F(ab)₂ fragments, and the CRP present in the sample reacts with microparticles. The resultant change in the turbidity of the solution is measured by the QuikRead go instrument. The CRP assay correlates well with results obtained from several other immunoturbidimetric methods.

The sample is added to the buffer in a prefilled cuvette and the blood cells in a whole blood sample are hemolyzed. The calibration data of the test are found on the barcode label of the cuvette, which the instrument reads automatically before the start of the test. The CRP value is automatically calculated by the instrument. When whole blood sample is used, the CRP value is automatically corrected according to the sample's hematocrit level. Hematocrit results are not displayed but are used in the calculations. The measuring range of the CRP assay is 1–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 1–120 mg/l with serum or plasma samples.

4. Reagents

Content of the kit (Table 1)

The content of the kit is listed in Table 1. The reagents contain preservatives, please see Section 5. "Warnings and precautions".

Storage

The shelf-life of the unopened QuikRead go easy CRP kit at 2–25°C is until the expiry date marked on the kit box.

Component	Shelf-life at 2–25°C	In-use stability at 2–8°C	In-use stability at 18–25°C
Reagent caps	Until the expiry date of the kit in an unopened reagent cap tube	6 months in a closed reagent cap tube	
Prefilled cuvettes	Until the expiry date of the kit in an intact foil pouch	6 months	3 months
Sample collectors	Until the expiry date of the sample collector box		
		Prefilled cuvette without foil cover: 2 hours	

Reagent preparation and storage conditions

All reagents are ready to be used. Keep the reagent caps away from moisture. Close the tube immediately after taking out the required number of reagent caps.

Reagent deterioration

Do not use a cuvette with visible dirt in the buffer.

5. Warnings and precautions

Health and safety information

- For *in vitro* diagnostic use only.
- QuikRead go Sample Collectors are intended for non-invasive sampling without skin contact.
- Do not smoke, eat or drink on premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- All patient samples and controls should be handled as potentially infectious material.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagent inside a reagent cap contains < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032). Contains synthetic polymer microparticles. Do not release the cap reagent outside the cuvette. Microparticles are contained within the closed cuvette when used according to instructions.
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration.
- Disposal: see Section 15.
- In case of a serious incident, please report it to the manufacturer or its representative and national authority. Summary of safety and performance is available on Eudamed.

Analytical precautions

- The QuikRead go easy CRP kit is intended solely to be used with the QuikRead go instrument.
- Do not use the product after the expiration date marked on the outer package.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. In addition, always ensure before using an individual cuvette that its foil covering is intact.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components already used for performing a test.
- Do not exceed the stated stability periods for opened reagents.
- Hold the cuvette from the sides. Do not touch the data matrix or the lower part of the cuvette (the optical part). Discard any cuvettes with fingerprints.
- QuikRead go easy CRP Reagent Caps are color-coded turquoise with a cap label to distinguish them from other QuikRead analytes.
- Keep the reagent caps away from moisture. Close the tube immediately after taking out the required number of reagent caps. Do not use reagent caps if they have not been stored according to the instructions.
- Do not attach labels to the cuvette. If needed, the identification information can be written on the opposite side of the label using a permanent pen.
- Do not splash liquid into the measuring well of the instrument.

6. Sample collection and preparation

Sample material, sample volume and collection

QuikRead go Sample Collector 10µl (Cat. no. 153211) included in the QuikRead go easy CRP kit (Cat. no. 153287) is required for transferring the sample into a cuvette. The capillary part of the sample collector is covered with EDTA and the appearance might vary depending on crystallization of the EDTA.

Sample material	Sample volume	Sample collection
Capillary blood	10 µl	Prick the clean and dry finger/heel with a lancet. Wipe off the first drop and collect 10 µl of blood from the second drop into EDTA covered capillary part of the sample collector.
Anticoagulated whole blood	10 µl	Use a venous blood sample collected in a tube containing heparin or EDTA. Mix whole blood by inverting the tube several times and take 10 µl with the sample collector.
Plasma	10 µl	Use EDTA/heparin plasma. Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Take 10 µl with a sample collector.
Serum	10 µl	Avoid clearly hemolyzed samples. Mild hemolysis will not affect the test results. Take 10 µl with a sample collector.

Sample dilution

QuikRead go easy CRP test does not require sample dilution. However, if the displayed result is marked with a greater than (>) symbol, it is possible to dilute the sample to obtain a result within the measuring range.

Sample material	Instructions
Whole blood	Do not dilute whole blood samples.
Plasma/ serum	Plasma or serum samples may be diluted with a 0.9% NaCl solution before the sample is added to the cuvette. The recommended dilution ratio is 1+3 (by volume, 1 part sample + 3 parts 0.9% NaCl). Take 10 µl of the diluted sample into a cuvette with the sample collector, analyse the sample and multiply the obtained result by 4.

Sample storage

Sample material	Short-term storage	Long-term storage
Capillary blood	Max. 2 minutes in the sample collector	Do not store
Anticoagulated whole blood	2–8°C for 3 days	Separate plasma and store as indicated below
Plasma	2–8°C for 7 days	Below –20°C for 3 years ¹⁵
Serum	2–8°C for 7 days	Below –20°C for 3 years ¹⁵
Diluted plasma/serum in 0.9% NaCl	18–25°C for 30 minutes	Do not store
Patient sample in a cuvette	18–25°C for 2 hours	Do not store

Allow the samples to reach room temperature (18–25°C) before testing. Frozen samples must be thawed completely, mixed thoroughly, and brought to room temperature prior to testing. If the thawed plasma samples contain clots, the samples should be centrifuged. Samples should not be repeatedly frozen and thawed.

7. Procedure

Materials required but not provided (Table 2)

The materials listed in Table 2 are not provided. The materials provided are listed in Table 1.

Assay procedure

Before starting an assay, read and follow the QuikRead go instrument and QuikRead go easy CRP test instructions for use. Open the foil pouch of the cuvette rack. Remember to mark the opening date on the cuvette rack and the reagent cap tube.

Sample collection (see figures 1–6)

QuikRead go Sample Collector 10µl (Cat. no. 153211) is required for transferring the sample into a prefilled cuvette. The sample collectors can be purchased separately, whereas the QuikRead go easy CRP kit (Cat. no. 153287) contains them already.

1. A prefilled cuvette must reach room temperature (18–25°C) before use. This takes 15 minutes for an individual refrigerated (2–8°C) prefilled cuvette. Do not touch the data matrix or the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on the results. The test must be performed within two (2) hours of opening the cuvette.
2. Open the sample collector box slightly to get one single the sample collector out of the box. The package filling can be discarded.
3. Hold the sample collector almost horizontally and take the sample by placing the tip of the sample collector just beneath the surface of the blood drop. Fill the capillary part of the sample collector completely.
4. Check that the sample collector is completely filled and there are no air bubbles in the capillary. Do not wipe the sample collector. If there is excess of sample on the sample collector, discard the sample collector and take a new sample with a new sample collector.
5. Within 2 minutes insert the sample collector into the cuvette. The sample collector remains in the cuvette during the whole test measurement procedure.
6. Close the cuvette with the reagent cap. Do not press down the turquoise inner part of the reagent cap. Close the reagent cap tube immediately after taking out the required number of reagent caps. Once the sample collector has been transferred into the cuvette the assay should be run within two (2) hours. Keep the cuvette upright and do not shake.

Note! If a venous sample is used, it is recommended to take the sample directly from the tube. Tilt the tube slightly and collect the sample from just beneath the surface of the sample with the sample collector. Do not wipe the sample collector.

Analysing the sample (see figures 7–9)

Begin the test by turning on the QuikRead go instrument. For full instructions, see the instrument instructions for use. The screen will guide you through the test.

7. Choose **Measure / Measure sample** on the screen of the QuikRead go instrument.
8. Place the cuvette into the instrument with the barcode facing you. The instrument detects the sample type. The screen shows how the measurement is proceeding.
9. The result is displayed on the screen, and the cuvette rises automatically when the measurement is completed.

Analysing the control

Please see the instructions for use provided with each control vial. Analyse the control samples as patient samples but choose **Quality control** on the screen of the QuikRead go instrument. The result will be stored in a separate result file for quality control results.

8. Quality control

Quality control testing must be performed in accordance with applicable national or local regulations or accreditation requirements and laboratory standard Quality Control Procedures. Regular use of QuikRead go CRP Control Low (Cat. no. 153765), QuikRead go CRP Control (Cat. no. 153764) and/ or QuikRead go CRP Control High (Cat. no. 153763) is recommended.

It is recommended to perform the control at least with each new shipment and new reagent lot and every time an unexpected result is obtained. The control should be used each time new QuikRead go system users are trained. Please consult Section 16. “Troubleshooting” and contact Aidian or your local distributor if you observe any changes in test functionality.

QuikRead go CRP controls: All QuikRead go CRP controls are ready-to-use. The control values have been determined for the QuikRead go instrument and the control is applied using the same procedure as for a plasma/serum sample. Perform the test as outlined in Section 7. “Procedure”, use the same control sample volume as for unknown blood specimens (10 µl). Tilt the control vial until the control solution reaches the bottleneck. Place the sample collector almost horizontally just below the control surface. Fill the capillary part completely. The control is stable for 15 minutes in the cuvette.

Other commercially available CRP controls: Handle and use the controls in accordance with the instructions for each control substance. Perform the test as outlined in Section 7. “Procedure”. Acceptable limits for the control must be determined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells, because these might not hemolyze normally.

9. Interpretation of results

Elevation of CRP value is unspecific and the results should be interpreted in the light of other clinical findings. Serial CRP measurements may be required to monitor a patient’s condition or response to therapy and when measuring CRP levels in newborns. Additionally, establishing the patient-specific baseline CRP level is recommended when assessing slightly elevated CRP levels.

CRP test result	Interpretation of result
< 10 mg/l	•Excludes many acute inflammatory diseases but does not specifically exclude inflammatory processes. ¹⁶ •Possible inflammations in newborn babies are not necessarily excluded. ^{10,17,18} •Possible inflammations or exacerbations in conditions with systemic inflammation are not necessarily excluded. ^{19–22}
10–50 mg/l	Elevated concentrations are found in acute disease that occur in the presence of a slight to moderate inflammatory process. ¹⁶
> 50 mg/l	Indicates high and extensive inflammatory activity. ¹⁶

If diluted plasma or serum samples have been used, the final results must be calculated in accordance with the instructions in Section 6. “Sample collection and preparation”. If necessary, please consult Section 16. “Troubleshooting”, and repeat the test.

10. Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. “Performance characteristics”. Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

11. Expected values

Reference value: < 3 mg/l^{22–24}, median value in apparently healthy adults: 0.8 mg/l²². CRP levels in newborns can be affected by e.g. prematurity, duration of delivery and complications related to the delivery. Immediately after birth the CRP levels in term neonates are generally low and will increase to some extent (median below 2 mg/l) over the next few days²³.

12. Performance characteristics

The performance characteristics of the QuikRead go easy CRP test have been demonstrated to be equivalent with both the QuikRead go Instrument and the QuikRead go Plus Instrument.

Measuring range

The measuring range of the CRP assay is 1–200 mg/l (at hematocrit level of 40%) when whole blood samples are used and 1–120 mg/l with serum or plasma samples. If the hematocrit is lower or higher than 40%, the measuring range changes in the manner shown in the table below (i.e. sliding result range). The CRP result is not displayed if the level of hematocrit is outside the range 15–75%.

Sample material	Hematocrit (%)	Sample volume	Measuring range (mg/l CRP)	
Whole blood	15–19	10 µl	1.0–140	In cases when the result is outside the measuring range, the result will be displayed as "CRP < 1.0 mg/l" or, for instance, as "CRP > 200 mg/l"
	20–24		1.0–150	
	25–29		1.0–160	
	30–33		1.0–170	
	34–36		1.0–180	
	37–39		1.0–190	
	40–42		1.0–200	
	43–45		1.1–210	
	46–50		1.2–220	
	51–53		1.3–240	
	54–57		1.4–260	
	58–60		1.5–280	
	61–64		1.7–300	
	65–68		1.9–340	
	69–71		2.1–380	
	72–75		2.4–420	
Plasma/serum	–	10 µl	1.0–120	

It is possible to turn off the sliding result range setting, in which case the measuring range for whole blood samples is fixed to 2.4–140 mg/l. In order to turn off the sliding result range, you have to change the test parameter settings of the instrument. Please see the instrument instructions for use.

Analytical sensitivity (limit of detection)

The detection limit of the assay is < 1 mg/l. Limit of detection was performed applying the Clinical and Laboratory Standards Institute (CLSI) guideline EP17-A2:2012. 4 blank and 5 low level samples were measured on 3 days with 5 replicates per day, to yield 15 results per sample.

Linearity

The recoveries in serial dilutions of CRP plasma samples have ranged from 85% to 104%. A linearity study was performed applying the CLSI guideline EP06-A:2003. 11 dilutions were tested as in 4 replicates.

Precision

A precision study was performed applying the CLSI guideline EP05-A3:2014. Whole blood samples were analyzed for 5 days with 3 instruments and 5 replicates per run. In total, 75 replicates per sample were measured. QuikRead go CRP Control Low, QuikRead go CRP Control, QuikRead go CRP Control High and plasma samples were analyzed for 20 days, 2 runs per day and 2 replicates per run. In total, 80 replicates per sample were measured.

Sample material	Specimen	Number of days	Mean CRP (mg/l)	Repeatability		Between-day precision		Reproducibility	
				SD	CV (%)	SD	CV (%)	SD	CV (%)
Whole blood	Low	5	14	0.4	2.9	0.1	0.9	0.4	3.2
	Medium		60	2.3	3.8	0.9	1.5	2.8	4.7
	High		175	14.0	8.0	3.9	2.2	14.5	8.3
Plasma	Low	20	1.1	0.08	6.8	0.02	1.7	0.09	7.8
	Medium		44	1.5	3.4	0.6	1.4	1.7	3.7
	High		96	8.5	8.8	0.5	0.5	8.5	8.9
Control Low	153765	20	2.5	0.11	4.5	0.02	0.9	0.12	4.9
Control	153764	20	24	0.9	3.7	0.2	0.9	0.9	3.8
Control High	153763	20	77	3.9	5.1	1.0	1.3	4.3	5.6

Antigen excess

CRP concentrations of less than 1000 mg/l do not give falsely low results.

Interference

Interference study was performed applying the CLSI guideline EP07-A2:2005. The interferents were compared to reference sample or results.

Interferent	No interference found up to concentration
Bilirubin	400 µmol/l
Vitamin C	200 µmol/l
Triglycerides	11.5 mmol/l
Cholesterol	10.0 mmol/l
Rheumatoid factors (RF)	525 IU/ml
Leucocytes	87x10 ⁹ /l
Anticoagulants (Li-heparin or EDTA)	No interference

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test because the assay antibodies lack the Fc part. In rare cases, interference by IgM myeloma protein has been observed.

13. Traceability

The calibrators used to calibrate the CRP assay of the QuikRead go easy CRP test are traceable to the ERM®-DA474/IFCC reference material. The values of the manufacturer’s calibrators are directly transferred to the test measurement. The QuikRead go instrument automatically reads the lot-specific calibration data from the cuvette. Hence, calibration by the user is not needed. The difference between lots was analyzed as average of 3 replicates of ERM®-DA474/IFCC reference material per lot.

Specimen	Certified value (uncertainty) mg/l	QuikRead go easy CRP results (mg/l)	
		Lot 1	Lot 2
ERM®-DA474/IFCC	41.2 (2.5)	40.3	41.0

14. Method comparison

Patient plasma samples were measured using one clinical laboratory method and one point of care method and QuikRead go easy CRP method. A summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method	Point of care test
y = 1.01x - 0.5 r = 0.988 N = 92	y = 0.93x + 1.2 r = 0.982 N = 89

Whole blood versus plasma comparability

In a comparison of 95 clinical samples, whole blood (y) and plasma (x) were found comparable. Regression equation for the comparability was y = 0.99x - 0.4 and r = 0.976.

15. Disposal

Dispose of contents according to national and local law. Treat and dispose of all patient samples, controls, used sample collectors, reagent caps, and cuvettes as potentially infectious materials. Discard unused reagent caps as hazardous waste containing synthetic polymer microparticles. Prevent release to the environment. For more information, please refer to the component-specific Safety Data Sheets (SDS).

Materials of the components of the QuikRead go easy CRP kit:

- Paper: Instructions for use
- Cardboard: Kit and sample collector boxes
- Plastic: Cuvette racks, and filling inside of sample collector box
- Metal: Reagent cap tubes, and cuvette foil covers
- Several: Reagent caps, lids of reagent cap tubes, foil pouches covering cuvette racks, unused cuvettes, and unused sample collectors

When used in accordance with Good Laboratory Practice, good occupational hygiene, and the instructions for use, the reagents supplied should not present a hazard to health.

16. Troubleshooting

Error messages

Error messages displayed by the QuikRead go instrument are listed below. For more detailed information on error messages, please consult the instrument instructions for use.

Error messages	Corrective actions
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner turquoise part of the cap is not pressed down.
Lot number of the reagent cap and cuvette do not match.	Discard the test. Perform a new test with new components. Make sure that the reagents are from the same kit lot.
Unable to read cap label.	Check from the reagent cap tube that the reagent cap is from the same kit lot than the cuvette. If the reagent cap and cuvette are from the same kit lot, you can continue the measurement. If the lots are not the same, cancel the measurement, discard the test and perform a new test with new correct components.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18–25°C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18–25°C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed, the cuvette may contain dirt or the sample interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. If this error message appears again, please make sure that the reagent caps have not been exposed to humidity. Also ensure that the reagent cap tube is properly closed.
No result. Hematocrit too low.	Whole blood sample's hematocrit is below < 15% and the test cannot be completed. In case of plasma/serum sample, the sample is clearly hemolyzed. Please collect a new venous sample.
No result. Hematocrit too high.	Whole blood sample's hematocrit is above > 75% and the test cannot be completed.

Unexpected low and high results

Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected low result.	Too small sample volume.	Perform a new test. Make sure that the sample collector is completely filled and there are no air bubbles in the capillary part. It is not possible to overfill it. Do not wipe the sample collector.
	Sample too long in the sample collector.	Perform a new test. Sample can be stored max. 2 minutes in the sample collector.
Unexpected high result.	Too large sample volume.	Perform a new test. Take the sample by placing the tip of the sample collector just beneath the surface of the sample and fill the capillary part completely. Do not wipe the sample collector. Make sure that there is no sample on the turquoise part of the sample collector.
	Too low buffer volume.	Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
	The cuvette is dirty.	Perform a new test. Do not touch the lower part of the cuvette.
Unexpected low/high result.	Use of reagents from different kit lots or tests.	Perform a new test. Make sure that the reagents are from the same kit lot.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored in the manner described in the instructions for use.
	Use of incorrect sample collector.	Perform a new test. Use QuikRead go Sample Collector 10µl for transferring the sample into the cuvette. Do not use pipette or other sampling apparatus.

References / Referenzen / Références / Riferimenti bibliografici

1. Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. A rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. Clinical chemistry and laboratory medicine 1999; 37(2): 109-113.

2. Pepys MB, Hirschfield GM. C-reactive protein: a critical update. The Journal of clinical investigation 2003; 111(12): 1805-1812.

3. Tonkin-Crane SKG, Tan PS, van Hecke O, et al. Clinician-targeted interventions to influence antibiotic prescribing behaviour for acute respiratory infections in primary care: an overview of systematic reviews (Review). Cochrane Database of Systematic Reviews 2017; 7(9): CD012252.

4. Butler CC, Gillespie D, White P, et al. C-Reactive Protein Testing to Guide Antibiotic Prescribing for COPD Exacerbations. New England Journal of Medicine 2019; 381(2): 111-120.

5. Bruns AH, Oosterheert JJ, Hak E, Hoepelman AI. Usefulness of consecutive C-reactive protein measurements in follow-up of severe community-acquired pneumonia. The European respiratory journal 2008; 32(3): 726-732.

6. Coelho L, Pova P, Almeida E, et al. Usefulness of C-reactive protein in monitoring the severe community-acquired pneumonia clinical course. Critical care (London, England) 2007; 11(4): R92.

7. Coelho LM, Salluh JI, Soares M, et al. Patterns of c-reactive protein RATIO response in severe community-acquired pneumonia: a cohort study. Critical care (London, England) 2012; 16(2): R53.

8. Pova P, Teixeira-Pinto AM, Carneiro AH. C-reactive protein, an early marker of community-acquired sepsis resolution: a multi-center prospective observational study. Critical care (London, England) 2011; 15(4): R169.

9. Hengst JM. The role of C-reactive protein in the evaluation and management of infants with suspected sepsis. Adv Neonatal Care 2003; 3(1): 3-13.

10. Benitz WE, Han MY, Madan A, Ramachandra P. Serial serum C-reactive protein levels in the diagnosis of neonatal infection. Pediatrics 1998; 102(4): E41.

11. Gans SL, Atema JJ, van Dieren S, Groot Koerkamp B, Boermeester MA. Diagnostic value of C-reactive protein to rule out infectious complications after major abdominal surgery: a systematic review and meta-analysis. International journal of colorectal disease 2015; 30(7): 861-873.

12. Straatman J, Harmsen AMK, Cuesta MA, Berkhof J, Jansma EP, van der Peet DL. Predictive Value of C-Reactive Protein for Major Complications after Major Abdominal Surgery: A Systematic Review and Pooled-Analysis. PloS one 2015; 10(7): e0132995.

13. Smolen JS, Aletaha D, Barton A, et al. Rheumatoid arthritis. Nature reviews Disease primers 2018; 4: 18001.

14. Chang S, Malter L, Hudesman D. Disease monitoring in inflammatory bowel disease. World journal of gastroenterology 2015; 21(40): 11246-11259.

15. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples WHO/DIL/LAB/99.1 Rev.2, 2002.

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17. Ehl S et al. C-reactive protein is a useful marker for guiding duration of antibiotic therapy in suspected neonatal bacterial infection. Pediatrics 1997; 99:216-221.

18. Meem M et al. Biomarkers for diagnosis of neonatal infections: A systematic analysis of their potential as a point-of-care diagnostics. J Glob Health 2011; 1:201-209.

19. Zhang Y et al. Association between C-reactive protein concentration and chronic obstructive pulmonary disease: a systematic review and meta-analysis. J Int Med Res 2012; 40:1629-1635.

20. van Durme YM et al. C-reactive protein levels, haplotypes, and the risk of incident chronic obstructive pulmonary disease. Am J Respir Crit Care Med 2009; 179:375-382.

21. Dessein PH et al. High sensitivity C-reactive protein as a disease activity marker in rheumatoid arthritis. J Rheumatol. 2004; 31 (6):1095-7.

22. Shine et al. Solid phase immunoassays for human C-reactive protein. Clin Chim Acta. 1981; 117: 13-23.

23. Chiesa C et al. C reactive protein and procalcitonin: Reference intervals for preterm and term newborns during the early neonatal period. Clin Chim Acta. 2011; 412 (11-12): 1053-9.

24. Schlebusch H et al. High Sensitive CRP and Creatinine: Reference intervals from Infancy to Childhood. J Lab Med 2002; 26 (5/6): 341-6.

Explanation of symbols / Zeichenerklärung / Explication des symboles / Spiegazione dei simboli

	English	Deutsch	Français	Italiano
	In vitro diagnostic medical device	In-vitro Diagnostikum	Dispositif médical de diagnostic <i>in vitro</i> .	Dispositivo medico diagnostico <i>in vitro</i>
	Device for near-patient testing	Produkt für patientennahe Tests	Dispositif de diagnostic près du patient	Dispositivo per analisi decentrate (near-patient testing)
	Catalogue number	Bestellnummer	Référence du catalogue	Numero di catalogo
	Batch code	Loscode	Code de lot	Codice di lotto
	Use by	Verwendbar bis	Utiliser jusqu'au	Utilizzare entro
	Shelf life after opening 3 months	Haltbarkeit nach Öffnung 3 Monate	Stabilité du flacon ouvert 3 mois	Validità del kit dopo l'apertura 3 mesi
	Shelf life after opening 6 months	Haltbarkeit nach Öffnung 6 Monate	Stabilité du flacon ouvert 6 mois	Validità del kit dopo l'apertura 6 mesi
	Temperature limitation	Temperaturbegrenzung	Limites de température	Limiti di temperatura
	Consult instructions for use	Nachschlagen in der Gebrauchsinformation	Consulter la notice d'utilisation	Consultare le istruzioni per l'uso
	Manufacturer	Hersteller	Fabricant	Fabbricante
	Sufficient for	Ausreichend für	Suffisant pour	Sufficiente per
	Do not reuse	Nicht wiederverwenden	Ne pas réutiliser	Non riutilizzare
	Reagent caps	Stopfen mit Reagenz	Bouchons réactifs	Capsule reagenti
	Buffer	Puffer	Tampon	Soluzione Tampone
	Sample collector	Probennehmer	Collecteur d'échantillons	Raccogliitore di campioni
	Ethylenediaminetetraacetic acid	Ethylendiamintetraessigsäure	Acide éthylènediaminetétraacétique	Acido etilendiamminotetraacetico
	Origin: sheep	Ursprung: Schaf	Origine: mouton	Origine: pecora
	Contains sodium azide	Enthält Natriumazid	Contient de l'azide de sodium	Contiene sodio azide
	Contains reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)	Enthält Reaction mass aus: 5-Chlor-2-methyl-2H-isothiazol-3-on [EG nr. 247-500-7] und 2-Methyl-2H-isothiazol-3-on [EG nr. 220-239-6] (3:1)	Contient un mélange de: 5-chloro-2-méthyl-2H-isothiazol-3-one [No. CE 247-500-7]; 2-méthyl-2H-isothiazol-3-one [No. CE 220-239-6] (3:1)	Contiene massa di reazione di: 5-cloro-2-metil-2H-isotiazol-3-one [EC no. 247-500-7] e 2-metil-2H-isotiazol-3-one [EC no. 220-239-6] (3:1)
	Software version	Softwareversion	Version du logiciel	Versione software
	This product fulfills the requirements of the European Parliament and of the Council on <i>in vitro</i> diagnostic medical devices	Dieses Produkt erfüllt die Anforderungen des Europäischen Parlaments und des Rates für <i>In-vitro</i> -Diagnostika	Ce produit répond aux exigences du Parlement européen et du Conseil sur les dispositifs médicaux de diagnostic <i>in vitro</i>	Questo prodotto soddisfa i requisiti imposti del Parlamento europeo e del Consiglio per i dispositivi medici diagnostici <i>in vitro</i>
	Authorised representative in Switzerland	Autorisierter Vertreter in der Schweiz	Représentant autorisé en Suisse	Rappresentante autorizzato in Svizzera

Revision history can be found from aidian.eu. / Die Revisionshistorie finden Sie auf aidian.eu. / L'historique des révisions peut être trouvé sur aidian.eu. / La cronologia delle revisioni è disponibile sul sito aidian.eu.

Aidian, QuikRead and QuikRead go are trademarks or registered trademarks of Aidian Oy.
Aidian, QuikRead und QuikRead go sind Marken oder eingetragene Marken von Aidian Oy.
Aidian, QuikRead et QuikRead go sont des marques enregistrées ou déposées d'Aidian Oy.
Aidian, QuikRead e QuikRead go sono marchi o marchi registrati di Aidian Oy.

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