

AIDIAN



New INR meter to improve your workflow

microINR[®] Expert

New INR meter to improve your workflow

Warfarin is a vitamin K antagonist (VKA) used in anticoagulant therapy to prevent thromboembolic events. The microINR Expert meter makes VKA therapy monitoring easy and efficient in healthcare settings. This meter is intended for professional use, and it can be connected to the laboratory information system ensuring accurate and fast operations with traceable results.

microINR Expert is reliable and easy to use

- Fully automatic meter with animated instructions
- Measuring range 0,8-8 INR
- Easy sample application
- Small sample volume (3 µl)
- Build in quality control system
- Results in less than 1 min
- Rechargeable battery
- Individually packaged test chips with room temperature storage

microINR Expert for smooth data transfer and traceability

- Integrated barcode scanner
- Automatic lot ID, no code chip needed
- Capacitive high resolution touch screen
- Wi-Fi and Bluetooth connections
- Ethernet connection
- POCT1-2A communication protocol

1. Chip



2. Sample



3. Result



For more detailed information, please consult the instructions for use.

Contact info

Aidian Oy

aidian@aidian.eu

+358 10 309 3000

For orders: orders@aidian.eu

Product	Prod. no	Order no
iLine microINR Expert Meter	KTP0001AA	155192
iLine microINR CHB Chips, 25 pcs	CHB0025AC	149438
iLine microINR EasyControl, 5x2 ml	EKA0005UU	144991

AIDIAN

www.aidian.eu