

Aidian – Your aid in diagnostics

Aidian is an *in vitro* diagnostics (IVD) company established in 1974 in Finland. We are dedicated to providing reliable solutions for healthcare and microbiological quality control worldwide. By combining quality and expertise, we deliver practical products that professionals can trust.

Our well-established Hygicult and Easicult product families are used by healthcare professionals, laboratories and industrial customers worldwide for routine hygiene monitoring, microbial screening and quality assurance applications. Designed to deliver reliable results with minimal effort, our dipslides offer a simple and cost-effective solution for microbiological monitoring. Their ease of use and versatility make them suitable for a wide range of applications across healthcare, food production, water testing, pharmaceutical manufacturing and industrial environments.

Besides our headquarters in Finland, we have local offices in 12 countries and collaborate with distributor partners in over 60 countries worldwide. Aidian is IVDR compliant and certified against ISO 13485:2016 Quality Management System standard, and complies with FDA quality system requirements.

References

- 1. Salo S, Laine A, Alanko T, Sjöberg A-M, Wirtanen G. Validation of the Microbiological Methods Hygicult Dipslide, Contact Plate, and Swabbing in Surface Hygiene Control: A Nordic Collaborative Study. Journal of AOAC International 2000; 83(6): 1357-1365.
- 2. NMKL 5: 2001. NordVal 018
- 3. Salo S, Alanko T, Sjöberg A-M, Wirtanen G. Validation of the Hygicult® E Dipslides Method in Surface Hygiene Control: A Nordic Collaborative Study. Journal of AOAC International 2002; 85(2): 388-394.

Contact and orders

Aidian Oy, Finland

aidian@aidian.eu / +358 10 309 3000
orders@aidian.eu

Product	Cat. no.
Hygicult TPC 10 tests	68010
Hygicult E 10 tests	68012
Hygicult E/6-GUR 10 tests	68267
Hygicult Y&F (Easicult M) 10 tests	68013
Cultura M 230 V	129997/05866
Cultura M 110/115 V	132754/06266
Cultura M 230 V (with a UK plug)	133192
Thermometer, accessory	129981/05724
Multirack, accessory	139248

8084-07EN, 10/2026

Simple surface hygiene monitoring

Hygicult



AIDIAN

Simple surface hygiene monitoring

Hygicult contact slides are intended for reliable and time saving on-site microbiological hygiene monitoring on surfaces. Hygicult slides can also be used as convenient transport media if further analysis of the samples is needed.

Available products

Hygicult TPC is for detection of bacteria and fungi on total plate count Agar for general hygiene monitoring. The method is AOAC validated¹ and NordVal certified².

Hygicult E is AOAC validated method³ for detection of faecal contaminations. Violet Red Bile Agar with added glucose (VRBG) is selective for the *Enterobacteriaceae*.

Hygicult E/B-GUR contains VRBG Agar for detecting faecal contamination, whereas β -GUR Agar is used for monitoring *Escherichia coli*.

Hygicult Y&F (Easicult M) contains malt agar for detection of fungal contaminations.

Reliable

- Established method with long experience
- AOAC validated methods^{1,3}
- Certified by NMKL/ NordVal (NordVal 018)²

Easy to use

- Three easy steps: sample-incubate-interpret
- Hinged slide for surface sampling
- On-site incubation and results interpretation
- No counting of colonies
- No microbiological skills needed

Flexible

- Long shelf-life
- Room temperature storage
- Incubation at ambient temperature or in incubator

For everyday use



Restaurants and catering services



Beverage manufacturing



Hospitals and healthcare centres



Domestic homes and day care centres



Laundries and textile services



Indoor swimming pools and spas



Gyms and sport centres



Beauty and tattoo studios



Cosmetics manufacturing



Food production plants: bakeries, dairies, slaughterhouses, fish processing factories

Only three steps to result



Dip, press or swab



Incubate



Interpret

Cultura M incubator is for easy and handy on-site culturing with temperature range of 25–45°C. Thermometer and multitrack for easy incubation of the tubes are also available.