

QuikRead go accessories.

QuikRead go system accessories make the system **more effective, more secure, and easier to use** in various testing situations.



View products and order [→](#)



QuikRead go Bag

For safely transporting the instrument and reagent kits to different locations.



QuikRead go workstation

For storing opened kit components. Suitable for all QuikRead go reagent kits.



QuikRead Cuvette Transport Box

For storing and transporting of cuvettes. The sample is stable in the cuvette for up to 2 hours in CRP/iFOBT tests and 4 hours in Strep A test.



Transportation caps

For closing cuvettes during transport.



QuikRead go sample collectors

Included in the QuikRead go easy CRP and QuikRead go HbA1c kits, also available separately.



Barcode reader

For reading patient and/or operator identification from a barcode.



QuikRead go WLAN Adapter

Connectivity through Wifi.



QuikRead go Accumulator

Internal battery for use without powercord. Sufficient for up to 50 CRP/iFOBT, 30 Strep A or 15 HbA1c measurements.

Contact and orders

Aidian Netherlands BV
+31 (0)317 35 18 38
Orders: cs@aidian.nl

Aidian Belgium BV
+31 (0)317 35 18 38
Orders: cs@aidian.be

Product	Product no.
QuikRead go Bag Torba TRM-42	OD-153902
QuikRead go workstation	OD-136206
QuikRead Cuvette Transport Box	OD-150000
Transportation caps (50 pcs)	OD-154666
QuikRead go Sample Collector 10µl (50 pcs)	OD-153211
QuikRead go Sample Collector 1µl (25 pcs)	OD-154457
Barcode reader L-46R Black USB with stand	JNC-14008
QuikRead go WLAN Adapter	OD-148045
QuikRead go Accumulator	OD-136205