

EntericBio realtime

- 🕒 Test Direct From Sample,
No DNA/RNA Extraction Or Purification
- 🕒 Freeze Dried, Room Temp Stable & Ready To Use Assays
- 🕒 Automated Setup With Pierceable Cap Technology
- 🕒 Easy To Use
- 🕒 Minimal Hands On Time
- 🕒 High Throughput

46 patient samples in under 3hrs with 30mins hands on time
276 patient samples in under 8hrs

fast, reliable & easy to use

entericbio.com

EntericBio® heatstation

EntericBio® workstation



Add Samples



Heat Treat



Automated Setup



Amplify & Detect

Revolutionising the standard test for faecal pathogens




Bacteria



Parasite



Viral



Distributie Nederland:



info@mediphos.nl
www.mediphos.com

Distribution Belgium - Luxembourg:



info@mediphos.be
www.mediphos.com Industrieweg 12-B, 6871 KA Renkum (NL), tel: +31 317 351838

Serosep Limited

Annacotty Business Park
Annacotty, Limerick,
V94FF83, Ireland

Telephone +353 (0)61 358 190
Email info@serosep.com

Serosep UK Limited

Worth Corner Business Centre
Turners Hill Road, Pound Hill
Crawley, West Sussex
RH10 7SL, UK

Telephone +44 (0) 129 388 9215