Aidian – Your aid in diagnostics

Aidian is an *in vitro* diagnostics (IVD) company established in 1974 in Finland. We are dedicated to providing reliable solutions for healthcare and microbiological quality control worldwide. By combining quality and expertise, we deliver practical products that professionals can trust.

Our well-established Hygicult and Easicult product families are used by healthcare professionals, laboratories and industrial customers worldwide for routine hygiene monitoring, microbial screening and quality assurance applications. Designed to deliver reliable results with minimal effort, our dipslides offer a simple and cost-effective solution for microbiological monitoring. Their ease of use and versatility make them suitable for a wide range of applications across healthcare, food production, water testing, pharmaceutical manufacturing and industrial environments.

Besides our headquarters in Finland, we have local offices in 12 countries and collaborate with distributor partners in over 60 countries worldwide. Aidian is IVDR compliant and certified against ISO 13485:2016 Quality Management System standard, and complies with FDA quality system requirements.

Contact and orders

Aidian Oy, Finland

aidian@aidian.eu / +358 10 309 3000
orders@aidian.eu

Product	Cat. no.
Easicult Combi 10 tests	67987
Easicult Combi PP 10 tests	05984
Hygicult Y&F (Easicult M) 10 tests	68013
Easicult TTC 10 tests	67683
Easicult TTC PP 10 tests	05988
Cultura M 230 V	129997/05866
Cultura M 110/115 V	132754/06266
Cultura M 230 V (with a UK plug)	133192
Thermometer, accessory	129981/05724
Multirack, accessory	139248

8086-07EN, 10/2026

Aidian and Easicult are registered trademarks of Aidian Oy. Aidian develops, manufactures, and markets diagnostic and microbiological testing solutions for professional use, with over 50 years of experience.

aidian.eu

Microbial control for various industries

Easicult



AIDIAN

Microbial control for various industries

Easicult is a dipslide culture method that is intended for monitoring microbial contaminations in industrial fluids. The test can be easily used on-site without any microbiological expertise and the slides can be used as convenient transport media if further analysis of the samples is needed.

Available products

Easicult Combi is for general hygiene monitoring. TTC medium is for total bacterial counts while Rose Bengal medium detects yeasts and moulds. The product with polypropylene (PP) tube offers longer stability and polystyrene (PS) tube is useful in varying temperature conditions.

Hygicult Y&F (Easicult M) contains malt agar for detection of fungal contaminations.

Easicult TTC contains TTC medium for estimation of total bacterial counts. The product is available in PP or PS tube.

Only three steps to result



Dip or swab



Incubate



Interpret

Cultura M incubator is for easy and handy on-site culturing with temperature range of 25–45°C. Thermometer and multirack for easy incubation of the tubes are also available.

Reliable

- High quality products with users all over the world
- Established method with long experience

Easy to use

- Three easy steps: sample-incubate-interpret
- No microbiological skills or special training needed
- Simple results interpretation – no counting of colonies

Flexible

- Long shelf-life
- Transportable
- Room temperature storage
- Incubation at ambient temperature or in incubator

For a wide range of applications

- Metal industry fluids
- Cooling towers, cooling waters
- Industrial cutting fluids
- Cosmetics industry
- Fuels, lubricating oils, diesel, biodiesel
- Fuel tanks: aircraft industry, ships and boats
- Paints, pigments, varnishes
- Manufacturers of cars and other motor vehicles

Simple microbial control, wherever you need it.